

Product datasheet for **SC330670**

C4ORF28 (PACRGL) (NM_001258345) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: C4ORF28 (PACRGL) (NM_001258345) Human Untagged Clone
Tag: Tag Free
Symbol: PACRGL
Synonyms: C4orf28
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330670 representing NM_001258345.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)


```
ATGCAGAAATCAGAGGGCTCTGGAGGTACACAGTTGAAAAACAGAGCAACAGGTAACATATGATCAAAGG
ACATCATCAAGCACACAGTTAAACACAGGAATGCAGTTCAGGGAAGCAAATCCTCATTGTCAACCGT
TCTCCAGAGTCTGCAAGAAAACCTCATCCTAGACCAAGTGATAAACTGAACCTAAAACAATTAATCCG
TTTGGTGAACAGTCACGAGTGCCTTCTGCATTTGCAGCTATTTACTCTAAAGGAGGTATTCCTTGCAGA
TTGGTACATGGTTCAGTAAAACACAGATTACAGTGGGAATGTCTCCTGAAAGTCTTTCATTTGATCCA
CTTCTTATTACTTTAGCTGAGGGTCTGAGAGAGACTAAGCATCCATACACTTTTGTGTCAAAGGAGGGT
TTTAGAGAATTACTTTTGGTCAAAGGTGCTCCTGAAAAAGCTATTCTTTGCTACCTAGACTGATTCCT
GTGCTAAAGGCAGCTCTGGTCCATTCGGATGATGAAGTGTGAAAGAGGATTGAATGCTCTAGTTCAG
CTAAGTGTGTTGTTGGTCTTCTCTAAACGACCATCTGAAGCATCTGCTTACAAGCCTTCCAAGAGA
TTAATGGACAAGAAATCAAAGAGCCAATCACCAGCGCATTACAAAAGCTAGAGCAACATGGTGGAAAT
GGGAGCCTTAGCATCATCAAATCTAAAATTCACATACTGCTCCATATGCTGTTGA
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001258345

Insert Size: 747 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001258345.1](#)

RefSeq Size: 2171 bp

RefSeq ORF: 747 bp

Locus ID: 133015

UniProt ID: [Q8N7B6](#)

Cytogenetics: 4p15.31

MW: 27.2 kDa