

Product datasheet for **SC330695**

UEVLD (NM_001261382) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: UEVLD (NM_001261382) Human Untagged Clone
Tag: Tag Free
Symbol: UEVLD
Synonyms: ATTP; UEV3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330695 representing NM_001261382.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGAGTTCGACTGCGAGGGCCTGAGACGGCTGCTTGGCAAGTACAAGTTCAGGGACCTAACTGTGGAA
GAACCTAAGGAATGTAATGTATTTTCCCACATTTCAAATATTCCATGGACACCTATGGTAATACATAT
AACATACCAATTCGTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCCCCTATTTGCTTCTTGAAGCCA
ACTGCAAATATGGGAATCTTAGTCGGAACATGTGGATGCTCAAGGCAGAATATATTTGCCCTATCTC
CAAACTGGAGCCATCCTAAATCTGTCTTGTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAG
GAACCTCCCATGTATTCTCTATCATCTGATGAGGCACGGCAGGTAGACTTGCTAGCCTATATTGCA
AAAATCACTGAAGGTGTTTCAGATACAAATTCAAAGAGCTGGGCAATCATGAGAATAAAACAGTCAAT
AAAATTACTGTGGTTGGAGGTGGAGAACTCGGTATTGCCTGCACATTAGCAATTTAGCAAAAGGGTATT
GCAGACAGGCTTGTCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCCACGATGGACCTTGAAATCTTC
AACCTTCTTAATGTGGAGATCAGCAAAGATTTGTCTGCCTCTGCTCATTCCAAGGTGGTATCTTCACA
GTCAACTCTTTGGGTAGTTCTCAGTCGTACCTTGATGTGGTACAGAGCAATGTGGATATGTTAGAGCC
CTTGTCACAGCTCTGGGACATTATAGTCAACACAGTGTCTGCTCGTTGCATCTCAACAGTGGAATC
ATGACCTATGTAACTGGAAGTGAAGTACATTTCTGCAATCGAGTATCGGAATGGATGTAATCTG
GATTCACAGAGATTACAGTATATTATTACAAATGTTTTGAAGGCACAGACTTCAGGCAAGAAGTATGG
GTTATTGGCGAGCAAGGAGAAGACAAAGTGTCTCAGTGGAGTGGCCAAGAAGAAGTATGAGTCATACC
TCTCAAGTGCAGCTGTCCAACAGAGCCATGGAAGTGTCTAAGAGTAAAGGTCAAAGATCCTGGTCTGTT
GGACTATCAGTAGCTGACATGGTTGACAGTATTGTAAACAATAAGAAGAAAGTGCATTCTGTATCAGCT
TTAGCAAAGGGATATTATGATATAAATAGTGAAGTGTGTTTAAAGTTTGCCTTGATCCTTGAACCAAT
GGAGTATCTGAAGTTATCAAAACCACTGAAAGAAGATACAGTTACTGAGAACTCAAAGCAGTGCA
TCCTCAATCCACAGTCTCCAACAACAGTTAAACTTTGA

Restriction Sites: SgfI-MluI
ACCN: NM_001261382
Insert Size: 1350 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001261382.1</u>
RefSeq Size:	2359 bp
RefSeq ORF:	1350 bp
Locus ID:	55293
UniProt ID:	<u>Q8IX04</u>
Cytogenetics:	11p15.1
MW:	49.8 kDa
Gene Summary:	Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (3) lacks an in-frame exon in the 5' coding region, compared to variant 1. The resulting isoform (c) lacks an internal segment, compared to isoform a. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.