

Product datasheet for **SC330696**

UEVLD (NM_001261383) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: UEVLD (NM_001261383) Human Untagged Clone
Tag: Tag Free
Symbol: UEVLD
Synonyms: ATTP; UEV3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330696 representing NM_001261383.
Blue=Insert sequence **Red**=Cloning site **Green**=Tag(s)

ATGGAGTTCGACTGCGAGGGCCTGAGACGGCTGCTTGGCAAGTACAAGTTCAGGGACCTAACTGTGGAA
GAACCTAAGGAATGTAAATGTATTTTCCCACATTTCAAATATTCCATGGACACCTATGGTAATACATAT
AACATACCAATTCGTTTCTGGATTTGGATTCTCACCTTTTCGCTCCCCCTATTTGCTTCTTGAAGCCA
ACTGCAAATATGGGAATCTTAGTCGGAACATGTGGATGCTCAAGGCAGAAATATATTTGCCCTATCTC
CAAACTGGAGCCATCTAAATCTGTCTTGTGGATTAATTAAGAAATGATTGCCAAGTTTCAAGAG
GAACTTCCCATGTATTCTCTATCATCTGATGAGGCACGGCAGGTAGACTTGCTAGCCTATATTGCA
AAAATCACTGAAGGTGTTTCAGATACAAATTCAGAGCTGGGCAATCATGAGAATAAAACAGTCAAT
AAAATTACTGTGGTTGGAGGTGGAGAACTCGGTATTGCCTGCACATTAGCAATTTAGCAAGGGTATT
GCAGACAGGCTTGTCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCCACGATGGACCTTGAAATCTTC
AACCTTCTTAATGTGGAGATCAGCAAAGATTTGTCTGCCTCTGCTCATTCCAAGGTGGTATCTTCACA
GTCAACTCTTTGGGTAGTTCTCAGTCGTACCTTGATGTGGTACAGAGCAATGTGGATATGTTAGAGCC
CTTGTCAGCTCTGGGACATTATAGTCAACACAGTGTCTGCTCGTTGCATCTCAACAGTGGAATC
ATGACCTATGTAACTGGAACCTGAGTACATTTCTGCAATCGAGTATCGGAATTGGATGTAATCTG
GATTCACAGAGATTACAGTATATTATTACAAATGTTTGAAGGCACAGACTTCAGGCAAGAAGTATGG
GTTATTGGCGAGCAAGGAGAAGACAAAGTGTCTCATGGAGTGGCCAAGAAGAAGTAGTGAGTCATACC
TCTCAAGTGCAGCTGTCCAACAGGGATATTATGATATAA

Restriction Sites: SgfI-MluI
ACCN: NM_001261383
Insert Size: 1074 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001261383.1</u>
RefSeq Size:	2235 bp
RefSeq ORF:	1074 bp
Locus ID:	55293
UniProt ID:	<u>Q8IX04</u>
Cytogenetics:	11p15.1
MW:	39.8 kDa
Gene Summary:	Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) lacks an in-frame exon in the 5' coding region and an exon in the 3' coding region, which results in a frameshift and an early stop codon, compared to variant 1. The resulting isoform (d) is shorter; it lacks an internal segment and has a distinct C-terminus, compared to isoform a. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.