

Product datasheet for SC330698

UEVLD (NM_001261385) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: UEVLD (NM_001261385) Human Untagged Clone
Tag: Tag Free
Symbol: UEVLD
Synonyms: ATTP; UEV3
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330698 representing NM_001261385.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGACACCTATGTTTTAAAGATAGTTCTCAGAAAGACCTGCTGAATTTTACTGGCACAATTCCTGTG
ATGTATCAGGTAATACATATAACATACCAATTCGTTTCTGGATTTTGGATTCTCACCTTTTCGCTCCC
CCTATTTGCTTCTTGAAGCCAAGTCAAATATGGGAATCTTAGTCGAAAAACATGTGGATGCTCAAGGC
AGAATATATTTGCCCTATCTCCAAAAGTGGAGCCATCCTAAATCTGTCTTGTGGATTAATTAAGAA
ATGATTGCCAAGTTTCAAGAGGAAGTCCCATGTATTCTCTATCATCTGATGAGGCACGGCAGGTA
GACTTGCTAGCCTATATTGCAAAAATCACTGAAGGTGTTTCAGATACAAATCAAAGAGCTGGGCAAAT
CATGAGAATAAAACAGTCAATAAAATTACTGTGGTTGGAGGTGGAGAACTCGGTATTGCCTGCACATTA
GCAATTTTCAGCAAAGGTATTGCAGACAGGCTTGCTCTTAGACCTCTCAGAAGGGACTAAAGGAGCC
ACGATGGACCTTGAAATCTTCAACCTTCTAATGTGGAGATCAGCAAAGATTTGTCTGCCTCTGCTCAT
TCCAAGGTGGTGATCTTCACAGTCAACTCTTGGGTAGTTCTCAGTCGTACCTTGATGTGGTACAGAGC
AATGTGGATATGTTTCAGAGCCCTTGTCAGCTCTGGGACATTATAGTCAACACAGTGTCTGCTCGTT
GCATCTCAACCAAGTGGAAATCATGACCTATGTAACATGGAAAGTGAAGTACATTTCTGCAAAATCGAGTG
ATCGGAATTGGATGTAATCTGGATTACAGAGATTACAGTATATTATTACAAATGTTTGAAGGCACAG
ACTTCAGGCAAAGAAGTATGGGTTATTGGCGAGCAAGGAGAAGACAAAGTGCTCACATGGAGTGGCCAA
GAAGAAGTAGTGAGTCATACCTCTCAAGTGCAGCTGTCCAACAGGGATATTATGATAAA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001261385

Insert Size: 1026 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001261385.1

RefSeq Size: 2258 bp

RefSeq ORF: 1026 bp

Locus ID: 55293

UniProt ID: Q8IX04

Cytogenetics: 11p15.1

MW: 37.6 kDa

Gene Summary: Possible negative regulator of polyubiquitination.[UniProtKB/Swiss-Prot Function]
 Transcript Variant: This variant (6) has an alternate 5' exon, which results in a downstream start codon, and also lacks an exon in the 3' coding region, which results in a frameshift and an early stop codon, compared to variant 1. The resulting isoform (f) is shorter and has a shorter N-terminus and a distinct C-terminus, compared to isoform a. Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.