

Product datasheet for **SC330701**

NOL10 (NM_001261394) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NOL10 (NM_001261394) Human Untagged Clone
Tag:	Tag Free
Symbol:	NOL10
Synonyms:	PQBP5
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330701 representing NM_001261394. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCAGGTCTCCAGCCTCAATGAGGTGAAGATTTACAGCCTCAGCTGCGGCAAGTCCCTTCTGAGTGG
CTTTCTGATAGGAAGAAGAGAGCGCTACAGAAGAAAGTGTAGATGTCCGTAGGAGAATTGAACCTATT
CAGGACTTTGAAATGCCTACTGTGTGTACCACTATTAAGGTGTCAAAGATGGACAGTACATTTAGCA
ACTGGAACATATAAACCTCGGGTTCGATGTTATGACACCTATCAATTATCCTTGAAGTTTGAAAGGTGT
TTAGATTGAGAAGTTGTACCTTTGAAATTTGTCTGATGACTACTCAAAGATTGTCTTCTACATAAT
GATAGATACATTGAATTTCAATCGCAATCAGGTTTTTACTACAAAACCAGAATACCAAAGTTTGGGAGA
GATTTCTCTTACCACTATCCATCCTGTGACTTGTACTTTGTTGGTGCAAGTTCTGAAGTTTATAGTTA
AACTTAGAACAAGGACGATACCTGAATCCTCTACAAACTGATGCTGCGGAGAATAATGTTTGTGACATA
AATTCACTGATGGCTTGTGTTGCCACAGGAACCATAGAGGGTAGAGTGGAGTGTCTGGGACCCAAGAACT
CGAAACAGAGTTGGCCTGTTAGACTGCGCCTTAAACAGTGTACAGCAGATTGAGAGATAAACAGTTTA
CCAACAATCTCTGCTTTGAAATTTAATGGTGCCTTGACCATGGCAGTTGGAACAACCACAGGGCAGGGA
AAAATATTTACTTCTTGGAGCCAGAGCATGACCTTAATGATGTTTGTCTCTACCCAACTCAGGCATG
CTTCTGACGGCCAATGAAACCCCAAGATGGGCATCTATTACATTCCAGTTTGGGTCTGCTCCTCGG
TGGTGTTCCTTCTAGACAACCTGACCGAAGAATTAGAAGAGAATCCAGAAAGCACAGTCTATGATGAT
TATAAATTTGTCAACAAGAAAGACCTTGAAAAATTTAGGGCTCACCCACCTCATTGGATCTCCTTTCCTC
CGGGCATATATGCATGGGTTTTTCATGGATATAAGACTCTATCACAAGGTGAAACTGATGGTAAATCCA
TTTGCTTATGAAGAATATAGGAAAGATAAAATACGACAGAAAAATAGAAGAAACACGTGCACAGAGAGTC
CAGTTAAAGAAATTGCCAAAAGTTAACAAAGAGCTGGCACTTAAATTAATTGAGGAAGAAGAGGAGAAG
CAGAAATCTACATGAAAAAGAAAGTTAAGAGTCTTCTAATATTCTCACCGATGATCGATTTAAAGTT
ATGTTTGAGAACCTGACTTCCAAGTAGATGAAGAGAGTGAAGAATTTAGGCTTCTGAATCCACTTGT
TCAAAAAATTAGTAAAAAAGGAAGAAGAACTAAGACTCTTAGAGCAACAAGAACTTCGTAAAAAGAA
GAGGAGGAAGAGCCGGAAGGAAAACCAAGTGATGCAGAAAGTTCGGAGAGTTCAGATGATAAAAAGCC
TGGGTTGAAGAGGTGAGGAAGCAACGCAGACTCCTCCAGCAGGAGGAAAAAGTGAAGCGGCAGGAACGA
CTCAAGGAGGACCAGCAGACAGTCTTAAAGCCCAAGTTTATGAGATCAAAGCAGGAGAAGAATTTAGA
AGCTTCAAAGATTCTGCCACAAAGCAGAACTGATGAACAAAACCTTGAAGATCGTTTGAAAAATTGAA
GCAAAAAATGGGACATTGAGTGTATCCGACACCACCGTTGGCAGCAACAATTGACATTCACGTTAAAG
AGGTCTGAACAGCAGAAGAAGCAACAGGAGGCTGAGAACTGCATCGACAAGAAAGGAAAAGACTCCGT
CGTTCGGCCGGACACCTGAAGTCAAGACACAAAAGAGGACGGTCGTTTCATGA

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001261394
Insert Size:	1917 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001261394.1</u>
RefSeq Size:	3365 bp
RefSeq ORF:	1917 bp
Locus ID:	79954
UniProt ID:	<u>Q9BSC4</u>
Cytogenetics:	2p25.1
MW:	74.5 kDa