

Product datasheet for SC330712

GNAL (NM_001261444) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	GNAL (NM_001261444) Human Untagged Clone
Tag:	Tag Free
Symbol:	GNAL
Synonyms:	DYT25
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330712 representing NM_001261444. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGTTTGTGTTGGTGGCCAGAGGGATGAGAGGAGAAAATGGATCCAGTGCTTTAACGATGTCACAGCT
ATCATTTCAGTCGCAGCCTGCAGTAGCTACAACATGGTGATTCGAGAAGATAACAACACCAACAGGCTG
AGAGAGTCCCTGGATCTTTTGAAGCATCTGGAACAACAGGTGGTTACGGACCATTCTATCATCTTG
TTCTTGAACAAACAAGATATGCTGGCAGAAAAAGTCTTGGCAGGGAAATCAAAAATTGAAGACTATTTCC
CAGAATATGCAAATTACTGTTCTGAAGACGCAACACCAGATGCAGGAGAAGATCCCAAAGTTACA
AGAGCCAAGTTCTTTATCCGGGACCTGTTTTGAGGATCAGCACGCCACCGGTGACGGCAAACATTAC
TGCTACCCGCACTTCACCTGCGCCGTGGACACAGAGAACATCCGAGGGTGTCAACGACTGCCCGGAC
ATCATCCAGCGGATGCACCTCAAGCAGTATGAGCTCTTGTA

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001261444
Insert Size:	525 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001261444.1
RefSeq Size:	5568 bp
RefSeq ORF:	525 bp
Locus ID:	2774
UniProt ID:	P38405
Cytogenetics:	18p11.21
Protein Families:	Druggable Genome
Protein Pathways:	Calcium signaling pathway, Olfactory transduction
MW:	20.5 kDa
Gene Summary:	<p>This gene encodes a stimulatory G protein alpha subunit which mediates odorant signaling in the olfactory epithelium. This protein couples dopamine type 1 receptors and adenosine A2A receptors and is widely expressed in the central nervous system. Mutations in this gene have been associated with dystonia 25 and this gene is located in a susceptibility region for bipolar disorder and schizophrenia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013]</p> <p>Transcript Variant: This variant (5) differs in the 5' UTR, lacks a large portion of the 5' coding region and initiates translation at a downstream, in-frame start codon, compared to variant 1. The encoded isoform (3) has a shorter N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>