

Product datasheet for SC330731

PTER (NM_001261837) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTER (NM_001261837) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTER
Synonyms:	HPHRP; RPR-1
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330731 representing NM_001261837. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCTTCCTTAAGTGGAAAAGTCCAAACCGTTTTGGGCCTTGTAGAGCCAAGCAAAGTGGGCCGTACC
CTGACCCATGAACACCTGGCCATGACCTTTGACTGCTGTTACTGTCCACCTCCCCGTGCCAGGAAGCT
ATTTCCAAAGAACCTATCGTGATGAAAAATTTATATTGGATTGAGAAAAACGCCTATTCCATAAAGAA
AACCTTCAATTAATCAGGAGACAGAAGCCATAAAGGAAGAACTGTTGTATTTAAAGCTAATGGTGGA
GGGGCTTTGGTGGAAAACACAACCACTGGGATTAGCCGAGACACAGACGTTGAAGAGGCTTGCAGAA
GAGACTGGCGTCCATATCATATCTGGAGCCGGGTTTTATGTGGATGCAACTCACTCCTCAGAGACCAGG
GCCATGTCAAGTGGAGCAGCTTACCGATGTCCTTATGAATGAAATTCTCCATGGAGCTGATGGAACCACT
ATCAAGTGTGGCATTATTGGAGAAATTGGTTGCTCCTGGCCTTTGACTGAGAGTGAAGAAAGGTTCTC
CAGGCCACAGCTCATGCCAGGCTCAGCTTGGTTGCTGTTATTATCCATCCTGGAGCGGAGCTCCAGG
GCACCATTTTCAAGATTATCCGAATATTGCAAGAAGCAGGCGCAGACATCTCCAAAACAGTCATGTCACAC
CTGGATAGGGTGGCTCTCCTGGTGAAGAGGGCTGTGAAGATCGAATTCTGGTAGCACATGACATACAT
ACGAAAACCCGGCTGATGAAATATGGAGGTCACGGCTATTCTCATATACTACCAATGTTGTTCTTAAA
ATGTTGCTGAGAGGCATAACTGAGAATGTGCTTGATAAGATTCTAATAGAGAACCCTAAGCAATGGCTA
ACTTTCAAATAG
  
```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001261837
Insert Size:	909 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001261837.1](#)

RefSeq Size: 3685 bp

RefSeq ORF: 909 bp

Locus ID: 9317

UniProt ID: [Q96BW5](#)

Cytogenetics: 10p13

MW: 33.4 kDa