

Product datasheet for **SC330747**

FRMD6 (NM_001267046) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FRMD6 (NM_001267046) Human Untagged Clone
Tag:	Tag Free
Symbol:	FRMD6
Synonyms:	C14orf31; c14_5320; EX1; Willin
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330747 representing NM_001267046. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAACAAATTGAATTTTCATAACAACAGAGTCATGCAAGACCGCCGAGTGTGTGCATTTTCCTTCCC
AACGATGAATCTCTGAACATCATCATAAATGTTAAGATTCTGTGTACCAGTTGCTGGTCCAGGTTTGT
GACCTGCTCAGGCTAAAGGACTGCCACCTCTTTGGACTCAGTGTATACAAAATAATGAACATGTGTAT
ATGGAGTTGTACAAAAGCTTTACAAATATTGTCCAAAAGAAATGGAAGAAAGAGGCCAGCAAGGTACGA
CAATACGAAGTCACTTGGGGTATCGACCAATTTGGGCTCCTATGATCATCCACTTCCGTGTGCAGTAC
TATGTGGAAAATGGCAGATTGATCAGTGACAGAGCAGCAAGATACTATTATTACTGGCACCTGAGAAAA
CAAGTTCTTCATTCTCAGTGTGTGCTCCGAGAGGAGGCCTACTTCTGCTGGCAGCCTTTGCCCTGCAG
GCTGATCTTGGGAACCTCAAAAGGAATAAGCACTATGGAAAATACCTCGAGCCAGAGGCTTACTTCCCA
TCTTGGGTTGTTTCCAAGAGGGGAAGGACTACATCCTGAAGCACATTCCAAACATGCACAAAGATCAG
TTTGCCTAACAGCTTCCGAAGCTCATCTTAAATATATCAAAGAGGCTGTCCGACTGGATGACGTCGCT
GTTCACTACTACAGATTGTATAAGGATAAAAGGAAAATGAAGCATCGTGACTCTTGGATTGACCATG
AGGGGAATACAGATTTTTCAGAAATTTAGATGAAGAGAAACAATTACTTTATGATTTCCCTGGACAAAT
GTTGGAATTTGGTGTGTTGTGGTAAGAAATTTGAGATTTGCCAGATGGCTTGCCTTCTGCCCGGAAG
CTCATATACTACACGGGTGCCCATGCGCTCCAGACACCTCTGCAACTTCTGAGCAACAGCCACCGC
CTCTATATGAATCTGCAGCCTGTCTGCGCCATATCCGGAAGCTGGAGGAAAACGAAGAGAAGAAGCAG
TACCGGAATCTTACATCAGTGACAACCTGGACCTCGACATGGACCAGCTGGAAAACGGTCGCGGGCC
AGCGGGAGCAGTGCGGGCAGCATGAAACACAAGCGCTGTCCCCTCATTCCACCGCCAGCCACAGCAGT
TCCACACCTCGGGCATTGAGGCAGACACCAAGCCCCGGGACACGGGGCCAGAAGACAGTACTCCAGC
AGTGCCATCCACCGCAAGCTGAAAACCTGCAGCTCAATGACCAGTCATGGCAGCTCCACACCTCAGGG
GTGGAGAGTGGCGGCAAGACCGGCTGGAAGAGGACTTACAGGACGATGAAATAGAGATGTTGGTTGAT
GACCCCCGGGATCTGGAGCAGATGAATGAAGAGTCTCTGGAAGTCAGCCAGACATGTGCATCTACATC
ACAGAGGACATGCTCATGTCGCGGAAGCTGAATGGACACTCTGGGTTGATTGTGAAAGAAATTTGGTCT
TCCACCTCGAGCTCTTCAGAAACAGTTGTTAAGCTTCGTGGCCAGAGTACTGATTCTCTCCACAGACT
ATATGTGCGAAACCAAAGACCTCCACTGATCGACACAGCTTGAGCCTCGATGACATCAGACTTTACCAG
AAAGACTTCTGCGCATTGCAGGTCTGTGTGAGGACACTGCTCAGAGTTACACCTTTGGATGTGGCCAT
GAAGTGGATGAGGAAGGCTCTATTGCAACAGTTGCTTGGCCAGCAGTGATCAACATCCAAGATGCT
TTTCCAGTCAAAAGAACAGCAAAATACTTTCTCTGGATCTCACTCATGATGAAGTTCCAGAGTTTGT
GTGTAA

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001267046
Insert Size:	1869 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001267046.1</u>
RefSeq Size:	4866 bp
RefSeq ORF:	1869 bp
Locus ID:	122786
UniProt ID:	<u>Q96NE9</u>
Cytogenetics:	14q22.1
MW:	72 kDa