

Product datasheet for **SC330810**

C10orf57 (TMEM254) (NM_001270367) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: C10orf57 (TMEM254) (NM_001270367) Human Untagged Clone
Tag: Tag Free
Symbol: C10orf57
Synonyms: bA369J21.6; C10orf57
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330810 representing NM_001270367.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGGGGCGGGACGGCTGGGGAGCTCCCAGCCAACTCTGTGTCTGGGTTTTCAAGAGGATGGTGAAG
 TCGGAGGCCAGCAATGCGGCTACCGCTATTTCCGTCCAGCAAGACCCTTGCCCTCTCTGGTCACAGTC
 CTGGGGCTGGGCTACTTCGCGTGGGTGTCTTCTGGCCTCAGAGTATCCCTATCAGAACCTTGGGCC
 CTGGGCCCTTCACTCAGTACTTGGTGGACCACCATCACACCCTCCTGTGCAATGGGTATTGGCTTGCC
 TGGCTGATTCATGTGGGAGAGTCCTTGTATGCCATAGTATTGTGAAGCATAAAGGCATCACAAGTGGT
 CGGGCTCAGCTACTCTGGTTCTACAGACTTTCTCTTTGGGATAGCGTCTCTCACCATCTTGATTGCT
 TACAAACGGAAGCGCCAAAACAACTTGA

Restriction Sites: SgfI-MluI
ACCN: NM_001270367
Insert Size: 444 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270367.1</u>
RefSeq Size:	2100 bp
RefSeq ORF:	444 bp
Locus ID:	80195
UniProt ID:	<u>Q8TBM7</u>
Cytogenetics:	10q22.3
Protein Families:	Druggable Genome, Transmembrane
MW:	16.4 kDa