

Product datasheet for SC330812

C10orf57 (TMEM254) (NM_001270369) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	C10orf57 (TMEM254) (NM_001270369) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM254
Synonyms:	bA369J21.6; C10orf57
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330812 representing NM_001270369. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGGGCGGGGACGGCTGGGGAGCTCCCAGCCAACCTCTGTGTCTGGGTTTTCAAGAGGATGGTGAAG
TCGGAGGCCAGCAATGCGGCTACCGCCTATTTCCGTCCAGCAAGACCCTTGCCCTCTCTGGTCACAGTC
CTGGGGCTGGGCTACTTCGCGTGGGTTGTCTTCTGGCCTCAGAGTATCCCTATCAGAACCTTGGGCC
CTGGGCCCTTCACTCAGTACTTGGTGGACCACCATCACACCCTCTGTGCAATGGGTATTGGCTTGC
TGGCTGATTCATGTGGGAGATCCTTGTATGCCATAGTATTGTGCAAATGA
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_001270369

Insert Size: 327 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270369.1</u>
RefSeq Size:	2152 bp
RefSeq ORF:	327 bp
Locus ID:	80195
UniProt ID:	<u>Q8TBM7</u>
Cytogenetics:	10q22.3
Protein Families:	Druggable Genome, Transmembrane
MW:	11.9 kDa