

Product datasheet for **SC330828**

LIN9 (NM_001270409) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LIN9 (NM_001270409) Human Untagged Clone
Tag:	Tag Free
Symbol:	LIN9
Synonyms:	BARA; BARPsv; Lin-9; TGS; TGS1; TGS2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330828 representing NM_001270409. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCACCGCGCGGGCAGCCTTTGAAAAAGCGCGCGGCTCGTTCAAGATGGCGGAGCTCGACCACTTG
CCTGACGAGAGCTCTTCAGCAAAAGCCCTTGTCAGTTTAAAGAAGGAAGCTTATCTAACACGTGGAAT
GAAAAGTACAGTTCTTTACAGAAAACACCTGTTTGAAAGGCAGGAATACAAGCTCTGCTGTGGAATG
AAATTTACAGCAACAATGTCAACACCAGATAAGAAAGCTTCACAGAAGATTGGTTTTCGATTACGTAAT
CTGCTCAAGCTTCTAAAGCACATAAATGGTGTATATACGAGTGGTTCTATTCAAATATAGATAAACCA
CTTTTTGAAGGTGATAATGACTTCTGTGTATGTCTAAAGGAATCTTTTCTAATTTGAAAACAAGAAAG
TTAACAAGAGTAGAATGGGGAAAAATTCGCGCGGCTTATGGGAAAACACGAGATGTTCTTCTGCATTT
TTTGAGGAAGAGAGATCAGCATTTAAACAGAAACGGCAGAAAATAAGGCTCTTACAACAAAGGAAAGTT
GCAGATGTTTCACAATTCAAAGATCTCCAGATGAAATTCCTTTGCCTCTGGTTATTGGAACGAAAGTT
ACAGCAGGATTACGTGGTGTTTCATGATGGTTTGTTCACCTGGACAAATAGATGCTGTGGTACTCTTAAT
GCTACTTATAGAGTAACTTTGTAGGACAGGGCTTGGAAACCATACCATCCCTGACTATGAAGTTCTC
AGTAATGAACCTCATGAGACAATGCCAATTGCTGCCTTTGGACAAAAACAGCGGCCTTCTCGATTTTTT
ATGACCCACACGCGTTACATTACTCCTCCTCCTCAGTCACCAATTATAGATAATGATCCTTTATTA
GGACAGTCGCGGTGGAGAAGTAAATTTCTGGCTCTGACACTGAAACATTAGGTGGTTTTCCAGTAGAA
TTTCTTATCCAAGTGACCAGATTATCAAAAATTCTCATGATTAAGGAACATATCAAGAAATTAAGG
GAAATGAACACAGAAGCAGAAAAATTGAAATCATATTCCATGCCATCAGCATTGAATTTAGCGGAGA
TATGCAACAATTGTTCTGGAGCTTGAACAGCTGAACAAGGACCTAAACAAAGTTTTGCATAAAGTTCAA
CAGTATTGCTATGAGCTTGCTCCAGACCAGGGGCTCCAGCCTGCAGATCAGCCAACAGATATGAGACGC
AGGTGTGAGGAAGAAGCACAGGAAATTTGTCGGCATGCAATTCCTCAACAGGACAGCCCTGCGTTGAA
AATGAAAATCTGACAGACTTAATTTCCAGGCTTACAGCTATTTTGTACAAATTAAGTGTCTAGCAGAA
GGAGGAGACCTGAATTCCTTTGAATTCAAATCACTTACAGACTCATTAAATGATATCAAGAGTACAATA
GACGCTTCTAATATCAGTTGCTTTCAGAAATATGTAGAAATCCATGTTGCACATATTCAGAGTGGCCTG
AGCCAGATGGGAACTTACATGCCTTTGCAGCAAATAACACCAACAGAGACTGA

```

Restriction Sites:	SgfI-MluI
ACCN:	NM_001270409
Insert Size:	1572 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270409.1</u>
RefSeq Size:	3240 bp
RefSeq ORF:	1572 bp
Locus ID:	286826
UniProt ID:	<u>Q5TKA1</u>
Cytogenetics:	1q42.12
Protein Families:	Druggable Genome
MW:	59.6 kDa
Gene Summary:	This gene encodes a tumor suppressor protein that inhibits DNA synthesis and oncogenic transformation through association with the retinoblastoma 1 protein. The encoded protein also interacts with a complex of other cell cycle regulators to repress cell cycle-dependent gene expression in non-dividing cells. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012] Transcript Variant: This variant (2) lacks an exon in the coding region but maintains the reading frame, compared to variant 1. The encoded protein (isoform 2) is shorter than isoform 1.