

Product datasheet for **SC330935**

CFAP410 (NM_001271441) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: CFAP410 (NM_001271441) Human Untagged Clone
Tag: Tag Free
Symbol: CFAP410
Synonyms: C21orf2; LRRC76; RDMS; SMDAX; YF5/A2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330935 representing NM_001271441.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```


ATGAAGCTGACGCGGAAGATGGTTCTGACCCGGGCCAAGGCCTCGGAGCTGCACAGCGTGCACAAGCTC
AACTGCTGGGGCAGCCGCTCACAGATATCTCCATTTGCCAGGAGATGCCAGCCTGGAGGTGATCACG
CTCAGTGTCAACAGCATCTCCACCCTGGAGCCTGTGAGCCGGTGCCAGCGCCTGAGTGTGCTGTACCTG
CGGAGGAACCGCATCCCCAGCCTGGCTGAGCTTCTACCTGAAGGGGCTGCCGCTGTGCGGGTGTG
TGGCTGGCCGAGAACCCGTGCTGCGGCACCAGCCCCACCGCTACCGCATGACCGTGTGCGCACCCCTG
CCGCGCTACAGAAGCTGGACAACCAGGCTGTGACGGAGGAGGAGCTGTCCCGTGCAGTGTGAGGGA
GAGGAGATCACTGCGGCCCCAGAGAGAGAGGGCACAGGCCACGGCGGCCCAAGCTATGCTGCACACTG
AGCTCCCTCAGCTCCGCTGCTGAGACTGGCCGGGACCCGCTGGACAGCGAGGAGGGAACCGGCCCC
CAGGATGAACGTGGCCTGAAGCCGCCTTCCCGGGGCCAGTTTCTTCCCTCTCAGCCAGGGATGCCTCG
AGCAGCCACAGGGGCAGGGTGTGAGTGGCGGGCCGCTAGGGGCCCGGCTGCCTCTGCCACTGCACCCAC
TGCACAGAAACCGTGGGGAGGGAGCATGGAGCCTCACAGGGCCCCGTGGGGAGGGAGCATGGAGCCTCA
CAGGGCCTTGAAGAGCTGTGCCCCAGGGGGAGCTGCGTGTGCGGGTCTGTGAATGCGCACACACGTGTA
ACACGTGCCCCGACGGAGCCGCTGGCCCCCTCAGCCTCTCCTGCTGTCTGCTGTGGAATGTGGG
CCCGGGCCCTGCTGGGCTGAGGGCAACAGGAGTACGTGGAAGAGTGCACACACACGCTCCACAGGGC
GGGCTCCTCTGCTCAGATTCTCCGAGTGTGCCAACGTCCTGACTGCCATCCTGCTGTGCTGCGGGAG
CTGGATGCAGAGGGGCTGGAGGCCGTGCAGCAGACTGTGGGCAGCCGGCTGCAGGCCCTGCTGGGGAA
GAGGTGCAGGAGCACGCCGAGTGA
  
```

Restriction Sites: SgfI-RsrII



[View online »](#)

Plasmid Map:



ACCN: NM_001271441

Insert Size: 1128 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271441.1](#)

RefSeq Size: 2611 bp

RefSeq ORF: 1128 bp

Locus ID: 755

UniProt ID: [O43822](#)

Cytogenetics: 21q22.3

MW: 40.4 kDa

Gene Summary: Four alternatively spliced transcript variants encoding four different isoforms have been found for this nuclear gene. All isoforms contain leucine-rich repeats. Three of these isoforms are mitochondrial proteins and one of them lacks the target peptide, so is not located in mitochondrion. This gene is down-regulated in Down syndrome (DS) brain, which may represent mitochondrial dysfunction in DS patients. [provided by RefSeq, Sep 2012]
Transcript Variant: This variant (3) uses an alternate in-frame splice site in the 3' coding region, compared to variant 1. The resulting isoform (3) has an additional segment in the C-terminal region, compared to isoform 1.