

Product datasheet for **SC330939**

SEL1L2 (NM_001271539) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SEL1L2 (NM_001271539) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEL1L2
Synonyms:	C20orf50; sel-1L2
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330939 representing NM_001271539. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGAAGCCCTGTCTCTGTTAATAGAGATATTGATAATTCTTGGGGTCACAATTAACCTATCAAAGCA
GAGGAACATAATAAAGACAAAAGGAAAGAAATGTCAACACACAGGTATCAGTGAACGAAATCAAACAA
TATTTATCACACATATTGGAACAAAGAACATCTAGTAATGTAATCAATAAAGAGAAAAATCTCTGGAG
AAAAAGAAGAATCAACGTAAAAAAGAATAAAGGAATTCAAAAAAGATATCTTGAAGAGAAATAAG
AATCATTTACAAAAGCAAGCAGAGAAAAATTTACAGATGAAGGAGACCAGCTATTTAAGATGGGCATC
AAGGTTCTCCAGCAGTCTAAAAGCCAAAAACAAAAGAAGAAGCCTACCTACTTTTTGCCAAAGCAGCT
GACATGGGAAACTTGAAAGCTATGGAGAAAATGGCTGACGCTTTGCTATTTGGAAATTTTGGCGTGCAA
AATATAACAGCAGCTATCCAATTATATGAGTCCTTGGCTAAAGAAGGATCATGTAAGGCCAAAAACGCA
TTAGGATTTTTGTCTTCTATGGAATAGGAATGGAATATGATCAAGCTAAGGCACTGATATATTACACC
TTTGGAAGTGCTGGAGGAAACATGATGTCCAGATGATTTGGGGTACAGATATTTGTCGGGAATCAAT
GTTCTACAGAATTGTGAAGTTGCCCTAAGTTATTACAAGAAAGTGGCAGATTATATTGCTGACACATTT
GAAAAAGTGAAGGTGTTCCAGTGGAAGAAAGTGAGACTAACGGAAGACCTGAAAATCTGAGTTCTAAC
AGTGAGATTTTGGATTGGGACATATACCAATACTATAAATTTTGGCAGAAAGAGGAGATGTTTCAGATA
CAAGTCTCTCTTGACAATTACATCTAATTGGCAGGAAAGGTCTAGATCAGGATTACTACAAAGCATT
CACTACTTCTTAAAGGCAGCAAGGCCGGGAGTGCAAAATGCCATGGCATTATAGGAAAGATGTATTTA
GAGGGGAATGCTGCCGTGCCGCAAAATAACGCTACTGCCTTCAAGTACTTTCCATGGCAGCCAGTAAG
GGCAATGCAATCGGCCTTCATGGCTTGGTCTTCTTTACTTTTCATGGAAGGAGTTCCCTGAATTAT
GCCGAAGCACTTAAATACTTTTCAGAAAGCTGCGGAAAAAGGGTGGCCGACGCACAGTTCCAGTTAGGC
TTCATGTACTACTCTGGCTCTGGAATATGGAAGGATTATAAACTTGCCTTCAAATATTTTACCTGGCA
TCTCAGAGTGGCAGCCCTTGCCATTTATTATCTGGCCAAGATGTATGCAACAGGAACAGGAGTAGTA
AGATCATGCAGAACTGCTGTGGAGAAAAGGCTAACATTCTTGAAGAGAGAAGATGTATCCAATGGCGC
TTCTCCTATGGAATCGAGCTGCCATTCAAGGACATTCACTTGGCCAGAAGATTGTACGACATGGCTGCT
CAAACGAGTCCAGATGCCACATACCTGTGCTCTTTGCCGTGATGAACTGGAACTACGCAATTTGCTC
CGGGATATCCTGTTTTTAATTTCAACAGAGATGGAAGTGGCTGAACTGGACAAACACCATTTGGACCA
CACTGGGACTTATTTGTGATTGGCTCATTGTTCTGGGCTGATTTGTGCTTAGAAATCACCATTGGG
TAG

Restriction Sites: SgfI-MluI



[View online »](#)

ACCN:	NM_001271539
Insert Size:	1728 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271539.1</u>
RefSeq Size:	1900 bp
RefSeq ORF:	1728 bp
Locus ID:	80343
UniProt ID:	<u>Q5TEA6</u>
Cytogenetics:	20p12.1
MW:	65.5 kDa