

Product datasheet for SC330955

ANKRD61 (NM_001271700) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: ANKRD61 (NM_001271700) Human Untagged Clone
Tag: Tag Free
Symbol: ANKRD61
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC330955 representing NM_001271700.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGGAATATAACCAGGAAAGGAAGCAGAGACCTGGTGGTTGACAGTGCCAAGTCCCTGGAAGATGGC
CCATCTGCAGCACTTCACTCGAACTCTATGAAGCCATCATGAGAGAAGACTGCACTACGATCGAGGTA
CTCCTGAGAAATCACCCTGTCAACCAGCCCATCACCATTCTGCCAACTCCGCCAGCAACAGATTACTT
CTGACCCAGCCGACAGAGTCTATCATCCCCATCCATCTGGCTGCCAAGTACCACAAGGCCAGAGTCTG
CTCTGCCTGTTACGGCACGGTGCTGACCCAGAAGTCAGGGACACGACAGGCCTCACCACACTCAACTTA
ATGCTACTGCAGTGGCCAGTCACTTCCACCACGTGGGCAAAACCAGGCAACAGAACGCACAGGATCCTG
ACAGACATTTCAGAATAGCAGCATCACATGTCTCCGCATCTTGTGTGCGCACGGAGCTCAAGTGAACACT
CAAGGGGAAATCAGCAACAAACGTTTACCACCTCCACCTGGCCATAGCATATGGTTGCTATCCAGTTCTC
TCCATTTTGACCCAAATGGTGCCGATGTCAATGCTATTAATGAAGCCAGCATGACACCCCTTCACATG
GCCGAAACATGCTGAATAAGGAGATGATGAAACGCTCATTGCCTATGGAGCAAACGTCAACTGTGCT
GTCTCTTCCACGGGGAACACGCCCTGAAGCTTGCACTGTGCACTGCATCAAGCAAAGCAGGCCGACTC
CTCGGGGCGGGGTCAGCTGCATCCGTCTGCTACTCACTACGGAGCCAAAGTCAACGCCAGGACTAC
AAGGGCCAAACAGCCATCCATGAGGCATGCTTTGGAGGCAGAGAGCAATCATCAATCTCCTGCTTGAA
TTTGAAGCAAATGTTAACATTTTAAACAAGAAACGGGGAATCTCCAATTTATATGTACCTTCAGCGCAGT
TGCAATGTAAGAGATACGGCACTTCTGGCCAGGCTACTTTATCACAATTATCCTCTGAGAATGACCAAT
AACCAAGGAATTCTACCTGCAGGAATCATGCTACCAGAATTCCGCCTCTTAAGGGACACCCTAATAAAG
CAATCGCAAAAACCTTTTACCTACAGGGTATCTGCAAAAGAAACATCAGGAATATTTATGGTGAGAAA
TACAAACAGCACTTGAAGCAATTCCTCCAGTGACAATATGGAATTCTGTCTACTGCTGTTATGACTTG
GCATATACCTCTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001271700
Insert Size: 1257 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001271700.1</u>
RefSeq Size:	1428 bp
RefSeq ORF:	1257 bp
Locus ID:	100310846
UniProt ID:	<u>A6NGH8</u>
Cytogenetics:	7p22.1
MW:	46.1 kDa