

Product datasheet for **SC330965**

MAGEB5 (NM_001271752) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAGEB5 (NM_001271752) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAGEB5
Synonyms:	CT3.3; MAGE-B5
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC330965 representing NM_001271752. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGACTTCTGCAGGTGTTTTTAATGCAGGATCTGACGAAAGGGCTAACAGTAGAGATGAGGAGTACCCA
TGTTCTCAGAGGTCTCACCCTCCACTGAGAGTTCATGCAGCAATTTCAAAATATTAAGGTGGGTTTG
TTGGAGCAGTTCCTGCTCTACAAGTTCAAATGAAACAGCGTATTTGAAGGAAGATATGCTGAAGATT
GTCAACCAAGATACCAAAACCAGTTTGTGAGATTACAGAAGAGCTTCTGAGCACATTGAGGTTGTC
TTTGCAGTTGACTTGAAGGAAGTCAACCAACTTGTCATTATATGACCTTGTGAGCAAGCTGAAACTC
CCCAACAATGGGAGGATTGATGTTGGCAAAGTGTACCAAGACTGGTCTCCTCATGACTTTCCTGGTT
GTGATCTTCTGAAAGGCAACTGTGCCAACAAGGAAGATACCTGGAAATTTCTGGATATGATGCAAATA
TATGATGGGAAGAAGTACTACATCTATGGAGAGCCAGGAAGCTCATCACTCAGGATTTCTGAGGCTA
ACGTACCTGGAGTACCACCAGGTGCCCTGCAGTTATCCTGCACACTATCAATTCCTTTGGGGTCCAAGA
GCCTATACTGAAACCAGCAAGATGAAAGTCTGGAATATTTGGCCAAGGTCAATGATATTGCTCCAGGT
GCCTTCTCATACAATATGAAGAGGCTTTGCAAGATGAGGAAGAGAGCCCAAGCCAGAGATGCAGCCGA
AACTGGCACTACTGCAGTGGCCAAGACTGTCTCAGGGCGAAGTTCAGCAGCTTCTCTCAACCTATTGA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001271752
Insert Size:	828 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001271752.1

RefSeq Size: 1216 bp

RefSeq ORF: 828 bp

Locus ID: 347541

UniProt ID: Q9BZ81

Cytogenetics: Xp21.3

MW: 31.9 kDa