

Product datasheet for SC331005

SEC11A (NM_001271919) Human Untagged Clone


Product data:

Product Type:	Expression Plasmids
Product Name:	SEC11A (NM_001271919) Human Untagged Clone
Tag:	Tag Free
Symbol:	SEC11A
Synonyms:	1810012E07Rik; SEC11L1; sid2895; SPC18; SPCS4A
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331005 representing NM_001271919. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGCTGTCTCTAGACTTTTTGGACGATGTGCGCGGATGAACAAGCGGCAGCTCTATTATCAAGTCCTA
AATTTTGGAAATGATTGTCTCATCGGCACTAATGATCTGGAAGGGTTAATGGTAATAACTGGAAGTGAA
AGTCCGATTGTAGTGGTCTCAGTGGCAGCATGGAACCTGCATTCATAGAGGAGATCTTCTTTCTA
ACAAATCGAGTTGAAGATCCCATACGAGTGGGAGAAAATTGTTGTTTTAGGATAGAAGGAAGAGAGATT
CCTATAGTTCACCGAGTCTTGAAGATTCATGAAAAATTTGTTCTTATATTGGAATTGTGACGATCCTC
ATGAATGACTATCCTAAATTTAAGTATGCAGTTCTTTTTGCTGGGTTATTCGTGCTGGTTCATCGT
GAGTAA

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001271919
Insert Size:	420 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271919.1
RefSeq Size:	1296 bp
RefSeq ORF:	420 bp
Locus ID:	23478
UniProt ID:	P67812
Cytogenetics:	15q25.2-q25.3
Protein Families:	Druggable Genome, Protease, Transmembrane
Protein Pathways:	Protein export
MW:	16.1 kDa
Gene Summary:	<p>This gene encodes a member of the peptidase S26B family. The encoded protein is an 18kDa subunit of the signal peptidase complex and has been linked to cell migration and invasion, gastric cancer and lymph node metastasis. Alternative splicing results in multiple transcript variants. A related pseudogene has been identified on chromosome 8. [provided by RefSeq, Dec 2012]</p> <p>Transcript Variant: This variant (6) uses an alternate in-frame exon and lacks an in-frame exon in the central coding region, compared to variant 1. The encoded isoform (6) is longer than isoform 1.</p>