

Product datasheet for **SC331077**

MFF (NM_001277061) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: MFF (NM_001277061) Human Untagged Clone
Tag: Tag Free
Symbol: MFF
Synonyms: C2orf33; EMPF2; GL004
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331077 representing NM_001277061.
Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
ATGAGTAAAGGAACAAGCAGTGACACATCACTAGGAAGGGTGAGCAGGGCAGCATTTCTTCTCCCACT
GCTGCTGAGATGGCAGAAATTAGTCGAATTCAGTACGAAATGGAATATACTGAAGGCATTAGTCAGCGA
ATGAGGGTCCCAGAAAAGTTAAAAGTAGCACCGCCAAACGCTGACCTGGAACAAGGATCCAAGAAGGA
GTTCCAAATGCTAGTGTGATAATGCAAGTCCGGAGAGGATTGTTGTAGCAGGAAATAATGAAGATGTT
TCATTTTCAAGACCAGCAGATCTTGACCTTATTCAGTCAACTCCCTTTAAACCCCTGGCACTGAAAAACA
CCACCTCGTGTACTTACGCTGAGTGAAAGACCACTAGATTTTCTGGATTTAGAAAGACCTCTACAACC
CCTCAAAATGAAGAAATCCGAGCAGTTGGCAGACTAAAAAGAGAGCGGTCTATGAGTGAAATGCTGTT
CGCCAAAATGGACAGCTGGTCAGAAATGATTCTCTGTGGCACAGATCAGATTCTGCCCAAGAAATAAA
ATTTCAAGGTTCCAGGCACCGATTTCTGCACCGGAGTACACTGTGACACCATCGCCACAACAGGCTCGG
GTCTGTCTCCCATATGTTACCTGAAGATGGAGCTAATCTTCTCTGCTCGTGGCATTGTTGTCGCTT
ATCCAGTCTTCTACTCGTAGGGCATACCAGCAGATCTGGATGTGCTGGATGAAATCGCAGACCTGTG
TTGCGTGGTGGTCTGCTGCCGCCACTTCTAATCCTCATCATGACAACGTCAGGTATGGCATTTCAAAT
ATAGATAACAACATTGAAGGAACGTCAGATGACCTGACTGTTGTAGATGCAGCTTCACTAAGACGACAG
ATAATCAAACCTAAATAGACGCTACAACCTCTGGAAGAGGAGAACAAGAACGTGCTAAAAGAGAAATG
GTCATGTATTCAATTACTGTAGCTTTCTGGCTGCTTAATAGCTGGCTCTGGTTTCGCCGCTAG
```

Restriction Sites: Sgfl-Mlul



[View online »](#)

Plasmid Map:



ACCN: NM_001277061

Insert Size: 1029 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001277061.1](#)

RefSeq Size: 2277 bp

RefSeq ORF: 1029 bp

Locus ID: 56947

UniProt ID: [Q9GZY8](#)

Cytogenetics: 2q36.3

Protein Families: Transmembrane

MW: 38.5 kDa

Gene Summary: This is a nuclear gene encoding a protein that functions in mitochondrial and peroxisomal fission. The encoded protein recruits dynamin-1-like protein (DNM1L) to mitochondria. There are multiple pseudogenes for this gene on chromosomes 1, 5, and X. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (a). Both variants 1 and 2 encode the same isoform (a).