

Product datasheet for **SC331148**

SLC33A1 (NM_001190992) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLC33A1 (NM_001190992) Human Untagged Clone
Tag: Tag Free
Symbol: SLC33A1
Synonyms: ACATN; AT-1; AT1; CCHLND; SPG42
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331148 representing NM_001190992.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCACCCACCATCTCCACAAGGACAGCAGCCGGCAACGGCGGCCAGGGAATTCAGTCACTCTCG
GATATGAAGAGCGGTCCCTGCCGCCAGCGGTTGGGATGACAGTCATTTGGACTCAGCGGGCCGGAA
GGGACAGAGAAGCTCTTCTGGGGATACCGGCACTGGCGACTTCTTAAAAGCCCCACAGAGCTCCGG
GCCGAATAAGCAGCATTTTGTACTACTCTTTCTTTACGTGCTTCAGGGTATCCCTGGGCTTGGCG
GGAAGCATCCCACTCATTTTCAAAGCAAAAATGTTAGCTATACAGACCAAGCTTTCTTCAGTTTGTG
TTTTGGCCCTTCAGTCTCAAATTAATCTGGGCCCGTTGGTTGATGCGGTCTACGTTAAGAACTTCGGT
CGTCGCAAATCTGGCTTGTCCGACACAGTATACTAGGACTTTCATGATCTATTTATCCACTCAG
GTGGACCGTTTGCTTGGGAATACCGATGACAGAACACCCGACGTGATTGCTCTCACTGTGGCGTTCTT
TTGTTTGAATCTTGGCCGCCACTCAGGACATTGCCGTGATGGTTGGCGTTAACTATGTTATCCAGG
GAAAATGTGGGTTATGCTTCTACTTGCAATTCGGTGGGCCAAACAGCGGTTACTTTTTGGGCAATGTT
TTGTTTTGGCCCTTGAATCTGCCGACTTTTGTAAACAATTTGCGGTTTCAGCCTCAACCCAGAGGA
ATCGTACTCTTTCAGATTTCTTTTTTCTGGGAACTGTATTTTTAATAACAACAACATTGGTTGCC
CTTCTGAAAAAAGAAAACGAAGTATCAGTAGTAAAAGAAGAAACACAAGGGATCACAGATACTTACAAG
CTGCTTTTTGCAATTATAAAAATGCCAGCAGTTCTGACATTTTGCCTTCTGATTCTAACTGCAAAGATT
GGTTTTTTCAGCAGCAGATGCTGTAAACAGGACTGAAATTGGTAGAAGAGGGAGTACCCAAAGAACATTTA
GCCTTATTGGCAGTTCCAATGGTTCCTTTGCAGATAAATACTGCCTCTGATTATCAGCAAATACACTGCA
GGTCCCGACGCATTAAACACATTTTACAAGCCATGCCCTACAGATTATTGCTTGGGTTAGAATATGCC
CTACTGGTTTGGTGGACTCCTAAAGTAGAACATCAAGGGGATTCCCTATATATTACTATATCGTAGTC
CTGCTGAGTTATGCTTTACATCAGGTTACAGTGTACAGCATGTATGTTTCTATAATGGCTTTCAATGCA
AAGGTTAGTGATCCACTTATTGGAGGAACATACATGACCCTTTTAAATACCGTGTCCAATCTGGGAGGA
AACTGGCCTTCTACAGTAGCTCTTTGGCTGTAGATCCCCTCACAGTAAAAGAGTGTGTAGGAGCATCA
AACCAGAATTGTGAACACCTGATGCTGTTGAGCTTTGCAAAAACCTGGGTGGCTCATGTGTTACAGCC
CTGGATGGTTATTATGTGGAGTCCATTATTTGTGTTTTTATTGGATTTGGTTGGTGGTCTTTCTTGGT
CCAAAATTTAAAAGTTACAGGATGAAGGATCATCTTCGTGAAAATGCAAAAGGAACAATTAA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001190992



[View online »](#)

Insert Size:	1650 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001190992.1</u>
RefSeq Size:	2449 bp
RefSeq ORF:	1650 bp
Locus ID:	9197
UniProt ID:	<u>O00400</u>
Cytogenetics:	3q25.31
Protein Families:	Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - ganglio series, Metabolic pathways
MW:	60.9 kDa
Gene Summary:	<p>The protein encoded by this gene is required for the formation of O-acetylated (Ac) gangliosides. The encoded protein is predicted to contain 6 to 10 transmembrane domains, and a leucine zipper motif in transmembrane domain III. Defects in this gene have been reported to cause spastic paraplegia autosomal dominant type 42 (SPG42) in one Chinese family, but not in similar patients of European descent. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2010]</p> <p>Transcript Variant: This variant (2) differs in the 3' UTR compared to variant 1. Variants 1 and 2 both encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>