

Product datasheet for **SC331171**

Matrin 3 (MATR3) (NM_001194955) Human Untagged Clone

Product data:

| | |
|---------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Matrin 3 (MATR3) (NM_001194955) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Matrin 3 |
| Synonyms: | ALS21; MPD2; VCPDM |
| Vector: | pCMV6-Entry (PS100001) |



[View online »](#)

Fully Sequenced ORF: >SC331171 representing NM_001194955.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCCAAGTCATTCCAGCAGTCATCTCTCAGTAGGGACTCACAGGGTCATGGGCGTGACCTGTCTGCG
GCAGGAATAGGCCTTCTTGCTGCTGCTACCCAGTCTTTAAGTATGCCAGCATCTTTGGAAGGATGAAC
CAGGGTACTGCACGCCTTGCTAGTTTAAATGAATCTTGAATGAGTCTTCATTGAATCAACAGGGAGCT
CATAGTGCACTGTCTTCTGCTAGTACTTCCCATAATTTGCAGTCTATATTTAACATTTGAAGTAGA
GGTCCACTCCCTTATCTTCTCAACACCCGTGGAGATGCAGACCAGGCCAGTAACATTTTGCCAGCTTT
GGTCTGTCTGCTAGAGACTTAGATGAACTGAGTCGTTATCCAGAGGACAAGATTACTCCTGAGAATTTG
CCCCAAATCCTTCTACAGCTTAAAAGGAGGAGAACTGAAGAAGGCCCTACCTTGAGTTATGGTAGAGAT
GGCAGATCTGCTACACGGGAGCCACCATACAGAGTACCTAGGGATGATTGGGAAGAAAAAGGCACTTT
AGAAGAGATAGTTTTGATGATCGTGGTCTAGTCTCAACCCAGTCTTGATTATGACCATGGAAGTCGT
TCTCAAGAATCTGGTTATTATGACAGAATGGATTATGAAGATGACAGATTAAGAGATGGAGAAAGGTGT
AGGGATGATTCTTTTTGGTGAGACCTCGCATAACTATCATAAATTTGACAGTGAATGAGAGAAATG
GGACGTGGTCTGGCCCTTACAAGAGAGATCTCTTTGAGAAAAAGAGAGGCGCTCCTCCAAGTAGC
AATATTGAAGACTCCATGGACTTTACCGAAGGGTTATCCCATCTGTGCTCTATATGTGATTTGCCA
GTTCAATTAATAAGGAGTGGAGTCAACATATCAATGGAGCAAGTCACAGTCGTCGATGCCAGCTTCTT
CTTGAAATCTACCCAGAAATGGAATCCTGACAATGATACAGGACACACAATGGGTGATCCATTGATTTG
CAGCAGTCTACAAATCCAGCACCAGGAATCTGGGACCTCCACCTCCCTCATTTCATCTTTGGGGACCA
GCAGTTGGACCAAGAGGAAATCTGGGTGCTGGAATGGAACTGCAAGGACCTAGACACATGCAGAAA
GGCAGAGTGGAACTAGCAGAGTTGTTCCACATCATGGATTTTCAACGAGGGAAAACTTGAGATACCAG
CTATTACAGCTGGTAGAACCATTTGGAGTCATTTCAAATCATCTGATTCTAAATAAAATTAATGAGGCA
TTTTATTGAAATGGCAACCACAGAGGATGCTCAGGCCGAGTGGATTATTACACAACCACCCAGCGTTA
GTATTTGGCAAGCCAGTGAGAGTTCAATTTATCCAGAAAGTATAAAAGAATAAAGAAACCTGAAGGAAAG
CCAGATCAGAAGTTTGATCAAAAAGCAAGAGCTTGACGTGTGATACATCTCAGCAATTTGCCGATTCT
GGCTATTCTGATAGTGCTGTTCTCAAGCTTGCTGAGCCTTATGGGAAAAATAAAGAATTACATATTGATG
AGGATGAAAAGTCAGGCTTTTATTGAGATGGAGACAAGAGAAGATGCAATGGCAATGGTTGACCATTGT
TTGAAAAAGCCCTTTGGTTTCAGGGGAGATGTGTGAAGTTGACCTGTCTGAGAAATATAAAAACTG
GTTCTGAGGATCCAAACAGAGGCATTGATTTACTGAAAAAGATAAATCCCGAAAAAGATCTTACTCT
CCAGATGGCAAAGATCTCCAAGTGATAAGAAATCCAAAATGATGGTTCCAGAAAGACTGAGAGTTCA
ACCGAAGGTAAAGAACAAGAAGAGAAGTCCGGTGAAGATGGTGAAGAACACAAAGGATGACCAGACA
GAGCAGGAACCTAATATGCTTCTTGAATCTGAAGATGAGCTACTTGTAGATGAAGAAGAAGCAGCAGCA
CTGCTAGAAAAGTGGCAGTTCAGTGGGAGACGAGACCGATCTTGCTAATTTAGGTGATGTGGCTTCTGAT
GGGAAAAAGGAACCATCAGATAAAGCTGTGAAAAAAGATGGAAGTCTTCCAGCAGCAGCAAGAAAAAG
CTTAAAAAGGTGGACAAGATCGAGGAACCTTGTCAAGAAAACGAAGCAGCGTTGGAAAAATGGAATTTAA
AATGAGGAAAACACAGAACCAGGTGCTGAATCTTCTGAGAACGCTGATGATCCCAACAAAGATACAAGT
GAAAACGCAGATGGTCAAAGTATGAGAACAAGGACGACTATACAATCCAGATGAGTATAGAATTGGA
CCATATCAGCCCAATGTTCTGTTGGTATAGACTATGTGATACCTAAAACAGGGTTTTACTGTAAGCTG
TGTTCACTCTTTTATACAAATGAAGAAGTTGCAAGAATACTCATTGCAGCAGCCTTCTCATTATCAG
AAATTAAGAAATTTCTGAATAAATTGGCAGAAGAACGCAGACAGAAGAAGGAACTTAA
  
```

Restriction Sites: Sgfl-Mlul

ACCN: NM_001194955

Insert Size: 2544 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

| | |
|-------------------------------|---|
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001194955.1 |
| RefSeq Size: | 5282 bp |
| RefSeq ORF: | 2544 bp |
| Locus ID: | 9782 |
| UniProt ID: | P43243 |
| Cytogenetics: | 5q31.2 |
| MW: | 94.6 kDa |
| Gene Summary: | <p>This gene encodes a nuclear matrix protein, which is proposed to stabilize certain messenger RNA species. Mutations of this gene are associated with distal myopathy 2, which often includes vocal cord and pharyngeal weakness. Alternatively spliced transcript variants, including read-through transcripts composed of the upstream small nucleolar RNA host gene 4 (non-protein coding) and matrin 3 gene sequence, have been identified. Pseudogenes of this gene are located on chromosomes 1 and X. [provided by RefSeq, Aug 2013]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. Variants 1-4 encode the same isoform (a). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |