

Product datasheet for SC331209

Ensconsin (MAP7) (NM_001198619) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ensconsin (MAP7) (NM_001198619) Human Untagged Clone
Tag: Tag Free
Symbol: Ensconsin
Synonyms: E-MAP-115; EMAP115
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331209 representing NM_001198619.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAAGGAGCCAGAAGCCAAAACAGAAGCATAACCGTTGGTCGTGGGAGGCTCTCCATGGGAGC
CCTAGCATCCACAGTGCAGATCCAGACAGGCGGTCAAGTTCCACCATGAATCTTTCGAAATATGTTGAT
CCCCTCATTAGCAAGCGGCTCTCCTCTTCATCTGCAACTTTACTAAATTCAGATAGAGCTCGCCGC
CTGCAGCTCAGCCATGGGAGAGCAGCGTTGTTAACAGACTCCTGACGCCACACATTGTTCTGGCC
AGAAGTAAAGCAGCAGCTGCCTTGTCTGGAGAAGCAGCATCTTGACGCCCATCATCATGCCCTACAAA
GCTGCACACTCTAGAAATTCGATGGATCGACAAAACCTTTGTAAACACCACCTGAGGGCTCTTCTCGC
AGGAGGATCATTATGGCACAGCGAGCTATAAAAAGAAAAGAGAGAGAGAAAATGTACTTCTCTACA
TCTGGCACCCGAAGGGCTGTATCTCCATCTAATCCCAAAGCAAGCAACCAGCTCGCTCCCGACTTTGG
CTTCCGTCCAAGTCTCTTCTCATTGTCCTGGCACACCCAGACCGACATCCTCCTTGCCACCCGGCTCA
GTCAAAGTGTCTCTGCTCAGGTCCGGCCCCATCCCCGGCAACATCCGCCCTGTCAAGAGGGAAGTC
AAAGTGGAGCCTGAGAAGAAAGTCTGAGAAGGAACCTCAGAAAGTTGCCAATGAGCCCTACTAAAG
GGCAGAGCACCTTTAGTGAAGGTAGAAGAAGCCACAGTTGAAGAGCGGACACCTGCTGAACCAGAAGTT
GGCCCTGCTGCTCCAGCCATGGCCCCAGCTCCAGCCTCGGCCAGCTCCAGCCTCGGCCAGCTCCA
GCCCGGTCCCCACCCAGCCATGGTCTCAGCCCCGTCACTCCACTGTGAATGCCAGTCTTCTGTTAAG
ACTTCTGCAGGCACCACCGACCCAGAGGAGGCCACAAGGCTTCTAGCTGAGAAGAGGCGGCTGGCCGA
GAGCAGAGAGAAAAGGAAGAAAGGGAGAGGAGGGAGCAGGAAGAGCTTGAAGACAAAAGAGAGAGGAA
TTGGCTCAACGTGTGGCTGAAGAGAGGACGACTCGCCGTGAGGAGGAGTGCAGCAGGCTGGAAGCCGAG
CAGGCCCGGAGAAAGGAGGAGCAGCTGCAGCGCAGGCGGAGGAGCGGGCGCTGCGCAGCGGGAGGAG
GCAGAGCGCGCCAGAGGCAGAAAGAAGAAGCTCGCGTTCTGTAAGAAGCAGAGAGGGTCCGGCAG
GAACGAGAGAAGCATTTCAGAGAGAAGAGCAAGAGCGCCTGGAGAGAAAGAAGCGACTTGAGGAGATT
ATGAAAAGAACCAGGAGAACAGAAGCTACAGATAAGAAAACAGTGATCAGAGAAACGGTATATAGCC
AAGGGAGCTCTACTGGAGGAACAGAGGTGTCTGCACTTCCATGTACAACAAACGCTCCGGGAAATGGA
AAGCCAGTTGGCAGCCACATGTGGTTACTCACACCAGTCAAAAAGTGACAGTGGAGAGCACTCCCGAT
TTGAAAAACAACCAATGAAAATGGTGTATCTGTTTCAGAAATGAAAATTTGAAGAAATTAATAACTTA
CCCATTGGATCTAAACATCCAGATTAGATGTCACCAACAGTGAGAGCCAGAAATTCCTTTGAATCCA
ATTTTGGCCTTTGATGATGAAGGGACACTTGGGCCCTGCCTCAGGTAGATGGTGTTCAGACACAGCAG
ACTGCAGAAGTTATATGA
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001198619
Insert Size:	1812 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001198619.1
RefSeq Size:	4248 bp
RefSeq ORF:	1812 bp
Locus ID:	9053
UniProt ID:	Q14244
Cytogenetics:	6q23.3
Protein Families:	Druggable Genome
MW:	66.9 kDa
Gene Summary:	<p>The product of this gene is a microtubule-associated protein that is predominantly expressed in cells of epithelial origin. Microtubule-associated proteins are thought to be involved in microtubule dynamics, which is essential for cell polarization and differentiation. This protein has been shown to be able to stabilize microtubules, and may serve to modulate microtubule functions. Studies of the related mouse protein also suggested an essential role in microtubule function required for spermatogenesis. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2010]</p> <p>Transcript Variant: This variant (10) has an alternate 5' exon, resulting in a downstream AUG start codon, as compared to variant 1. Variants 9 and 10 encode the same isoform (8), which is shorter at the N-terminus, as compared to isoform 1.</p>