

## Product datasheet for **SC331218**

### Leptin Receptor (LEPR) (NM\_001198687) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Leptin Receptor (LEPR) (NM_001198687) Human Untagged Clone
Tag:	Tag Free
Symbol:	Leptin Receptor
Synonyms:	CD295; LEP-R; LEPRD; OB-R; OBR
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC331218 representing NM\_001198687.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGATTTGTCAAAAATTCTGTGTGGTTTTGTTACATTGGGAATTTATTTATGTGATAACTGCGTTTAAAC
TTGTCATATCCAATTACTCCTTGGAGATTTAAGTTGTCTTGCATGCCACCAAATCAACCTATGACTAC
TTCTTTTTGCCTGCTGGACTCTCAAAGAATACTCAAATTCGAATGGACATTATGAGACAGCTGTTGAA
CCTAAGTTTAATTCAAGTGGTACTCCTTTTCTAACTTATCCAAAACAACCTTCCACTGTTGCTTTCCG
AGTGAGCAAGATAGAAACTGCTCCTTATGTGCAGACAACATTGAAGGAAAAGACATTTGTTTCAACAGTA
AATTCTTTAGTTTTTCAACAAATAGATGCAAACTGGAACATACAGTGCTGGCTAAAAGGAGACTTAAAA
TTATTCATCTGTTATGTGGAGTCATTATTTAAGAATCTATTCAGGAATTATAACTATAAGGTCCATCTT
TTATATGTTCTGCCTGAAGTGTTAGAAGATTCACCTCTGGTTCCCAAAAAGGCAGTTTTTCAGATGGTT
CACTGCAATTGCAGTGTTTCAATGTTGTGAATGTTGTGCCTGTGCCAACAGCCAACTCAACGAC
ACTCTCCTTATGTGTTTAAAAATCACATCTGGTGGAGTAAATTTCCAGTCACCTCTAATGTCAGTTCAG
CCCATAAATATGGTGAAGCCTGATCCACCATTAGGTTTGCATATGAAAATCACAGATGATGGTAATTTA
AAGATTTCTTGGTCCAGCCACCATTGGTACCATTCCACTTCAATATCAAGTGAATATTCAGAGAAT
TCTACAACAGTTATCAGAGAAGCTGACAAGATTGTCTCAGCTACATCCCTGCTAGTAGACAGTATACTT
CCTGGGTCTTCGATGAGGTTCCAGTGAGGGGCAAGAGACTGGATGGCCAGGAATCTGGAGTGACTGG
AGTACTCCTCGTGTCTTTACCACACAAGATGTCATATACTTTCCACCTAAAATTTCTGACAAGTGTGGG
TCTAATGTTTCTTTTCACTGCATCTATAAGAAGGAAAACAAGATTGTTCCCTCAAAGAGATTGTTTGG
TGGATGAATTTAGCTGAGAAAATTCCTCAAAGCCAGTATGATGTTGTGAGTGATCATGTTAGCAAAGTT
ACTTTTTCAATCTGAATGAAACCAAACCTCGAGGAAAGTTTACCTATGATGCAGTGTACTGCTGCAAT
GAACATGAATGCCATCATCGCTATGCTGAATTATATGTGATTGATGTCATATCAATATCTCATGTGAA
ACTTGGGACTTAACTAAAATGACTTGCAGATGGTCAACCAGTACAATCCAGTCACTTCCGGAAGC
ACTTTGCAATTGAGGTATCATAGGAGCAGCCTTTACTGTTCTGATATCCATCTATTCATCCCATATCT
GAGCCCAAAGATTGCTATTTGCAGAGTGATGGTTTTTATGAATGCATTTTCCAGCCAATCTTCTATTA
TCTGGCTACACAATGTGGATTAGGATCAATCACTCTCTAGGTTCACTTGACTCTCCACCAACATGTGTC
CTTCTGATTCTGTGGTGAAGCCACTGCCTCCATCCAGTGTGAAAGCAGAAATTAATAAACATTGGA
TTATTGAAAATATCTTGGGAAAAGCCAGTCTTCCAGAGAATAACCTTCAATCCAGATTCGCTATGGT
TTAAGTGGAAAAGAAGTACAATGGAAGATGTATGAGGTTTATGATGCAAAATCAAATCTGTCAGTCTC
CCAGTTCAGACTTGTGTGCAGTCTATGCTGTTCCAGTGCGCTGTAAGAGGCTAGATGGACTGGGATAT
TGGAGTAATTGGAGCAATCCAGCCTACACAGTTGTCATGGATATAAAAGTTCTATGAGAGGACCTGAA
TTTTGGAGAATAATTAATGGAGATACTATGAAAAGGAGAAAAATGTCACCTTACTTTGGAAGCCCTG
ATGAAAATGACTCATTGTGCAGTGTTCCAGAGATATGTGATAAACCATCATACTTCTGCAATGGAACA
TGGTCAGAAGATGTGGGAAATCACACGAAATTCACCTTCTGTGGACAGAGCAAGCACATACTGTTACG
GTTCTGGCCATCAATTCATTTGGTGTCTTCTGTTGCAAAATTTAATTTAACCTTTTTCATGGCCTATGAGC
AAAGTAAATATCGTGCAGTCACTCAGTGCTTATCCTTTAAACAGCAGTTGTGATTGTTTCTGGATA
CTATCACCCAGTGATTACAAGCTAATGTATTTTATTATTGAGTGGAAAAATCTTAATGAAGATGGTGAA
ATAAAATGGCTTAGAATCTTTCATCTGTTAAGAAGTATTATATCCATGATCATTTTATCCCCATTGAG
AAGTACCAGTTCAGTCTTACCCAATTTTATGGAAGGAGTGGGAAAACCAAAGATAATTAATAGTTTC
ACTCAAGATGATATTGAAAACACCAGAGTGATGCAGGTTTATATGTAATTGTGCCAGTAATTTATTC
TCTTCCATCTTATTGCTTGAACATTATTAATATCACACCAAAGAATGAAAAAGCTATTTTGGGAAGAT
GTTCCGAACCCCAAGAATTGTTCTGGGCAACAAGGACTTAATTTTCCAGAGATGCTTGAAGGCAGCATG
TTCGTTAAGAGTCATCACCCTCCCTAATCTCAAGTACCAGGGACACAAACACTGCGGAAGGCCACAG
GGTCTCTGCATAGGAAAACCAGAGACCTTTGTTCACTTGTATCTGCTGACCCCTCCCTCACTATTG
TCCTATGACCCTGCCAAATCCCCCTCTGTGAGAAACACCCAAGAATGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001198687  
**Insert Size:** 2877 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001198687.1</a>
<b>RefSeq Size:</b>	3044 bp
<b>RefSeq ORF:</b>	2877 bp
<b>Locus ID:</b>	3953
<b>UniProt ID:</b>	<a href="#">P48357</a>
<b>Cytogenetics:</b>	1p31.3
<b>Protein Families:</b>	Druggable Genome, Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Adipocytokine signaling pathway, Cytokine-cytokine receptor interaction, Jak-STAT signaling pathway, Neuroactive ligand-receptor interaction
<b>MW:</b>	109.4 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene belongs to the gp130 family of cytokine receptors that are known to stimulate gene transcription via activation of cytosolic STAT proteins. This protein is a receptor for leptin (an adipocyte-specific hormone that regulates body weight), and is involved in the regulation of fat metabolism, as well as in a novel hematopoietic pathway that is required for normal lymphopoiesis. Mutations in this gene have been associated with obesity and pituitary dysfunction. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. It is noteworthy that this gene and LEPROT gene (GeneID:54741) share the same promoter and the first 2 exons, however, encode distinct proteins (PMID:9207021).[provided by RefSeq, Nov 2010]</p> <p>Transcript Variant: This variant (4) contains alternate 5' UTR and 3' terminal exon compared to variant 1, resulting in a shorter isoform (2) with a distinct C-terminus compared to isoform 1. Variants 2 and 4 encode the same isoform.</p>