

## Product datasheet for **SC331270**

### **BAT4 (GPANK1) (NM\_001199237) Human Untagged Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** BAT4 (GPANK1) (NM\_001199237) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** BAT4  
**Synonyms:** ANKRD59; BAT4; D6S54E; G5; GPATCH10  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331270 representing NM\_001199237.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCCCGCCCTTGCTCATCACCTTCACCCAGCCACTGACCCAGCGACCTCTGGAAGGATGGGCAG
CAGCAGCCACAGCCCAGAAAGCCAGAGTCCACCCTGGATGGGGCTGCAGCCCAGCTTTCTATGAGGCC
CTGATTGGGGATGAGAGCAGCGCTCCTGACTCCCAGAGATCTCAGACTGAACCTGCCAGAGAAAGAAAG
AGAAAGAAAGAAGAATAATGAAGGCACCAGCAGCAGAAGCAGTGGCAGAAGGAGCATCAGGAAGACAT
GGACAAGGGAGATCCCTTGAGGCTGAGGATAAGATGACTCACCGATACTGAGGGCAGCCAGGAGGGG
GACCTGCCAGAACTTAGGAGACTGCTGGAACCGCATGAGGCAGGAGGAGCTGGGGGAATATCAACGCC
CGGGATGCCTTCTGGTGGACCCACTGATGTGTGCTGCTCGAGCGGGCCAGGGGGCAGCTGTGAGCTAT
CTCCTGGGCCGTGGGGCTGCCTGGTGGGGTCTGTGAGCTGAGTGGCAGGGATGCGGCTCAGCTCGCT
GAAGAAGCTGGCTTCCCTGAGGTAGCCCGCATGGTCAGGGAGAGCCATGGAGAGACAAGGAGCCCGAA
AACCGGTCTCCTACTCCCTCCCTCCAGTACTGCGAGAAGTGTGACACCCACTTCCAAGATTTCAACCAC
CGCACATCCACTGCTCACCTGCTGTCACCTGTGCGAGGGTCTCAGCCTCCAACCTTCCACTTGGGGTG
CCCATCTCCAGCCCGGGCTTCAAAGTGTGCTGAGGGGGGGCTGGGAGCCAGGAATGGGGCTGGGACCC
CGGGGTGAGGGCCGTGCCAATCCCATCCCAGTGTCTCAAGAGGGACCAGGAAGGACTAGGCTACAGA
TCAGCACCCAGCCCGAGTGACACATTTCCAGCTTGGGATACCCGAGCTGTGGCTGGGAGGGAGAGA
CCCCCTCGGGTGGCCACACTGAGCTGGAGGGAGGAGAAGGAGGGAGGAAAGACAGGGCTTGGGAG
CGGGATCTAAGGACTTACATGAACCTCGAGTTCGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001199237  
**Insert Size:** 1071 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001199237.1</a></u>
<b>RefSeq Size:</b>	2793 bp
<b>RefSeq ORF:</b>	1071 bp
<b>Locus ID:</b>	7918
<b>UniProt ID:</b>	<u><a href="#">O95872</a></u>
<b>Cytogenetics:</b>	6p21.33
<b>MW:</b>	39.3 kDa
<b>Gene Summary:</b>	<p>This gene is located in a cluster of HLA-B-associated transcripts, which is included in the human major histocompatibility complex III region. This gene encodes a protein which is thought to play a role in immunity. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Nov 2010]</p> <p>Transcript Variant: This variant (1) represents the longest transcript.</p>