

Product datasheet for **SC331298**

ZNF615 (NM_001199324) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF615 (NM_001199324) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF615
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331298 representing NM_001199324.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGATGCAGGCCAGGAATCCCTAACACTGGAGGATGTGGCTGTGGACTTCACCTGGGAGGAGTGGCAG
 TTCCTGAGCCCTGCTCAGAAGGACCTGTACCGGACGTGATGTTGGAGAACTACAGCAACCTGGTGGCA
 GTGGGGTATCAAGCCAGCAAAACAGATGCACTCTCCAAATTGGAACGAGGAGAAGAACTGCACAACA
 GAAGATGAAATCTACTCTCGAATCTGTTCTGATTGAGGAGTGCATCAGGAGGTGCATATGCAGAAATC
 AGGAAAATTGATGATCCTCTGCAGCATCACTTGCAAAATCAAAGTATTCAGAAAGAGTGTGAAACAGTGC
 CATGAACAGAATATGTTTGGAAATATTGTTAATCAGAACAAGGTCATTTCTGCTGAAGCAAGATTGT
 GATACGTTTGACTTACATGAAAAACCTTTAAAATCAAATTTAAGTTTGGAAAACAGAAAAGGAGCTCT
 GGCCTAAAGAACTCTGCTGAGTTTAATAGAGATGGGAAATCCCTTTTTCATGCTAACATAAACAATTT
 TATACTGAAATGAAGTTTCTGCAATTGCAAAACCTATTAATAAGTCCCAGTTTCATTAAGCAACAGAGA
 ACTCACAACATAGAGAATGCCATGTATGCAGTGAATGTGGGAAAGCCTTCCTCAAGTTGTCTCAGTTT
 ATTGATCATCAGAGAGTTCACACTGGAGAAAAACCTCATGTATGCAGTATGTGTGGGAAAGCTTTCTCC
 AGAAAATCCAGACTAATGGACCATCAGAGAACTCATAAGAACTGAAACATTATGAATGCACTGAATGT
 GACAAAACCTTCCTCAAGAAATCACAGCTCAATATACATCAGAAAACCTATATGGGAGGGAAACCTTAC
 ACATGTAGCCAATGTGGGAAAGCCTTCATCAAGAAGTGTGCGCTCATTTATCATCAACGAACTCATACA
 GGAGAGAAACCCCATGGATGCAGTGTATGTGGGAAAGCCTTCTCTACAAAGTTCAGTCTCACTACACAT
 CAGAAAACCTATACAGGAGAAAAACCTTATATATGTAGTGAATGTGGAAAAGGCTTCATTGAGAAGAGG
 CGTCTTACTGCACATCATCGAACTCATACTGGTGAGAAACCTTTATATGCAATAAATGTGGGAAAGGC
 TTCACCTTGAAGAAGCTTTTCAATGAAGCACTGTCTCATGGTACATCAACGAACTCATACTGGAGAGAAA
 CCTTATAAATGCAATGAGTGTGGAAAGGCTTCGCTTTGAAGAGCCCACTCATCAGACATCAGCGAACA
 CATACTGGAGAGAAACCTATGTATGCACCGAATGTGCGAAAAGGTTTCACCATGAAGAGTGACCTCATT
 GTACATCAGCGAACTCATACTGCAGAGAAGCCATATATATGCAATGATTGTGAAAAGGCTTCAGTGTG
 AAGAGCCGCTTATTGTGCATCAGCGAACTCATACTGGAGAAAAACCTATGTATGTGGTGAAGTGTGGA
 AAAGGCTTTCCAGCAAAGATCCGGCTAATGGGACATCAACGAACTCATACTGGAGAGAGAAACCTTATATT
 TGCAATGAGTGTGAAAAGGCTTCACTGAGAAGAGTCACTCAATGTACATCGGCGCACTCATACAGGA
 GAGAAACCTATGTATGCAGTGAATGTGGCAAAGGCTTAACTGGGAAAAGCATGCTCATTGCACATCAG
 CGAACTCATACTGGGAGAAACCTTATATATGCAATGAATGTGAAAAGGCTTCACCATGAAGAGTACT
 CTCAGTATACATCAGCAAACCTCATACTGGAGAGAAGCCATACAAATGCAATGAATGTGATAAACCTTC
 AGGAAGAAGACATGCCTCATACAACATCAGCGATTTACACAGGAAAGACTTCCTTTGCATGTACTGAA
 TGTGAAAATTCTTTTGCAGCAAAATGATCTTATTACACATCAGAGAATTACACAGGAGAGAAACCG
 TACAAATGCAGTGACTGCGGAAAGCCTTCACTACAAAATCAGGGCTCAATGTTTCATCAAAGAAAACAT
 ACAGGAGAGAGGCCCTATGGATGTAGTGATTGTGGGAAAGCTTTTGCACACTGTCTATCCTTGTAAAC
 CACAGGAGAATTACAGGTAG

Restriction Sites: Sgfl-MluI

ACCN: NM_001199324

Insert Size: 2229 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199324.1</u>
RefSeq Size:	4110 bp
RefSeq ORF:	2229 bp
Locus ID:	284370
UniProt ID:	<u>Q8N8J6</u>
Cytogenetics:	19q13.41
Protein Families:	Transcription Factors
MW:	84.6 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]