

Product datasheet for **SC331306**

EPB4IL2 (EPB41L2) (NM_001199388) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EPB4IL2 (EPB41L2) (NM_001199388) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPB4IL2
Synonyms:	4.1-G; 4.1G
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331306 representing NM_001199388.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGACTACTGAAGTAGGCTCTGTGTCTGAAGTGAAGAAGGACTCTAGCCAGTTAGGAACAGATGCAACC
AAGGAAAAACCTAAAGAAGTAGCAGAAAAATCAGCAGAATCAGTCTTCCGATCCAGAGGAGGAAAAAGGT
TCCCAGCCACCTCCTGCAGCTGAAAGCCAAAGTAGTCTACGCCCCAGAAGAGAGAGAAGGAAACATCG
GAGAGCAGGGGTATTTCTCGGTTTCATACGCCATGGCTTAAAGAAGCAAAAGTCATATACCTTAGTAGTG
GCCAAAGATGGAGGAGATAAAAAAGAGCCTACCCAAGCTGTTGTTGAAGAACAGGTCTTAGATAAAGAG
GAACCCCTTCCAGAAGAACAGAGACAGGCTAAGGGTGATGCTGAAGAAATGGCTCAGAAGAAACAAGAG
ATTAAAGTTGAAGTCAAGGAAGAAAAACCTCAGTGAGCAAGGAAGAAAAACCTCAGTGAGCAAAGTG
GAGATGCAGCCTACTGAATTAGTAAGTAAGGAGAGAGAAGAGAAGGTAAAAGAAACACAGGAAGACAAA
TTAGAAGGAGGAGCAGCAAAAAGGGAGACCAAGGAAGTGCAGACCAATGAGCTGAAAGCAGAGAAGGCA
TCTCAAAAAGTCACCAAGAAGACCAAACTGTCCAGTGTAAAGTGACCCTCTTAGATGGCACCGAATAC
AGCTGTGACCTGGAGAAACATGCCAAGGGACAAGTGTATTTGACAAAGTGTGTGAACACCTCAATCTC
TTGGAGAAAGACTACTTTGGACTTTTGTTCAGGAAAGCCCTGAGCAGAAAAACTGGTTAGATCCTGCT
AAAGAAATAAAGAGACAACCTGAGAAACCTTCCATGGCTATTCACTTTTAATGTGAAGTTTTATCCTCCT
GATCCTTCTCAATTGACTGAAGATATCACCAGATACTTCTTGTGCCTTCAGCTCCGGCAGGACATTGCC
TCTGGCCGCTGCCCTGCTCTTTTGTGACTCATGCTCTCCTGGGATCCTACACCTGCAGGCTGAACTT
GGTGACTATGACCCAGAAGAACATGGCAGCATCGACCTCAGTGAATTCAGTTTGGCCCTACTCAGACT
AAGGAGCTGGAAGAGAAGGTGGCAGAGCTGCACAAAACCCACAGGGGCTTATCGCCAGCACAAAGCTGAT
TCCCAGTCTTAGAAAAATGCAAAGAGGCTTTCCATGTATGGTGTGACCTACATCATGCCAAGGACTCA
GAAGGTGTGGACATCAAGCTGGGCGTGTGTCTAATGGACTTCTCATTTACAAGACAGACTGCGAATC
AATCGTTTTGCTTGGCCGAAAATCTTAAAAATTTCCATAAAACGCAGTAACTTCTACATTAAGCTCAGA
CGGCAGAGCTGGAACAGTTTGAGAGTACCATTGGATTCAAACCTGCAAAACCCAGCCGAGCGAAAAGA
CTATGGAAGTGTGCGTGGAGCATCATACTTTCTACAGGCTTGTTCCTCAGAGCAGCCACCAAAAGCC
AAGTTCCTGACCTTGGGGTCCAAATTTTCGCTATAGTGGCCGACCCAAAGCACAGACCCGCCAGGCCAGC
ACCCTCATAGATAGGCCAGCACCACTTTGAGCGCACTTCTAGTAAACGGGTCTCCAGGAGTCTAGAT
GGAGCTCCGATTGGTGTGATGGACCAAAGTCTTATGAAGGATTTTCTGGCGTCTGGGGAGATTTCA
GCCTATGGACCTGGACTGTGACGATTGCCGTGTACAAGATGGGGACGGCAGGAGGGAAGTGAGAAGC
CCAATAAAGCCACATTTGCAGCTCATTGAAGGAAAGATTACATGAGACTCTGAATATAGTGGAG
GAGAAGAAGCGGGCAGAGGTTGGGAAAGACGAAAGAGTAATCACAGAAGAAATGAATGGTAAAGAGATA
TCACCTGGGAGTGGTCTGGGGAGATTCTGAAGTGGAGCCTGTGACACAAAAGACTCCACCTCCCTG
TCTTCTGAGAGCAGCAGCAGCAGTGTGAGAGTGAAGGAGGAAGACGTGGGAGAGTACCGTCCCCACCAC
CGAGTGACCGAGGGCACCATCAGGGAGGAACAGGAGTATGAAGAAGAGGTGGAGGAAGAACCCCGCCCG
GCAGCCAAGCCACCAAGTGGTAAAAACAGAGATGGTAACAATTTCTGATGCCTCACAAAGGACAGAAATC
TCCACCAAGGAAGTCCCCATTGTCCAACTGAGACCAAAACCATCACATATGAGTCTCCACAGATTGAT
GGCGGGCTGGTGGTGTATCGGGCAGTTACTGACCGCACAAACCATCACATCTGAGTCCGTGTCAACA
ACGACAACCACACATCACCAAGACTGTAAAAGGTGGAATTTCTGAAACAAGAATTGAGAAACGCATT
GTGATCACAGGAGATGGAGATATTGATCATGACCAGGCACTGGCTCAGGCGATCAGGGAAGCCAGAGAG
CAGCACCTGACATGTCCGTACAAGAGTGGTGGTACACAAAGAAACAGAGTTGGCTGAGGAAGGGGAA
GATTAA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001199388
Insert Size: 2559 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001199388.2
RefSeq Size:	4023 bp
RefSeq ORF:	2559 bp
Locus ID:	2037
Cytogenetics:	6q23.1-q23.2
Protein Pathways:	Tight junction
MW:	95.5 kDa
Gene Summary:	Required for dynein-dynactin complex and NUMA1 recruitment at the mitotic cell cortex during anaphase (PubMed:23870127).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR, uses an alternate splice site in the coding region and lacks two consecutive coding exons, but maintains the reading frame, compared to variant 1. The encoded isoform (c) is shorter than isoform a. Variants 4 and 15-17 all encode the same isoform (c).