

## Product datasheet for **SC331307**

### **EPB4IL2 (EPB41L2) (NM\_001199389) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	EPB4IL2 (EPB41L2) (NM_001199389) Human Untagged Clone
Tag:	Tag Free
Symbol:	EPB4IL2
Synonyms:	4.1-G; 4.1G
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC331307 representing NM\_001199389.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGACTACTGAAGTAGGCTCTGTGTCTGAAGTGAAGAAGGACTCTAGCCAGTTAGGAACAGATGCAACC
AAGGAAAAACCTAAAGAAGTAGCAGAAAAATCAGCAGAATCAGTCTTCCGATCCAGAGGAGGAAAAAGGT
TCCCAGCCACCTCCTGCAGCTGAAAGCCAAAGTAGTCTACGCCCCAGAGAGAGAGAAGGAAACATCG
GAGAGCAGGGGTATTTCTCGGTTACATACGCCATGGCTTAAAGAAGCAAAAGTCATATACCTTAGTAGTG
GCCAAAGATGGAGGAGATAAAAAAGAGCCTACCCAAGCTGTTGTTGAAGAACAGGTCTTAGATAAAGAG
GAACCCCTTCCAGAAGAACAGAGACAGGCTAAGGGTGATGCTGAAGAAATGGCTCAGAAGAAACAAGAG
ATTAAAGTTGAAGTCAAGGAAGAAAAACCTCAGTGAGCAAGGAAGAAAAACCTCAGTGAGCAAAGTG
GAGATGCAGCCTACTGAATTAGTAAGTAAGGAGAGAGAAGAGAAGGTAAAAGAAACACAGGAAGACAAA
TTAGAAGGAGGAGCAGCAAAAAGGGAGACCAAGGAAGTGACACCAATGAGCTGAAAGCAGAGAAGGCA
TCTCAAAAAGTCACCAAGAAGACCAAACTGTCCAGTGTAAAGTGACCCTCTTAGATGGCACC GAATAC
AGCTGTGACCTGGAGAAACATGCCAAGGGACAAGTGTTATTTGACAAAGTGTGTGAACACCTCAATCTC
TTGGAGAAAGACTACTTTGGACTTTTGTTCAGGAAAGCCCTGAGCAGAAAAACTGGTTAGATCCTGCT
AAAGAAATAAAGAGACAACCTGAGAAACCTTCCATGGCTATTCACTTTTAATGTGAAGTTTTATCCTCCT
GATCCTTCTCAATTGACTGAAGATATCACCAGATACTTCTTGTGCCTTCAGCTCCGGCAGGACATTGCC
TCTGGCCGCTGCCCTGCTCTTTTGTGACTCATGCTCTCCTGGGATCCTACACCTGCAGGCTGAACTT
GGTGACTATGACCAGAAGAACATGGCAGCATCGACCTCAGTGAATTCAGTTTGGCCCTACTCAGACT
AAGGAGCTGGAAGAGAAGGTGGCAGAGCTGCACAAAACCCACAGGGGCTTATCGCCAGCACAAAGCTGAT
TCCCAGTCTTAGAAAAATGCAAAGAGGCTTTCCATGTATGGTGTGACCTACATCATGCCAAGGACTCA
GAAGGTGTGGACATCAAGCTGGGCGTGTGTCTAATGGACTTCTCATTTACAAGACAGACTGCGAATC
AATCGTTTTGCTTGGCCGAAAATCTTAAAAATTTCTATAAAACGCAGTAACTTCTACATTAAGGTGAGA
CGGCAGAGCTGGAACAGTTTGAGAGTACCATTGGATTCAAACCTGCAAAACCCAGGACGCGAAAAGA
CTATGGAAGTGTGCGTGGAGCATCATACTTTCTACAGGCTTGTTCCTCAGAGCAGCCACCAAAAGCC
AAGTTCCTGACCTTGGGGTCCAAATTTTCGCTATAGTGGCCGACCCAAGCACAGACCCGCCAGGCCAGC
ACCCTCATAGATAGGCCAGCACCACTTTGAGCGCACTTCTAGTAAACGGGTCTCCAGGAGTCTAGAT
GGAGCTCCGATTGGTGTGATGGACCAAGTCTTATGAAGGATTTTCTGGCGTCTGGGGAGATTTCA
GCCTATGGACCTGGACTGTGACGATTGCCGTGTACAAGATGGGACGGCAGGAGGGAAGTGAGAAGC
CCAATAAGCCACATTTGCAGCTCATTGAAGGAAAGATTACATGAGACTCTGAATATAGTGGAG
GAGAAGAAGCGGGCAGAGGTTGGGAAAGACGAAAGAGTAATCACAGAAGAAATGAATGGTAAAGAGATA
TCACCTGGGAGTGGTCTGGGGAGATTCTAAGGTGGAGCCTGTGACACAAAAGACTCCACCTCCCTG
TCTTCTGAGAGCAGCAGCAGCAGTGTGAGAGTGAAGGAGGAAAGACGTGGGAGAGTACCGTCCCCACCAC
CGAGTGACCGAGGGCACCATCAGGGAGGAACAGGAGTATGAAGAAGAGGTGGAGGAAGAACCCCGCCCG
GCAGCCAAGCCACAGTGGTAAAAACAGAGATGGTAACAATTTCTGATGCCTCACAAAGGACAGAAATC
TCCACCAAGGAAGTCCCCATTGTCCAACTGAGACCAAAACCATCACATATGAGTCTCCACAGATTGAT
GGCGGGCTGGTGGTATTTCGGGCACGTTACTGACCGCACAAACCATCACATCTGAGTCCGTGTCAACA
ACGACAACCACACATCACCAAGACTGTAAAAGGTGGAATTTCTGAAACAAGAATTGAGAAACGCATT
GTGATCACAGGAGATGGAGATATTGATCATGACCAGGCACTGGCTCAGGCGATCAGGGAAGCCAGAGAG
CAGCACCTGACATGTCCGTACAAGAGTGGTGGTACACAAAGAAACAGAGTTGGCTGAGGAAGGGGAA
GATTAA
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_001199389

**Insert Size:** 2559 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001199389.1</u>
<b>RefSeq Size:</b>	4051 bp
<b>RefSeq ORF:</b>	2559 bp
<b>Locus ID:</b>	2037
<b>Cytogenetics:</b>	6q23.1-q23.2
<b>Protein Pathways:</b>	Tight junction
<b>MW:</b>	95.5 kDa
<b>Gene Summary:</b>	Required for dynein-dynactin complex and NUMA1 recruitment at the mitotic cell cortex during anaphase (PubMed:23870127).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (5) has an alternate 5' UTR exon and lacks two consecutive in-frame exons and an in-frame segment in the CDS, compared to variant 1. The resulting isoform (c) lacks two internal segments, compared to isoform a. Variants 4 and 5 encode the same isoform.