

Product datasheet for **SC331308**

NEK1 (NM_001199397) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: NEK1 (NM_001199397) Human Untagged Clone
Tag: Tag Free
Symbol: NEK1
Synonyms: ALS24; NY-REN-55; SRPS2; SRPS2A; SRTD6
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331308 representing NM_001199397.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGAGAAGTATGTTAGACTACAGAAGATTGGAGAAGGTTCAATTTGGAAGCCATTCTTGTAAATCT
ACAGAAGATGGCAGACAGTATGTTATCAAGGAAATTAACATCTCAAGAATGTCCAGTAAAGAAAGAGAA
GAATCAAGGAGAGAAGTTGCAGTATTGGCAAACATGAAGCATCCAATATTGTCCAGTATAGAGAATCA
TTTGAAGAAAATGGCTCTCTACATAGTAATGGATTACTGTGAGGGAGGGGATCTGTTTAAAGCGAATA
AATGCTCAGAAAGGCGTTTTGTTTCAAGAGGATCAGATTTGGACTGGTTGTACAGATATGTTTGGCC
CTGAAACATGTACATGATAGAAAATTCTTCATCGAGACATTAATCTCAGAACATATTTTAACTAAA
GATGGAACAGTACAACCTGGAGATTTTGGAAATGCTAGAGTTCTTAATAGTACTGTAGAGCTGGCTCGA
ACTTGCATAGGGACCCATACTACTTGTACCTGAAATCTGTGAAAACAAACCTTACAATAATAAAAGT
GACATTTGGGCTCTGGGGTGTGTCCTTTATGAGCTGTGTACACTTAAACATGCTTTTGAAGCTGGCAGT
ATGAAAACCTGGTACTGAAGATAATATCTGGATCTTTCCACCTGTGTCTTTGCATTATTCCTATGAT
CTCCGCAGTTTGGTGTCTCAGTTATTTAAAGAAATCCTAGGGATAGACCATCAGTCACTCCATATTG
GAGAAAGGTTTTATAGCCAAACGCATTGAAAAGTTTCTCTCTCCTCAGCTTATTGCAGAAGAAATTTGT
CTAAAAACATTTTCGAAGTTTGGATCACAGCCTATACCAGCTAAAAGACCAGCTTCAGGACAAAACCTCG
ATTTCTGTTATGCCTGCTCAGAAAATTACAAAGCCTGCCGTAATATGGAATACCTTTAGCATATAAG
AAATATGGAGATAAAAAATTACACGAAAAGAAACCACTGCAAAAACATAAACAGGCCATCAAACCTCCA
GAGAAGAGAGTGAATACTGGAGAAGAAAGGAGGAAAATATCTGAGGAAGCAGCAAGAAAGAGAAGGCTG
GAATTTATTGAAAAGAAAAGAAAACAAAAGGATCAGATTATTAGTTTAAATGAAGGCTGAACAAATGAAA
AGGCAAGAAAAGGAAAGGTTGGAAAGAATAAATAGGGCCAGGGAACAAGGATGGAGAAATGTGCTAAGT
GCTGGTGAAGTGGTGAAGTAAAGGCTCCTTTTCTGGGCAGTGGAGGACTATAGCTCCATCATCTTTT
TCTTCTCGAGGACAGTATGAACATTACCATGCCATTTTGGACCAATGCAGCAACAAAGAGCAGAAGAT
AATGAAGCTAAAAGGAAAAGAGAAATATATGGTCGAGGCTTCCAGAAAGAGGAATTCGCTGGAGTT
CGTCCAGGATTTCTTATGGGGCTGCAGGTCATCACCATTTTCTGATGCTGATGATATTAGAAAACCT
TTGAAAAGATTGAAGGCGGTGTCTAAACAAGCCAATGCAAACAGGCAAAAAGGGCAGCTAGCTGTAGAA
AGAGCTAAACAAGTAGAAGAGTTCTGCAGCGAAAACGGGAAGCTATGCAGAATAAAGCTCGAGCCGAA
GGCATAATGGGAATCCTGCAAAAACCTGGCAGCTATGATGGAGGCAGGCCAGCTCTTCAAGAGGAGGG
AAGCCAAGAAAACAAGAGGAAGAGGTTTATCTGGCAAGACTGAGGCAAAATAAGACTACAGAATTTCAAT
GAGCGCCAACAGATTAAGCCAACTTCGTGGTAAAAGAAAAGAGCTAATCATTCTGAAGGACAAGAA
GGAAGTGAAGAGGCTGACATGAGGCGCAAAAAATCGAATCACTGAAGGCCCATGCAATGCACGTGCT
  
```




[View online »](#)

GCTG TACTAAA GAACA AACTAG AACGAA AGAGAA AGGAGG CTTATG AGAGAG AAAAAA AGTGTGGG AA
 GAGCAT TTTGGT GGCTAA AGGAGT TAAGATT CTGATG TTTCTCC ACCTTTGGG ACAGCATG AAACAGGT
 GGCTCTCC ATCAAAG CAACAGATG AGATCTG TTTCTG TAACTTCAG CTTTGA AAGAAGTTGG CGTG
 GACAGTAG TTTAACTG ATACCCGGG AAACCTT CAGAAG AGATGCA AAAAGACCA ACAATGCTAT TTTCAAGT
 AAGCGAGA AACTTTCG TAGATTA ATGAAA ATCTTA AGCTCA GAAGAATG AAAAGG AAAAGC AGAAT
 CTCTCTG ATACTTTT GAGATAA ATGTTCA TGAAGATG CCAAAG AGCATG AAAAAG AAAAATCAG TTTCA
 TCTGATCG CAAGAAG TGGGAG GCAGGAGT CAACTTGT GATTCC CTGGATG AGTTA AACTAG ATACA
 TCCTTCTC TACAAC TGAAGAC ATACAGT GGGAGA AGTTAT TAAATTAG GTCTA ATGGATCTCCA AGA
 AGAGCCTG GGGGAAA AGTCCG ACAGATTCT GTTCTAA AGATACTT GGAGAAG CTGAACTACA ACTTCAG
 ACAGAACT ATTAGAAA TACAAC TATTAGA AGTGAG ATTTCTCC CGAAGGGG AAAAGTACA AACCTTA
 ATTACTGG AGAAAAA AAGTACA ATGTAT TTTACATG AAATAA ACCCATCAG CTATTG TTGATTCTCCT
 GTTGAGAC AAAAAGTCC CGAGTTCAG TGAGGC ATCTCC ACAGATG TCATTG AAAC TGAAGG AAAATTTA
 GAAGAAC TGTGAT TTTGGA AACAGAA TTCTACA AGAGCCA AGTGAAC AAAACA AAGATG AGAGCTTG
 CCATGCA CTATTACTG ATGTGTGG ATTAGT GAGAAAA AAGAAA CAAGGAA ACTCAG TCGGCAG ATAGG
 ATCACC ATTCAGG AAAATGA AGTTTCTG AAGATG GAGTCT CGAGTACT GTGGACCA ACTTAGT GACATT
 CATATAG AGCCTG GAACCA ATGATTCT CAGCACTCTA AATGTG ATGTAG ATAAGTCTGTG CAACCGGAA
 CCATTTTTC CAT AAGGTGG TTTCTTCTG AACACTT GAACTTAGTCC CTCAAGT TCAATCAG TTCAAGTGT
 TCACCAGA AAGATCCT TTTGCAT TTTGATCT CACTCGC ATTTACC ACCAAAAA AAAAAACA GAATTCC
 TTGCTGAT TGGACTT TCAACTG GTCTG TTTGATG CAAACA ACCCAAAG ATGTTA AGGACATG TTCACTT
 CCAGATCTCTCAA AGCTG TTTCA GAACCTT ATGGATG TTTCC ACCCGT AGGAGATG TTTGTC AAGACAAT
 CTTGAAATAG ATGAAATG AAGATG AAAACATTA AAGAAGG ACCTTCTG ATTTCTG AAGACATTGTG TTT
 GAAGAACTG ACACAG ATTTACA AGAGCTG CAGGCCTCG ATGGAACAG TTA CTAGGGAACA ACCTGGT
 GAAGAATAC AGTGAAGA GAAGAGT CAGTCTT GAAGAACAG TGAATG TGGAGCCA ACTGCA AATGGGACA
 GATGTGGC AGATGAAG ATGACA ATCCCAGCAGT GAAAGT GCCCTG AACGAAGA ATGGCACTCAG ATAAC
 AGTGATGG TGAATTGCTAGTGA ATGTGAATG CGATAGTGTCTTTA ACCATTTAGAGGAA CTGAGACTT
 CATCTGGAGCAGGAAA TGGGCTTTGAAAA ATCTTTGAGG TTTATGAGAAAATAA AGGCTATTCATGAA
 GATGAAGATG AAAATATTGAAATTTGTTCA AAAATAGTTC AAAATATTTTGGGAAATGAACATCAGCAT
 CTTTATGCCAAGATTTCTCATTTAGT CATGGCAGATGGAGCTACCAAGAAGATAATGATGAA TAA

Restriction Sites:

SgfI-MluI

Plasmid Map:



ACCN: NM_001199397

Insert Size: 3861 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199397.1</u>
RefSeq Size:	5695 bp
RefSeq ORF:	3861 bp
Locus ID:	4750
UniProt ID:	<u>Q96PY6</u>
Cytogenetics:	4q33
Protein Families:	Druggable Genome, Protein Kinase
MW:	145.8 kDa
Gene Summary:	<p>The protein encoded by this gene is a serine/threonine kinase involved in cell cycle regulation. The encoded protein is found in a centrosomal complex with FEZ1, a neuronal protein that plays a role in axonal development. Defects in this gene are a cause of polycystic kidney disease (PKD). Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2010]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>