

Product datasheet for **SC331371**

SBK3 (NM_001199824) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SBK3 (NM_001199824) Human Untagged Clone
Tag:	Tag Free
Symbol:	SBK3
Synonyms:	SGK110
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331371 representing NM_001199824. Blue=Insert sequence Red=Cloning site Green=Tag(s)

ATGGAGCGCAGGGCCTCCGAGACCCCTGAGGATGGGGACCCAGAGGAGGACACAGCCACAGCCCTCCAA
CGGCTGGTGGAGCTGACGACCAGCAGGGTGACCCCGGTGAGGAGCCTTCGGGACCACTACACCTCATC
CGGAAGCTGGGCTCCGGCTCCTACGGCCGCGTGCTCCTTGCCAGCCTCACCAGGGGGTCCAGCTGTG
GCTCTGAAGCTCCTGCGTCGGGATTTGGTCCTGAGAAGCACCTTCTGAGGGAGTTCTGTGTGGGCCGC
TGCGTCTCTGCACACCCAGGCCCTGCTGCAGACCCCTGCGAGGACCCCTACAGACCCCGCTATTTTGCC
TTCGCCAGGAGTACGCGCCCTGTGGGACCTCAGCGGATGCTGCAGGAAAGGGCCTCCAGAACTG
CTGGTGAAGCGGGTGGTGGCCAGTTGGCAGGAGCTCTGGAATTCTCCACAGCCGGGGCTGGTCCAC
GCAGATGTCAAACCTGACAACGTGCTGGTCTTCGACCCGGTCTGCAGCCGTGTGGCCTGGGAGACCTG
GGTCTGACCCGGCCAGAGGGCAGCCCGACCCCGCCCACTGACCTCTGCCACGGCACCAGCTGAG
CTCTGTCTCTGTCTACCGCCGACACCTGCCTCTGCGGCCAGCCGTGGACTCCTGGGGCTGGGGTG
CTTCTCTCTGTGCTGCCACTGCCTGTTCCCTTGGGACGTGGCACTGGCCCCAACCTGAGTTCGAG
GCCTTCGCTGGCTGGGTGACCACCAAGCCTCAGCCACCTCAGCCACCACCACCTGGGACCAGTTTGCG
CCCCAGCCCTGGCCTTGTCCAGGGGCTTCTGGACCTGGATCCCGAGACTAGGAGCCCCCACTGGCT
GTCCTGGAATTCCTGGGGGACGACTGGGGGTTGCAAGGGAACAGAGAGGACCTGGGGTTTGGGGAGC
GCCGTGTCCTATGAGGACAGGGAGGAGGGAGGCTCAAGCCTGGAGGAGTGGACAGATGAGGGTGATGAC
AGCAAAAGTGGTGGGAGGACGGGGACAGATGGGGAGCTCCCTGA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001199824
Insert Size:	1080 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001199824.2</u>
RefSeq Size:	1270 bp
RefSeq ORF:	1080 bp
Locus ID:	100130827
UniProt ID:	<u>P0C264</u>
Cytogenetics:	19q13.42
MW:	38.5 kDa