

Product datasheet for **SC331410**

SAFB (NM_001201338) Human Untagged Clone

Product data:

| | |
|---------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | SAFB (NM_001201338) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | SAFB |
| Synonyms: | HAP; HET; SAB-B1; SAF-B; SAF-B1; SAFB1 |
| Vector: | pCMV6-Entry (PS100001) |



[View online »](#)

Fully Sequenced ORF: >SC331410 representing NM_001201338.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCGGAGACTCTGTCAGGCCTAGGTGATTCTGGAGCGGCGGGCGGGCGGCTCTGAGCTCCGCTCG
TCAGAGACCGGGACCGGGCGCCTCAGCGACCTGCGAGTGATCGATCTGCGGGCGGAGCTGAGGAAACGG
AATGTGGACTCGAGCGGCAACAAGAGCGTTTTGATGGAGCGGCTGAAGAAGGCAATTGAAGATGAAGGT
GGTAATCTGACGAAATTGAAATTACCTCCGAGGAAACAAGAAAACATCAAAGAGGTCTAGCAAAGGG
CGCAAACCAAGAAGAAGAGGGTGTGGAAGATAACGGGCTGGAGGAAAACCTTGGGGATGGACAGGAGGAT
GTTGAGACCAGTCTGGAGAACTGCAGGACATCGACATCATGGATATCAGTGTGTTGGATGAAGCAGAA
ATTGATAATGGAAGCGTTGCAGATTGTGTGCAAGACGATGATGCTGATAACCTCCAGGAGTCCCTGTGCG
GATAGTAGAGAGCTAGTCGAGGGGAAATGAAAGAGCTTCCGAGCAGCTTCAGGAAACATGCTATAGAG
GACAAAGAACTATAACAATTTAGATACTTCATCATCTGACTTCACTATATTACAGGAAATTGAAGAG
CCATCCCTGGAGCCAGAAAATGAGAAAATACTCGACATTTGGGGGAAACTTGTAATCTGAGCCAGTA
AAAGAAGAAAGTCCGAGCTGGAGCAGCATTTCACAGGACACAAGTAGCGTGGGGCCAGACAGAAAG
CTTGCGGAGGAAGAGGACCTATTTGACAGCGCCATCCGGAAGAGGGTGTATTTAGATTTGGCCAGCGAG
TCAACAGCACACGCTCAGTCGAGCAAGGACAGACGCTGTAGCGGTAGTGAAAAGGGAGCCCGCGGAG
CAGCCAGGCGATGGCGAGAGGACGGACTGTGAGCCTGTAGGGCTAGAGCCGGCAGTTGAGCAGAGTAGT
GCGGCCTCCGAGCTCGCGGAGGCTCTAGCGAGGAGCTCGCAGAAGCACCACGGAAGCCCAAGCCCA
GAAGCCAGAGATAGCAAAGAAGACGGGAGGAAGTTGATTTTGACGCTTGTAAATGAAGTCCCTCCGCT
CCTAAAGAGTCTCAACCAGTGAGGGCGCTGATCAGAAAATGAGTTCTCCCGAAGATGACTCGGATACA
AAAAGGCTTTCCAAGAGGAAAAGGGTTCGACGAGTGTGGTAGAAAATTTCTGGGTTAGTGGACTCTCT
TCTACAACCAGAGCTACAGATTTGAAGAATCTTTTCAGCAAATATGGGAAGGTGGTGGCCCAAGGTT
GTGAAAATGCCGGAGTCTGGAGCTCGCTGTTACGGTTTTGTCACGATGTCCACAGCAGAGAAGGCC
ACAAAATGCATTAACCACCTGCACAAGACGGAGCTCCACGAAAAGATGATCTCCGTGGAGAAAGCCAAA
AATGAACCTGTGGGAAAGAAAACCTCTGACAAAAGAGACAGTGACGGGAAAAGGAGAAGTCGAGCAAC
AGTGACAGATCTACAAACCTTAAGAGGGATGATAAATGTGACAGAAAAGATGATGCTAAGAAGGGTGAC
GACGGAAGTGGAGAAAAGAGTAAGGACCAAGATGATCAGAAAACCTGGCCCTCAGAGCGATCTCGAGCC
ACAAAGTCAGGAAGTTCGAGGGACCGAACGGACTGTAGTAAATGGATAAATCCAAAGGGTGCCTGTGATT
AGTGTA AAAACGTCGGGTCCAAGAGAGAGCTTCCAAAAGCCAGGATCGCAAATCAGCCAGCAGAGAG
AAGCGGTCGTCGTGTCCTTTGATAAAGTCAAGGAGCCTCGGAAGTCAAGAGACTCAGAGTCCCATAGC
AGGGTGCGTGAACGCAGTGAACGCGAACACGCATGCAGGCGCAGTGGGAGCGCAGGAGCGTGAGCGG
CTGGAGATTGCCCGAGAGAGGCTGGCCTTCCAGCGCCAGCGGCTGGAGCGGGAGCGCATGGAGCGGGAA
CGGCTGGAGCGCAACGCATGCACGTGGAGCACGAGCGCAGGCGGAGCAGGAGCGCATCCACCGTGAG
CGCGAGGAGCTGAGGCGCCAGCAGGAACTGCGCTATGAGCAGGAGCGGGCCCGCGGTGCGGGCGCCC
TACGACCTGGACCGGCGAGATGATGCCTATTGGCCGGAAGCCAAGCGGGCCCGCTGGATGAGCGCTAC
CATTCTGACTTTAACCGCCAGGACCGCTTCCACGACTTTGACCACAGGGACCGCGGCCGCTACCCCGAC
CACTCGGTGGACAGGAGAGAAGGTTCAAGGTCAATGATGGGAGAACGAGAAGGACAGCATTACCCAGAA
CGCCATGGAGGACCAGAGCGCCACGGCCGGGACTCCCGGATGGCTGGGGGGCTATGGCTCTGACAAG
AGGATGAGCGAGGGCCGGGGCTGCCTCCTCCCCAGGGGCAGACGTGACTGGGGGGACCATGGCCGA
AGAGAGGATGACCGGTCAATGGCAGGGCACGGCCGACGGGGCATGATGGACAGGGATCACAAGAGGTGG
CAAGGTGGCGAGAGAAGCATGTCCGGTCACTCCGGCCTGGCCACATGATGAACCAGGAGGAATGTCA
GGGCGCGGAGCTTTGCCCGAGGCGGGCTCCCGGGCCACCCCATCCACACGGTGGCATGCAGGGC
GGGTTTGGAGGCCAGAGCCGGGGAGCAGGCCAGCGATGCCCGTTCACTCGCCGCTACTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001201338
Insert Size: 2754 bp

| | |
|-------------------------------|--|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001201338.1 |
| RefSeq Size: | 3121 bp |
| RefSeq ORF: | 2754 bp |
| Locus ID: | 6294 |
| UniProt ID: | Q15424 |
| Cytogenetics: | 19p13.3 |
| Protein Families: | Druggable Genome, Transcription Factors |
| MW: | 102.9 kDa |
| Gene Summary: | <p>This gene encodes a DNA-binding protein which has high specificity for scaffold or matrix attachment region DNA elements (S/MAR DNA). This protein is thought to be involved in attaching the base of chromatin loops to the nuclear matrix but there is conflicting evidence as to whether this protein is a component of chromatin or a nuclear matrix protein. Scaffold attachment factors are a specific subset of nuclear matrix proteins (NMP) that specifically bind to S/MAR. The encoded protein is thought to serve as a molecular base to assemble a 'transcriptosome complex' in the vicinity of actively transcribed genes. It is involved in the regulation of heat shock protein 27 transcription, can act as an estrogen receptor co-repressor and is a candidate for breast tumorigenesis. This gene is arranged head-to-head with a similar gene whose product has the same functions. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2011]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest protein (isoform 1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |