

Product datasheet for **SC331411**

SAFB (NM_001201339) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SAFB (NM_001201339) Human Untagged Clone
Tag:	Tag Free
Symbol:	SAFB
Synonyms:	HAP; HET; SAB-B1; SAF-B; SAF-B1; SAFB1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331411 representing NM_001201339.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCGGAGACTCTGTCTAGGCCTAGGTGATTCTGGAGCGGCGGGCGGGCGGCTCTGAGCTCCGCTCG
TCAGAGACCGGGACCGGGCGCCTCAGCGACCTGCGAGTGATCGATCTGCGGGCGGAGCTGAGGAAACGG
AATGTGGACTCGAGCGGCAACAAGAGCGTTTGGATGGAGCGGCTGAAGAAGGCAATTGAAGTGAAGGT
GGTAATCTGACGAAATTGAAATTACCTCCGAGGAAACAAGAAAACATCAAAGAGGTCTAGCAAAGGG
CGCAAACCAGAAGAAGAGGGTGTGGAAGATAACGGGCTGGAGGAAAACCTGGGGATGGACAGGAGGAT
GTTGAGACCAGTCTGGAGAACTGCAGGACATCGACATCATGGATATCAGTGTGTTGGATGAAGCAGAA
ATTGATAATGGAAGCGTTGCAGATTGTGTGCAAGACGATGATGCTGATAACCTCCAGGAGTCCCTGTCTG
GATAGTAGAGAGCTAGTCGAGGGGAAATGAAAGAGCTCCGGAGCAGCTTCAGGAAACATGCTATAGAG
GACAAAGAACTATAACAATTTAGATACTTCATCATCTGACTTCACTATATTACAGGAAATTGAAGAG
CCATCCCTGGAGCCGAAATGAGAAAATACTCGACATTTGGGGGAACTTGTAAATCTGAGCCAGTA
AAAGAAGAAAGTCCGAGCTGGAGCAGCATTGACAGGACACAAGTAGCGTGGGGCCAGACAGAAAG
CTTGCGGAGGAAGAGGACCTATTTGACAGCGCCATCCGGAAGAGGGTGTATTTAGATTTGGCCAGCGAG
TCAACAGCACACGCTCAGTCGAGCAAGGACAGACGCTGTAGCGGTAGTGAAAAGGGAGCCCGCGGAG
CAGCCAGGCGATGGCGAGAGGACGGACTGTGAGCCTGTAGGGCTAGAGCCGGCAGTTGAGCAGAGTAGT
GCGGCCTCCGAGCTCGCGGAGGCTCTAGCGAGGAGCTCGCAGAAGCACCACGGAAGCCCAAGCCCA
GAAGCCAGAGATAGCAAAGAAGACGGGAGGAAGTTGATTTTACGCTTGTAAATGAAGTCCCTCCGCT
CCTAAAGAGTCTCAACCAGTGAGGGCGCTGATCAGAAAATGAGTTCTCCCGAAGATGACTCGGATACA
AAAAGGCTTTCCAAGAGGAAAAGGGTTCGACGAGTGTGGTAGAAAATTTCTGGTTAGTGGACTCTCT
TCTACAACAGAGCTACAGATTTGAAGAATCTTTTTCAGCAAATATGGGAAGGTGGTGGCCCAAGTT
GTGAAAATGCCGAGTCTGGAGCTCGCTGTTACGGTTTTGTCACGATGTCCACAGCAGAGAGGCC
ACAAAATGCATTAACCACCTGCACAAGACGAGCTCCACGAAAAGATGATCTCCGTGGAGAAAGCCAAA
AATGAACCTGTGGGAAAAGAAACCTCTGACAAAAGAGACAGTGCAGGAAAAGGAGAAAGTCGAGCAAC
AGTGACAGATCTACAAACCTTAAGAGGGATGATAAATGTGACAGAAAAGATGATGCTAAGAAGGGTAC
GACGGAAGTGGAGAAAAGAGTAAGGACCAAGATGATCAGAAAACCTGGCCCTCAGAGCGATCTCGAGCC
ACAAAGTCAGGAAGTCCGAGGACCGAACGGACTGTAGTAAATGGATAAATCCAAAGGGTGCCTGTGATT
AGTGTAAAACGTCGGGTCCAAGAGAGAGCTTCCAAAAGCCAGGATCGCAAATCAGCCAGCAGAGAG
AAGCGGTCCGTCTGCTTTGATAAAGTCAAGGAGCCTCGGAAGTCAAGAGACTCAGAGTCCCATAGG
GTGCGTGAACGAGTGAACGCGAACAACGCATGCAGGCGCAGTGGGAGCGGAGGAGCGTGAGCGGCTG
GAGATTGCCCGAGAGAGGCTGGCCTTCCAGCGCCAGCGGCTGGAGCGGGAGCGCATGGAGCGGGAACGG
CTGGAGCGCAACGCATGCACGTGGAGCACGAGCGCAGGCGGAGCAGGAGCGCATCCACCGTGAGCGC
GAGGAGCTGAGGCGCCAGCAGGAACTGCCTATGAGCAGGAGCGGCGGCCCGCGGTGCGGCGGCCCTAC
GACCTGGACCGCGGAGATGATGCCTATTGGCCGGAAGCAAGCGGGCCCGCTGGATGAGCGCTACCAT
TCTGACTTTAACCGCCAGGACCGCTTCCACGACTTTGACCACAGGGACCGCGGCCCTACCCCGACCAC
TCGGTGGACAGGAGAGAAGTTCAAGGTCAATGATGGGAGAACGAGAAGGACAGCATTACCCAGAACGC
CATGGAGGACAGAGCGCCACGGCCGGGACTCCCGCATGGCTGGGGGGGCTATGGCTCTGACAAGAGG
ATGAGCGAGGGCCGGGGGCTGCCTCCTCCCCAGGGGCAGACGTGACTGGGGGACCATGGCCGAAGA
GAGGATGACCGGTGATGGCAGGGCACGGCCGACGGGGCATGATGGACAGGGATCACAAGAGGTGGCAA
GGTGGCGAGAGAAGCATGTCGGTCACTCCGGGCTGGCCACATGATGAACCGAGGAGGAATGTCAGGG
CGCGGCAGCTTTGCCAGGCGGGGCTCCCGGGCCACCCATCCACACGGTGGCATGCAGGGCGGG
TTTGGAGGCCAGAGCCGGGGGAGCAGGCCAGCGATGCCCGTTCCTACTCGCCGCTACTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001201339
Insert Size: 2751 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001201339.1
RefSeq Size:	3118 bp
RefSeq ORF:	2751 bp
Locus ID:	6294
UniProt ID:	Q15424
Cytogenetics:	19p13.3
Protein Families:	Druggable Genome, Transcription Factors
MW:	102.8 kDa
Gene Summary:	<p>This gene encodes a DNA-binding protein which has high specificity for scaffold or matrix attachment region DNA elements (S/MAR DNA). This protein is thought to be involved in attaching the base of chromatin loops to the nuclear matrix but there is conflicting evidence as to whether this protein is a component of chromatin or a nuclear matrix protein. Scaffold attachment factors are a specific subset of nuclear matrix proteins (NMP) that specifically bind to S/MAR. The encoded protein is thought to serve as a molecular base to assemble a 'transcriptosome complex' in the vicinity of actively transcribed genes. It is involved in the regulation of heat shock protein 27 transcription, can act as an estrogen receptor co-repressor and is a candidate for breast tumorigenesis. This gene is arranged head-to-head with a similar gene whose product has the same functions. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2011]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the 3' coding region compared to variant 1. The resulting protein (isoform 2) is shorter compared to isoform 1.</p> <p>Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>