

Product datasheet for **SC331455**

ZNF816A (ZNF816) (NM_001202457) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF816A (ZNF816) (NM_001202457) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF816A
Synonyms:	ZNF816A
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331455 representing NM_001202457.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTTACGTGAGGAAGCCACCAAGAAGAGCAAAGAAAAGGAGCCAGGGATGGCTCTTCTCAGGGACGC
TTGACTTTCAGGGATGTGGCTATAGAGTTCTCTTTGGAGGAGTGGAAATGCCTGAACCTGCACAGAGG
GCTTTATACAGGGCTGTGATGTTGGAGAAGTACAGGAACCTGGAGTTTGTGGATAGCTCTTTAAAAATCC
ATGATGGAGTTCTCATCAACCAGGCACAGTATTACAGGAGAAGTATCCACACAGGGACGTTGCAAAGA
CATAAAAGTCATCACATTGGAGATTTTTGCTTCCAGAAAATGAAGAAAGATATTCATCACTTTGAGTTT
CAGTGGCAAGAAGTTGAAAGAAATGGCCATGAAGCACCCATGACAAAAATCAAAAAGTTGACTGGTAGT
ACAGACCGAAGTGATCACAGGCATGCTGGAAACAAGCCTATTAAGATCAGCTTGGATTAAGCTTTCAT
TCGCATCTGCCTGAACCTCCACATGTTTCAGACTAAAGGGAAAAATTAGTAACCAATTGGACAAGTCTATC
GGTGCTTCTCAGCTTCCAGAAATCCCAAAGAATTTCTTGTAGGCTCAAACTCATATTTCTAATAAGTAT
GGGAAGAATTTCTCCATTCTTCATTACACAAAATACAGGAAATATGCATGAGAGAAAAACCTTGCCAA
AGTAATGAGTGTGGCAAAGCCTTAATTATAGCTCACTCTTAAGGAGACACCACATAACCCATTCAAGA
GAGAGAGAATATAAATGTGATGTATGTGGCAAGATCTTTAATCAGAAGCAATACATTGTATATCATCAT
AGATGTCACACTGGTGGAGAAAACCTTACAAGTGTAAATGAGTGTGGGAAAGACCTTCACTCAGATGTCATCC
CTTGATGCCATCGTAGACTTCATACTGGAGAGAAACCTTACAAGTGTAAATGAGTGTGGCAAGACCTTC
AGTGAGAAGTCATCCCTTAGATGCCATCGTAGACTTCATACTGGAGAGAAACCTTACAAGTGTAAATGAG
TGTGGCAAGACTTTTGGTCGAAATTCAGCCCTTGTAAATCATAAGGCAATTCATACTGGAGAGAAACCT
TACAAGTGTAAATGAGTGTGGCAAGACCTTCAGTCAGAAATCATCCCTTCAATGCCATCATATACTTAC
ACTGGAGAGAAACCTTACAAATGTGAAGAATGTGACAAATGTTTACATTGCAGATCACACCTTGAAGA
CATAGGAAGATTCATACTGGAGAGGGATCATACAAATGTAAGGTTTGTGACAAGGTTTTCCGGAGTGT
TCATACCTTGCAGAACATCAGAGAGTTCATACTGGAGAGAAACCATACAAGTGTAAATAAATGTGGCAGG
AGTTTTAGTGTGGCAAGCTTCACTTCAATACCATCATACTTCACTGGAGAGAAACCTTACACATGT
AATGAATGTGGCAAGGTTTTTAGTGAAGAGAAAACCTTGCACGTCATCATAGACTTCATGCTGGAGAG
AAACCTTACAAATGTGAAGAATGTGACAAATGTTTTCAGTCGCAGATCACACCTTGAAGACATAGGAGA
ATTCATACTGGGAAAAACCATACAAATGTAAGGTTTGTGACAAGGCTTTCCGGAGTGATTCATGCCTT
GCAAACCATACGAGAGTTCATACTGGAGAGAAACCTTACAAGTGTAAATAAATGTGCGAAGGTTTTAAT
CAAAAAGGAATCCTTGCACAACATCAGAGAGTTCATACTGGAGAGAAACCTTACAAGTGTAAATGAATGT
GGCAAGGTTTTAATCAAAAAGCAAGCCTTGCAAAACATCAGAGAGTTCATACTGCAGAGAAACCTTAC
AAGTGTAAATGAGTGTGGCAAAGCCTTACTGGACAGTCAACACTTATTCACCATCAAGCAATCCATGGG
TGTAGGAAACCTTACAAATGTAA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001202457

Insert Size: 1956 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001202457.1](#)

RefSeq Size: 2609 bp

RefSeq ORF: 1956 bp

Locus ID: 125893

UniProt ID: [Q0VGE8](#)

Cytogenetics: 19q13.41

MW: 75.7 kDa

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (3) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein.