

## Product datasheet for SC331500

### MBD1 (NM\_001204141) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** MBD1 (NM\_001204141) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** MBD1  
**Synonyms:** CXXC3; PCM1; RFT  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331500 representing NM\_001204141.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTGAGGACTGGCTGGACTGCCCGGCCCTGGGCCCTGGCTGGAAGCGCCGCGAAGTCTTTTCGCAAG
TCAGGGGCCACCTGTGGACGCTCAGACACCTATTACCAGAGCCCCACAGGAGACAGGATCCGAAGCAAA
GTTGAGCTGACTCGATACCTGGGCCCTGCGTGTGATCTCACCTCTTCGACTTCAAACAAGGCATCTTG
TGCTATCCAGCCCCAAGGCCATCCCGTGGCGGTTGCCAGCAAGAAGCGAAAGAAGCCTTCAAGGCCA
GCCAAGACTCGGAAACGTCAGGTTGGACCCAGAGTGGTGAAGTCAAGGAGGAGCCCGAGGGATGAG
ACCAAGGTGACTGACACAGCCAGCTTCATTCCCTGCTCCTGGGTGCTGTGAGAATGTGGAATC
AGCTTCTCAGGGGATGGCACCCAAAGGCAGCGGCTCAAAACGTTGTGCAAGACTGTCGAGCACAGAGA
ATTGCCTTCAACCGGGAACAGAGAATGTTAAGAGCCGAGGGTGTGGAGTATGCCGGGGCTGTCAGACC
CAAGAGGATTGTGGCCATTGCCCATCTGCCTTCGCCCTCCCCGCCCTGGTCTCAGGCGCCAGTGGAAA
TGTGTCCAGCGACGTTGCTACGGGGTAAACATGCCCGCCGCAAGGGAGGCTGTGACTCCAAGATGGCT
GCCAGGCGGCGCCCCGGAGCCAGCCACTGCCTCCACCACCCCATCACAGTCCCAGAGCCACAGAG
CCGACCCAGAGCCCTGGCCCCCTCGCCACCTGCCGAGTTCATCTATTACTGTGTAGACGAGGACGAG
CTACCTTACACGAACCGCCGCGAGAACCAGCAAGTGCAGGCGCTGTGCAGCCTGCCTACGGCGGATGGAC
TGTGGCCGCTGCGACTTCTGCTGCGACAAGCCAAATTCGGGGGAGCAACCAGAAGCGCCAGAAGTGT
CGTTGGCGCCAATGCCTGCAGTTTGCCATGAAGCGGCTGCTGCCAGTGTCTGGTCAGAGTCTGAGGAT
GGGGCAGGATCGCCCCACCTTACCCTGCTCGAAAGAGGCCAGCTCTGCCCGACGGCACCATCTTGGC
CCTACCTTGAAGCCACCTTGGCTACACGCACAGCCCAACCAGACCATAACCAGGCTCCAACGAAGCAG
GAAGCAGGTGGTGGCTTTGTGCTGCCCGCCCTGGCACTGACCTTGTGTTTTACGGGAAGGCGCAAGC
AGTCTGTGCAGGTGCCGGGCCCTGTTGCAGCTTCCACAGAAGCCCTGTTGCAGGAGGCCAGTGTCT
GGCCTGAGTTGGTTGTGGCCTTACCCAGGTGAAGCAAGAGAAGGCGGATACCCAGGACGAGTGGACA
CCAGGCACAGCTGCTGACTTCTCCCGTATTGGTGCCTGGCTGCCCTAGCAAGGACGATAGCCAGGC
CTGCCTTCTGTGAAGCAAGAGCCACCTGACCCAGAGGAGGACAAGGAGGAGAACAAGGATGATTCTGCC
TCCAAATTGGCCCCAGAGGAAGGCGAGGAGGGCTGGCACCCGATGATCACGGAGATTTTCAGCCTG
GGTGAACCCGCTTCCGAGATACAGCAGTCTGGTTGCCAAGTATTACCACCTTGCTCTTACTGAGGAAA
TGCAACTGTGGTTACCACCTGTGCTGCAGGTCCGTCCTGGTTCCCTGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001204141



[View online »](#)

<b>Insert Size:</b>	1704 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001204141.1</u>
<b>RefSeq Size:</b>	2860 bp
<b>RefSeq ORF:</b>	1704 bp
<b>Locus ID:</b>	4152
<b>UniProt ID:</b>	<u>Q9UIS9</u>
<b>Cytogenetics:</b>	18q21.1
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	62.6 kDa
<b>Gene Summary:</b>	<p>The protein encoded by this gene is a member of a family of nuclear proteins related by the presence of a methyl-CpG binding domain (MBD). These proteins are capable of binding specifically to methylated DNA, and some members can also repress transcription from methylated gene promoters. This protein contains multiple domains: MBD at the N-terminus that functions both in binding to methylated DNA and in protein interactions; several CXXC-type zinc finger domains that mediate binding to non-methylated CpG dinucleotides; transcriptional repression domain (TRD) at the C-terminus that is involved in transcription repression and in protein interactions. Numerous alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Feb 2011]</p> <p>Transcript Variant: This variant (10) lacks 2 non-consecutive coding exons, and contains an alternate 3' terminal exon compared to variant 1. This results in a frame-shift in the 3' coding region, and a shorter isoform (10) with a distinct C-terminus compared to isoform 1.</p>