

## Product datasheet for **SC331504**

### HDAC9 (NM\_001204145) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** HDAC9 (NM\_001204145) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** HDAC9  
**Synonyms:** HD7; HD7b; HD9; HDAC; HDAC7; HDAC7B; HDAC9B; HDAC9FL; HDRP; MITR  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331504 representing NM\_001204145.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCACAGTATGATCAGCTCAGTGGATGTGAAGTCAGAAGTTCCTGTGGGCTGGAGCCCATCTCACCT
TTAGACCTAAGGACAGACCTCAGGATGATGATGCCCGTGGTGGACCCTGTTGTCCGTGAGAAGCAATTG
CAGCAGGAATTACTTCTTATCCAGCAGCAGCAACAAATCCAGAAGCAGCTTCTGATAGCAGAGTTTCAG
AAACAGCATGAGAACTTGACACGGCAGCACCAGGCTCAGCTTCAGGAGCATATCAAGGAACTTCTAGCC
ATAAAACAGCAACAAGAACTCCTAGAAAAGGAGCAGAACTGGAGCAGCAGAGGCAAGAAGCAAGTA
GAGAGGCATCGCAGAGAACAGCAGCTTCTCCTCTCAGAGGCAAAGATAGAGGACGAGAAAGGGCAGTG
GCAAGTACAGAAGTAAAGCAGAAGCTTCAAGAGTTCCTACTGAGTAAATCAGCAACGAAAGACACTCCA
ACTAATGGAAAAAATCATTCCGTGAGCCGCCATCCCAAGCTCTGGTACACGGCTGCCACCACACATCA
TTGGATCAAAGCTCTCCACCCCTTAGTGGAAACATCTCCATCCTACAAGTACACATTACCAGGAGCACA
GATGCAAAGGATGATTTCCCCCTTCGAAAAACTGAATCCTCAGTCAGTAGCAGTTCTCCAGGCTCTGGT
CCCAGTTCACCAAACAATGGGCCAAGTGGAAAGTGTACTGAAAATGAGACTTCGGTTTTGCCCCCTACC
CCTCATGCCGAGCAAATGGTTTTACAGCAACGCATTCTAATTCATGAAGATTCATGAACCTGCTAAGT
CTTTATACCTCTCCTTTGCCCAACATTACCTTGGGGCTTCCCGCAGTGCCATCCCAGCTCAATGCT
TCGAATTCACCTCAAAGAAAAGCAGAAGTGTGAGACGCAGACGCTTAGGCAAGGTGTTCTCCTGCCG
CAGTATGGAGGCAGCATCCCGGCATCTCCAGCCACCCTCATGTTACTTTAGAGGAAAAGCCACCCAAC
AGCAGCCACCAGGCTCTCCTGCAGCATTATTATTGAAAGAACAAATGCGACAGCAAAAGCTTCTTGTA
GCTGGTGGAGTTCCTTACATCCTCAGTCTCCCTTGGCAACAAAAGAGAGAATTTACCTGGCATTAGA
GGTACCCACAAAATGGCCCGTACAGACCCCTGAACCGAACCCAGTCTGCACCTTTGCCTCAGAGCAG
TTGGCTCAGCTGGTCATTCAACAGCAACACCAGCAATTCTGGAGAAGCAGAAGCAATACCAGCAGCAG
ATCCACATGAACAACTGCTTTGAAAATCTATTGAACAACTGAAGCAACCCAGGAGTACCTTGAGGAA
GCAGAGGAAGAGCTTCAGGGGGACCAGGCGATGCAGGAAGACAGAGCGCCCTCTAGTGGCAACAGCACT
AGGAGCGACAGCAGTGTGTGTGGATGACACACTGGGACAAGTTGGGGCTGTGAAGGTCAAGGAGGAA
CCAGTGGACAGTATGAAGATGCTCAGATCCAGGAAATGGAATCTGGGGAGCAGGCTGCTTTTATGCAA
CAGGTAATAGGCAAGATTTAGCTCCAGGATTTGTAATTAAGTCATTATCTGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001204145



[View online »](#)

<b>Insert Size:</b>	1641 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001204145.1</a>
<b>RefSeq Size:</b>	4126 bp
<b>RefSeq ORF:</b>	1641 bp
<b>Locus ID:</b>	9734
<b>UniProt ID:</b>	<a href="#">Q9UKV0</a>
<b>Cytogenetics:</b>	7p21.1
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	60.8 kDa
<b>Gene Summary:</b>	<p>Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by this gene has sequence homology to members of the histone deacetylase family. This gene is orthologous to the Xenopus and mouse MITR genes. The MITR protein lacks the histone deacetylase catalytic domain. It represses MEF2 activity through recruitment of multicomponent corepressor complexes that include CtBP and HDACs. This encoded protein may play a role in hematopoiesis. Multiple alternatively spliced transcripts have been described for this gene but the full-length nature of some of them has not been determined. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (7) has an extended 5' UTR, differs in the 3' UTR and coding sequence, and lacks an alternate in-frame exon compared to variant 1. The resulting isoform (7) has a shorter and distinct C-terminus and lacks an internal segment compared to isoform 1. Variants 7, 30 and 31 encode the same isoform (7).</p>