

## Product datasheet for SC331521

### p73 (TP73) (NM\_001204188) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** p73 (TP73) (NM\_001204188) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** p73  
**Synonyms:** P73  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331521 representing NM\_001204188.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCCCAGTCCACCGCCACCTCCCCTGATGGGGCACCACGTTTGAGCACCTCTGGAGCTCTCTGGAA
CCAGACAGCACCTACTTCGACCTTCCCCAGTCAAGCCGGGGGAATAATGAGGTGGTGGCGGAACGGAT
TCCAGCATGGACGTCTTCCACCTGGAGGGCATGACTACATCTGTCATGGCCAGTTCAATCTGCTGAGC
AGCACCATGGACCAGATGAGCAGCCGCGGCCCTCGGCCAGCCCCACACCCAGAGCAGCCGCCAGC
GTGCCACCCACTCGCCCTACGCACAACCCAGCTCCACCTTCGACACCATGTCGCCGCGCCTGTCATC
CCCTCCAACCCGACTACCCCGACCCACCACCTTTGAGGTCACCTTCCAGCAGTCCAGCAGGCCAAG
TCAGCCACTGGACGTACTCCCCGCTCTTGAAGAACTCTACTGCCAGATCGCCAAGACATGCCCCATC
CAGATCAAGGTGTCCACCCCGCCACCCAGCCACCGCCATCCGGGCCATGCCTGTTACAAGAAAGCG
GAGCACGTGACCGACGTGTAACAGCAGCTGCCCCAACACGAGCTCGGGAGGGACTTCAACGAAGGACAG
TCTGCTCCAGCCAGCCACCTCATCCGCGTGAAGGCAATAATCTCTCGCAGTATGTGGATGACCCCTGC
ACCGGCAGGCAGAGCGTGTGGTGCCTATGAGCCACCACAGGTGGGGACGGAATACACCACCATCCTG
TACAACTTCATGTGTAACAGCAGCTGTGTAGGGGGCATGAACCGCGGCCCATCCTCATCATCACC
CTGGAGATGCGGGATGGCAGGTGCTGGCCCGCGGTCTTTGAGGGCCGCATCTGCGCCTGTCTGGC
CGCGACCGAAAAGCTGATGAGGACCACTACCGGGAGCAGCAGGCCCTGAACGAGAGCTCCGCCAAGAAC
GGGGCCGCCAGCAAGCGTGCCTTCAAGCAGAGCCCCCTGCCGTCCCCGCCCTTGGTGCCGGTGTGAAG
AAGCGGGGCATGGAGACGAGGACAGTACTACCTTCAGGTGCGAGGCCGGGAGAAGTTTGAATCCTG
ATGAAGCTGAAAGAGAGCCTGGAGCTGATGGAGTTGGTGCAGCCACTGGTGGACTCCTATCGGCAG
CAGCAGCAGTCTACAGAGGCCTTTTTAAACAGGATTGGGGTGTCCAACTGCATCGAGTATTTACC
TCCAAGGTTACAGAGCATTACACCTGCAGAACCTGACCATTGAGGACCTGGGGCCCTGAAGATC
CCCAGCAGTACCGCATGACCATCTGGCGGGCCTGCAGGACCTGAAGCAGGGCCAGACTACAGCACC
GCGCAGCAGTGTCCGCTCTAGCAACGCGGCCACCATCTCCATCGCGGCTCAGGGAACTGCAGCGC
CAGCGGGTTCATGGAGCCGTGCACTTCCGCGTGCGCCACACCATCACCATCCCCAACCGCGCGGCCCA
GGCGCGGCCCTGACAGTGGGGGACTTCGGCTTCGACCTGCCGACTGCAAGGCCCGCAAGCAGCC
ATCAAGGAGGAGTTCACGGAGGCCGAGATCCACTGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001204188



[View online »](#)

<b>Insert Size:</b>	1623 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001204188.1</a></u>
<b>RefSeq Size:</b>	4866 bp
<b>RefSeq ORF:</b>	1623 bp
<b>Locus ID:</b>	7161
<b>UniProt ID:</b>	<u><a href="#">O15350</a></u>
<b>Cytogenetics:</b>	1p36.32
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>Protein Pathways:</b>	Neurotrophin signaling pathway, p53 signaling pathway
<b>MW:</b>	59.9 kDa
<b>Gene Summary:</b>	<p>This gene encodes a member of the p53 family of transcription factors involved in cellular responses to stress and development. It maps to a region on chromosome 1p36 that is frequently deleted in neuroblastoma and other tumors, and thought to contain multiple tumor suppressor genes. The demonstration that this gene is monoallelically expressed (likely from the maternal allele), supports the notion that it is a candidate gene for neuroblastoma. Many transcript variants resulting from alternative splicing and/or use of alternate promoters have been found for this gene, but the biological validity and the full-length nature of some variants have not been determined. [provided by RefSeq, Feb 2011]</p> <p>Transcript Variant: This variant (12) lacks two alternate coding exons compared to variant 1. The resulting isoform (l, also known as TA p73 zeta) lacks an alternate internal segment compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments. There are no full-length transcripts supporting this RefSeq in human; it is predicted based on partial transcript alignments and on full-length transcript support reported in PMID:12154353.</p>