

Product datasheet for SC331549

EMA (MUC1) (NM_001204294) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	EMA (MUC1) (NM_001204294) Human Untagged Clone
Tag:	Tag Free
Symbol:	EMA
Synonyms:	ADMCKD; ADMCKD1; ADTKD2; CA 15-3; CD227; EMA; H23AG; KL-6; MAM6; MCD; MCKD; MCKD1; MUC-1; MUC-1/SEC; MUC-1/X; MUC1/ZD; PEM; PEMT; PUM
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331549 representing NM_001204294. Blue=Insert sequence Red=Cloning site Green=Tag(s)
	ATGACACCGGGCACCCAGTCTCCTTTCTTCCTGCTGCTGCTCCTCACAGTGCTTACAGTTGTTACGGGT TCTGGTCATGCAAGCTCTACCCCAGGTGGAGAAAAGGAGACTTCGGCTACCCAGAGAAGTTCAGTGCCC AGCTCTACTGAGAAGAATGCTATCCCAGCACCGACTACTACCAAGAGCTGCAGAGAGACATTTCTGAAA TGGCCAGGATCTGTGGTGGTACAATTGACTCTGGCCTTCCGAGAAGGTACCATCAATGTCCACGACGTG GAGACACAGTTCAATCAGTATAAAACGGAAGCAGCCTCTCGATATAACCTGACGATCTCAGACGTCAGC GTGAGTGATGTGCCATTTCCTTTCTCTGCCCAGTCTGGGGCTGGGGCTGCCAGGCTGGGGCATCGCGCTG CTGGTGCTGGTCTGTGTTCTGGTTGCGCTGGCCATTGTCTATCTCATTGCCTTGGCTGTCCGGCGCG CGCCGAAAGAACTACGGGCAGCTGGACATCTTTCCAGCCCGGGATACCTACC
Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001204294
Insert Size:	693 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

CRIGENE EMA (MUC1) (NM_001204294) Human Untagged Clone – SC331549

Reconstitution Method:	 Centrifuge at 5,000xg for 5min. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. Close the tube and incubate for 10 minutes at room temperature. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM 001204294.1</u>
RefSeq Size:	1091 bp
RefSeq ORF:	693 bp
Locus ID:	4582
UniProt ID:	<u>P15941</u>
Cytogenetics:	1q22
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	24.6 kDa
Gene Summary:	This gene encodes a membrane-bound protein that is a member of the mucin family. Mucins are O-glycosylated proteins that play an essential role in forming protective mucous barriers on epithelial surfaces. These proteins also play a role in intracellular signaling. This protein is expressed on the apical surface of epithelial cells that line the mucosal surfaces of many different tissues including lung, breast stomach and pancreas. This protein is proteolytically cleaved into alpha and beta subunits that form a heterodimeric complex. The N-terminal alpha subunit functions in cell-adhesion and the C-terminal beta subunit is involved in cell signaling. Overexpression, aberrant intracellular localization, and changes in glycosylation of

this protein have been associated with carcinomas. This gene is known to contain a highly polymorphic variable number tandem repeats (VNTR) domain. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Feb 2011]

Transcript Variant: This variant (18) has multiple differences compared to variant 1, one of which results in a translational frameshift in the central coding region. The encoded isoform (18) is shorter than isoform 1.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US