

Product datasheet for **SC331556**

Amyloid Precursor Protein (APP) (NM_001204302) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Amyloid Precursor Protein (APP) (NM_001204302) Human Untagged Clone
Tag:	Tag Free
Symbol:	Amyloid Precursor Protein
Synonyms:	AAA; ABETA; ABPP; AD1; alpha-sAPP; APPI; CTFgamma; CVAP; PN-II; PN2; preA4
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331556 representing NM_001204302.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTGCCCGGTTTGGCACTGCTCCTGCTGGCCGCCTGGACGGCTCGGGCGCTGGAGGTACCCACTGAT
GGTAATGCTGGCTGCTGGCTGAACCCAGATTGCCATGTTCTGTGGCAGACTGAACATGCACATGAAT
GTCCAGAATGGGAAGTGGGATTCAGATCCATCAGGGACCAAAACCTGCATTGATACCAAGGAAGGCATC
CTGCAGTATTGCCAAGAAGTCTACCCTGAACCTGCAGATCAACCAATGTGGTAGAAGCCAACCAACCATG
ACCATCCAGAAGCTGGTGAAGCGGGCCGCAAGCAGTGAAGACCCATCCCCACTTTGTGATTCCTTAC
CGCTGCTTAGTTGGTGTGTTTGAAGTGTGCCCCTTCTCGTTTCTGACAAGTGAATTTCTTACACCAG
GAGAGGATGGATGTTTGCGAAACTCATCTTCACTGGCACACCGTCGCAAGAGACATGCAGTGAGAAG
AGTACCAACTTGCATGACTACGGCATGTTGCTGCCCTGCGGAATTGACAAGTTCGAGGGGTAGAGTTT
GTGTGTTGCCACTGGCTGAAGAAAGTGAATGTGGATTCTGCTGATGCGGAGGAGGATGACTCGGAT
GTCTGGTGGGCGGAGCAGACACAGACTATGCAGATGGGAGTGAAGACAAAGTAGTAGAAGTAGCAGAG
GAGGAAGAAGTGGCTGAGGTGGAAGAAGAAGAAGCCGATGATGACGAGGACGATGAGGATGGTGTGAG
GTAGAGGAAGAGGCTGAGGAACCTACGAAGAAGCCACAGAGAGAACCACCAGCATTGCCACCACCACC
ACCACCACCACAGAGTCTGTGGAAGAGGTGGTTGAGAGAGTGTGCTCTGAACAAGCCGAGACGGGGCCG
TGCCGAGCAATGATCTCCCGCTGGTACTTTGATGTGACTGAAGGGAAGTGTGCCCCATTCTTTACGGC
GGATGTGGCGGCAACCGGAACAACCTTTGACACAGAAGAGTACTGCATGGCCGTGTGTGGCAGCGCCATT
CCTACAACAGCAGCCAGTACCCCTGATGCCGTTGACAAGTATCTCGAGACACCTGGGGATGAGAATGAA
CATGCCATTTCCAGAAAGCCAAAGAGAGGCTTGAGGCCAAGCACCGAGAGAGAAATGTCCAGGTCATG
AGAGAATGGGAAGAGGCAGAACGTCAAGCAAAGAACTTGCCCTAAAGCTGATAAGAAGGCAGTTATCCAG
CATTTCCAGGAGAAAGTGAATCTTTGGAACAGGAAGCAGCCAACGAGAGACAGCAGCTGGTGGAGACA
CACATGGCCAGAGTGAAGCCATGCTCAATGACCGCCCGCCCTGGCCCTGGAGAACTACATCACCCT
CTGCAGGCTGTTCTCCTCGGCCCTCGTCACGTGTTCAATATGCTAAAGAAGTATGTCCGCGCAGAACAG
AAGGACAGACAGCACACCCTAAAGCATTTCGAGCATGTGCGCATGGTGGATCCCAAGAAAGCCGCTCAG
ATCCGGTCCCAGGTTATGACACACCTCCGTGTGATTTATGAGCGCATGAATCAGTCTCTCTCCCTGCTC
TACAACGTGCCTGCAGTGGCCGAGGAGATTGAGGATGAAGTTGATGAGCTGCTTCAGAAAGAGCAAAAC
TATTCAGATGACGCTTTGGCCAACATGATTAGTGAACCAAGGATCAGTTACGAAACGATGCTCTCATG
CCATCTTTGACCGAAACGAAAACACCCGTGGAGCTCCTCCCGTGAATGGAGAGTTCAGCCTGGACGAT
CTCCAGCCGTGGCATTCTTTGGGGCTGACTCTGTGCCAGCCAACACAGAAAACGAAGGTTCTGGGTTG
ACAAATATCAAGACGGAGGAGATCTCTGAAGTGAAGATGGATGCAGAATCCGACATGACTCAGGATAT
GAAGTTCATCATCAAAAATTGGTGTCTTTGCAGAAGATGTGGGTCAAACAAAGGTGCAATCATTGGA
CTCATGGTGGGCGGTGTTGTCATAGCGACAGTATCGTCATCACCTTGGTGTGCTGAAGAAGAAACAG
TACACATCCATTCATCATGGTGTGGTGGAGTTGACGCCGCTGTCAACCCAGAGGAGCGCCACCTGTCC
AAGATGCAGCAGAACGGCTACGAAAATCCAACCTACAAGTTCTTTGAGCAGATGCAGAACTAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001204302

Insert Size: 2202 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204302.1](#)

RefSeq Size: 3537 bp

RefSeq ORF: 2202 bp

Locus ID: 351

UniProt ID: [P05067](#)

Cytogenetics: 21q21.3

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Alzheimer's disease

MW: 82.9 kDa

Gene Summary: This gene encodes a cell surface receptor and transmembrane precursor protein that is cleaved by secretases to form a number of peptides. Some of these peptides are secreted and can bind to the acetyltransferase complex APBB1/TIP60 to promote transcriptional activation, while others form the protein basis of the amyloid plaques found in the brains of patients with Alzheimer disease. In addition, two of the peptides are antimicrobial peptides, having been shown to have bacteriocidal and antifungal activities. Mutations in this gene have been implicated in autosomal dominant Alzheimer disease and cerebroarterial amyloidosis (cerebral amyloid angiopathy). Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Aug 2014]

Transcript Variant: This variant (9) lacks two alternate in-frame exons compared to variant 1. The resulting isoform (i, also known as L-APP733) has the same N- and C-termini but is shorter compared to isoform a. No full-length transcript is available for this transcript; however, it is supported by PMID:1429732 and other publications.