

Product datasheet for **SC331580**

RGS6 (NM_001204420) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: RGS6 (NM_001204420) Human Untagged Clone
Tag: Tag Free
Symbol: RGS6
Synonyms: GAP; HA117; S914
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331580 representing NM_001204420.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCTCAAGGATCCGGGGATCAAAGAGCAGTGGGGGTTGCTGACCCAGAGGAGAGTTCTCCAAACATG
ATCGTTTACTGCAAAATTGAAGACATCATTACAAAGATGCAAGATGACAAGACAGGGGTGTGCCCATC
AGAACAGTCAAGAGCTTTCTCTCCAAAATCCCCAGTGTGTCACAGGACTGACATTGTGCAGTGGCTT
ATGAAGAACCTTTCCATTGAGGACCCAGTTGAAGCAATACACTTGGGGAGCCTTATCGTGGCCAGGGC
TACATCTTTCCAATCTCAGACCATGTTCTCACCATGAAGGATGATGGCACCTTTTATCGTTTCCAGGCT
CCGTACTTCTGGCCTTCGAACTGCTGGGAACCTGAAAACACTGACTATGCCATCTATCTCTGTAAGAGG
ACAATGCAAAATAAAGCAAGGCTGGAACCTGGCAGATTATGAAGCAGAAAACCTTAGCAAGACTCCAGAGG
GCCTTTGCGAGGAAGTGGGAATTCATCTTTATGCAAGCAGAAGCACAAGTAAAGATTGACCGGAAAAAA
GACAAGACAGAAAAGGAAAAATTTGGATAGTCAAGAACGAGCCTTTTGGGATGTCCACAGGCCTGTGCCA
GGCTGTGTGAACACAACAGAAATGGATATCCGAAAATGTCGACGTTTGAAGAATCCACAAAAGGTTAAA
AAGTCCGTGTATGGCGTGACTGAAGAGTCCCAGGCACAGAGCCCGGTGCATGTACTCAGCCAACCAATC
AGGAAAACAACAAAAGAGGACATCCGAAAACAGATAACATTTTGAACGCACAGATCGACAGACATTGT
TTGAAAATGTCCAAAGTGGCTGAAAGCAAAGAGCCCAGCCAACAGCGAGTAAAAAGATGGGGCTTCTCT
TTCGATGAGATATTGAAGGACCAGTGGGGCGGGACCAGTTTCTACGATTCTGGAGTCCGAATTCAGT
TCAGAAAACCTCAGGTTCTGGCTGGCTGTCCAAGATCTTAAGAAAACAACCCCTACAGGATGTGGCCAAG
AGGGTAGAAGAAATCTGGCAAGAGTTTCTGGCTCCAGGGGCTCCAAGTGAATCAACCTGGATTCTCAC
AGCTATGAGATAACAGTCAAAATGTCAAAGATGGAGGGAGATATACATTTGAAGACGCCAGGAGCAC
ATCTACAAGCTGATGAAGAGTGACAGCTATGCCCGCTTCTCCGGTCAAATGCTTACCAGGATTTGCTG
CTGGCCAAGAAGAAGGAAAGTCGCTGGCGGGCAAGCGCCTCACGGGCCTGATGCAGTCTCTCTGA
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001204420
Insert Size: 1308 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001204420.1
RefSeq Size:	5639 bp
RefSeq ORF:	1308 bp
Locus ID:	9628
UniProt ID:	P49758
Cytogenetics:	14q24.2
Protein Families:	Druggable Genome
MW:	50.2 kDa
Gene Summary:	<p>This gene encodes a member of the RGS (regulator of G protein signaling) family of proteins, which are defined by the presence of a RGS domain that confers the GTPase-activating activity of these proteins toward certain G alpha subunits. This protein also belongs to a subfamily of RGS proteins characterized by the presence of DEP and GGL domains, the latter a G beta 5-interacting domain. The RGS proteins negatively regulate G protein signaling, and may modulate neuronal, cardiovascular, lymphocytic activities, and cancer risk. Many alternatively spliced transcript variants encoding different isoforms with long or short N-terminal domains, complete or incomplete GGL domains, and distinct C-terminal domains, have been described for this gene, however, the full-length nature of some of these variants is not known.[provided by RefSeq, Mar 2011]</p> <p>Transcript Variant: This variant (6) lacks 2 non-consecutive in-frame coding exons compared to variant 1, resulting in a shorter isoform (6, also known as RGS6Lalpha2(-GGL)) missing 2 internal protein segments, compared to isoform 1.</p>