

Product datasheet for **SC331598**

KIFAP3 (NM_001204514) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIFAP3 (NM_001204514) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIFAP3
Synonyms:	dj190116.1; FLA3; KAP-1; KAP-3; KAP3; SMAP; Smg-GDS
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331598 representing NM_001204514.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCGATTGCCCTTGGCAACACTTACCGGAGAAGTCTTTCAGATTGAGAACTGTTTCTGACTCCAA
GCAAGCTTATGTGAAGAGAAAAAGAAAAATCAAGCAAGCCTAAAGATCCACCTCCTTTTGAAGGAATG
GAGATTGATGAAGTTGCTAACATTAATGACATGGATGAATATATTGAGTTATTATGAAGATATTCCT
GACAAAGTTCGGGGTTCTGCTTTGATCCTGCAGCTTGTCTGAAAATCCTGATAACTTGAAGAAGCTACTA
TTGAATGAAACTGCCCTTGGTGCAATAGCAAGGTCCTGAGAGAAGACTGGAAGCAAAGTGTGAGTTA
GCTACAAACATAATTTACATCTTTTTTTGTTTCTCCAGCTTTTCTCAATTCATGGACTTATTACTCAC
TATAAAATTGGAGCTCTGTGTATGAATATTATTGATCATGAGTTAAAAAGACATGAGCTTTGGCAAGAA
GAACTCTCAAAGAAGAAGAAAGCTGTTGATGAAGACCCTGAAAACCAAACCTTGAGAAAGGATTATGAA
AAAACCTTTAAAAAGTACCAGGGGCTTGTGGTAAAACAGGAACAGCTATTACGAGTTGCTCTTTATTTG
CTTCTGAATCTTGTGAGGATACTCGTACCGAACTGAAAATGAGGAACAAGAACATAGTTCACATGTTG
GTGAAAGCCCTTGATCGGACAATTTGAGCTGCTAATTTTAGTTGTGTCATTCTTGAAGAACTCAGC
ATTTTTATGGAGAATAAAAATGATATGGTGAAAATGGATATTGTTGAAAACTGGTGAAAATGATACCT
TGTGAGCATGAAGACCTGCTGAATATCACCTCCGACTTTTACTAACCTATCCTTTGACACAGGACTG
AGGAATAAGATGGTACAAGTTGGACTGCTTCCCAAGCTCACTGCCTCCTAGGCAATGACAACTACAAA
CAAATAGCAATGTGTGTTCTTTACCACATAAGCATGGATGACCGCTTTAAATCAATGTTTGACATACCT
GACTGTATACCACAGTTAATGAAGATGCTGTTGAATGTTGAGATGAACGAATTGACTTGGAACTCATT
TCTTTCTGCATTAATCTTGTCTGCTAACAAAAGAAATGTACAGCTTATCTGTGAAGGAAATGGGCTGAAG
ATGCTCATGAAGAGGGCTCTGAAGTTTAAAGGATCCATTGCTGATGAAAATGATTAGAAACATTTCTCAG
CATGATGGACCAACTAAAAATCTGTTTATTGATTATGTTGGGGACCTTGCAGCCCAGATCTCTAATGAT
GAAGAAGAGGAGTTTGTGATTGAATGTTTGGGAACTCTTGCAAACCTTGACCATTCCAGACTTAGACTGG
GAATTGGTTCTTAAAGAATATAAGTTGGTCCATACCTCAAGGATAAACTAAAACCAAGGTGCTGCAGAA
GATGATCTTGTTTAGAAGTGGTTATAATGATTGGAAGTGTATCCATGGATGACTCTTGTGCTGCATTG
CTAGCCAAATCTGGCATAATCCCTGCCTCATTGAATTGCTAAATGCTCAACAAGAAGATGATGAATTT
GTGTGTCAGATAATTTATGTCTTCTACCAGATGGTTTTCCACCAAGCCACAAGAGACGTCAATCAAG
GAAACACAGGCTCCAGCATATCTCATAGACCTAATGCATGATAAGAATAATGAAATCCGAAAGGTCTGT
GATAATACATTAGATATTATAGCGGAATATGATGAAGAATGGGCTAAGAAAATTCAGAGTGAAAAGTTT
CGCTGGCATAACTCTCAGTGGCTGGAGATGGTAGAGAGTCGTGAGATGGATGAGAGTGAGCAGTACTTG
TATGGTGTGATGCAATTGAGCCATACATTCATGAAGGAGATATTCTCGAAAGACCTGACCTTTTCTAC
AACTCAGATGGATTAATTGCCTCTGAAGGAGCCATAAGTCCCATTCTTCAATGATTACCACCTCAA
AATGGAGATGTTGTTGGGAGCATTCTTTCTGGCAGCCTTGGAAATGGATGGCTTTGGCCAACAGTT
GGCATTCTTGGACGCCCTGCCACAGCATATGGATTCCGCCTGATGAACCTTACTACTATGGCTATGGA
TCTTGA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001204514

Insert Size: 2145 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204514.1](#)

RefSeq Size: 3175 bp

RefSeq ORF: 2145 bp

Locus ID: 22920

UniProt ID: [Q92845](#)

Cytogenetics: 1q24.2

Protein Families: Druggable Genome

MW: 82.1 kDa

Gene Summary: The small G protein GDP dissociation stimulator (smg GDS) is a regulator protein having two activities on a group of small G proteins including the Rho and Rap1 family members and Ki-Ras; one is to stimulate their GDP/GTP exchange reactions, and the other is to inhibit their interactions with membranes. The protein encoded by this gene contains 9 'Armadillo' repeats and interacts with the smg GDS protein through these repeats. This protein, which is highly concentrated around the endoplasmic reticulum, is phosphorylated by v-src, and this phosphorylation reduces the affinity of the protein for smg GDS. It is thought that this protein serves as a linker between human chromosome-associated polypeptide (HCAP) and KIF3A/B, a kinesin superfamily protein in the nucleus, and that it plays a role in the interaction of chromosomes with an ATPase motor protein. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]

Transcript Variant: This variant (2) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (2) is longer at the N-terminus compared to isoform 1.