

Product datasheet for **SC331600**

KIFAP3 (NM_001204517) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KIFAP3 (NM_001204517) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIFAP3
Synonyms:	dj190116.1; FLA3; KAP-1; KAP-3; KAP3; SMAP; Smg-GDS
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331600 representing NM_001204517.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGGGGACCCCATGTTGGGAGAACGAAAAGAATGTCAAAAAATCATTGACTTAAGAGTCTCAATGCC
AACACAGATATAACTCCCTGGCAAGGAAGGTGGTTGAAGAATGTAACTCATTTCATCCTTCAAAACTA
AATGAGGTAGAACAGCTGTTGTACTATCTACAGAACCCTGATTGATTGTCAGGAAAAGAGAAAAA
GAAAAATCAAGCAAGCCTAAAGATCCACCTCCTTTTGAAGGAATGGAGATTGATGAAGTTGCTAACATT
AATGACATGGATGAATATATTGAGTTATTATGAAGATATTCCTGACAAAGTTCGGGGTTCTGCTTTG
ATCCTGCAGCTTCTCGAAATCCTGATAACTTGAAGAAGTACTATTGAATGAACTGCCCTTGGTGCA
TTAGCAAGGGTCTGAGAGAAGACTGGAAGCAAAGTGTGAGTTAGCTACAAACATAATTTACATCTTT
TTTTGTTTCTCCAGCTTTTCTCAATTTTCTGACTTATTACTCACTATAAAATTGGAGCTCTGTGTATG
AATATTATTGATCATGAGTTAAAAAGACATGAGCTTTGGCAAGAAGAACTCTCAAAGAAGAAGAAAGCT
GTTGATGAAGACCCTGAAAACCAAACCTTGAGAAAGGATTATGAAAAACCTTTAAAAAGTACCAGGGG
CTTGTGGTAAAACAGGAACAGCTATTACGAGTTGCTCTTTATTTGCTTCTGAATCTTGTGAGGACT
CGTACCGAACTGAAAATGAGGAACAAGAACATAGTTTACATGTTGGTGAAGCCCTTGATCGGGACAAT
TTTGAGCTGCTAATTTTGTGTTGTCATTCTTGAAGAACTCAGCATTTTATGGAGAATAAAAAATGAT
ATGGTGGAAATGGATATTGTTGAAAACTGGTGAAGATGACCTTGTGAGCATGAAGACCTGCTGAAT
ATCACCTCCGACTTTTACTAAACCTATCCTTTGACACAGGACTGAGGAATAAGATGGTACAAGTTGGA
CTGCTTCCCAAGCTCACTGCACTCCTAGGCAATGACAACATAAAACAAATAGCAATGTGTGTTCTTTAC
CACATAAGCATGGATGACCGCTTAAATCAATGTTTGCATACACTGACTGTATACCACAGTTAATGAAG
ATGCTGTTTGAATGTTTCCAGATGAACGAATTGACTTGGAACTCATTCTTTCTGCATTAATCTTGTGCT
AACAAAAGAAATGTACAGCTTATCTGTGAAGGAAATGGGCTGAAGATGCTCATGAAGAGGGCTCTGAAG
TTTAAGGATCCATTGCTGATGAAAATGATTAGAAACATTTCTCAGCATGATGGACCACTAAAAATCTG
TTTATTGATTATGTTGGGGACCTTGACGCCAGATCTCTAATGATGAAGAAGAGGAGTTTGTGATTGAA
TGTTTGGAACTCTTGCAAACTTGACCATTCCAGACTTAGACTGGGAATTGGTTCTTAAAGAATATAAG
TTGTTCCATACCTCAAGGATAAACTAAAACAGGTGCTGCAGAAGATGATCTTGTTTTAGAAGTGTT
ATAATGATTGAACTGTATCCATGGATGACTCTTGTGCTGATTGCTAGCCAAATCTGGCATAATCCCT
GCACTCATTGAATTGCTAAATGCTCAACAAGAAGATGATGAATTTGTGTGTCAGATAATTTATGTCTTC
TACCAGATGGTTTTCCACCAAGCCACAAGAGACGTCATAATCAAGGAAACACAGGCTCCAGCATATCTC
ATAGACCTAATGCATGATAAGAATAATGAAATCCGAAAGGTCTGTGATAATACATTAGATATTATAGCG
GAATATGATGAAGATGGGCTAAGAAAATTCAGAGTGAAGTTTCGCTGGCATAACTCTCAGTGGCTG
GAGATGGTAGAGAGTCGTCAGATGGATGAGAGTGAAGTACTTGTATGGTGTGATGATCGAATTGAGCCA
TACATTATGAAGGAGATATTCTCGAAAGACCTGACCTTTTCTCAAACTCAGATGGATTAATTGCCTCT
GAAGGAGCCATAAGTCCCGATTTCTTCAATGATTACCACCTTCAAAATGGAGATGTTGTTGGGCAGCAT
TCATTTCTGGCAGCCTTGGAAATGGATGGCTTTGGCCAACAGTTGGCATTCTTGGACGCCCTGCCACA
GCATATGGATTCCGCCCTGATGAACCTTACTACTATGGCTATGGATCTTGA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001204517

Insert Size: 2259 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204517.1](#)

RefSeq Size: 4120 bp

RefSeq ORF: 2259 bp

Locus ID: 22920

UniProt ID: [Q92845](#)

Cytogenetics: 1q24.2

Protein Families: Druggable Genome

MW: 86.7 kDa

Gene Summary: The small G protein GDP dissociation stimulator (smg GDS) is a regulator protein having two activities on a group of small G proteins including the Rho and Rap1 family members and Ki-Ras; one is to stimulate their GDP/GTP exchange reactions, and the other is to inhibit their interactions with membranes. The protein encoded by this gene contains 9 'Armadillo' repeats and interacts with the smg GDS protein through these repeats. This protein, which is highly concentrated around the endoplasmic reticulum, is phosphorylated by v-src, and this phosphorylation reduces the affinity of the protein for smg GDS. It is thought that this protein serves as a linker between human chromosome-associated polypeptide (HCAP) and KIF3A/B, a kinesin superfamily protein in the nucleus, and that it plays a role in the interaction of chromosomes with an ATPase motor protein. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]

Transcript Variant: This variant (4) differs in the 5' UTR and coding sequence compared to variant 1. The resulting isoform (4) is shorter at the N-terminus compared to isoform 1.