

Product datasheet for **SC331605**

Cadherin 16 (CDH16) (NM_001204746) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cadherin 16 (CDH16) (NM_001204746) Human Untagged Clone
Tag:	Tag Free
Symbol:	Cadherin 16
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331605 representing NM_001204746.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTCCCTGCCTGGCTGTGGCTGCTTTGTGTCTCCGTCCCCAGGCTCTCCCCAAGGCCAGCCTGCA
GAGCTGTCTGTGGAAGTCCAGAAAATATGGTGGAAATTTCCCTTTATACCTGACCAAGTTGCCGCTG
CCCCGTGAGGGGGCTGAAGGCCAGATCGTGCTGTACAGGGGACTCAGGCAAGGCAACTGAGGGCCATTT
GCTATGGATCCAGATTCTGGCTTCCTGCTGGTGACCAAGGCCCTGGACCGAGAGGAGCAGGCAGAGTAC
CAGCTACAGGTACCCCTGGAGATGCAGGATGGACATGTCTTGTGGGGTCCACAGCCTGTGCTTGTGCAC
GTGAAGGATGAGAATGACCAGGCCTCAGGCCACCAGGCCACTGCCACCGTGAAGTCTCCATCATAGAG
AGCACCTGGGTGTCCCTAGAGCCTATCCACCTGGCAGAGAATCTCAAAGTCTATACCCGCACCACATG
GCCAGGTACACTGGAGTGGGGGTGATGTGCACTATCACCTGGAGAGCCATCCCCGGGACCCCTTTGAA
GTGAATGCAGAGGGAAACCTCTACGTGACCAGAGAGCTGGACAGAGAAGCCAGGCTGAGTACCTGCTC
CAGGTGCGGGCTCAGAATCCCATGGCGAGGACTATGCGGCCCTCTGGAGCTGCACGTGCTGGTATG
GATGAGAATGACAACGTGCCTATCTGCCCTCCCCGTGACCCACAGTCAGCATCCCTGAGCTCAGTCCA
CCAGGTACTGAAGTGACTAGACTGTCAGCAGAGGATGCAGATGCCCCGGCTCCCCAATCCCACGTT
GTGTATACAGTCCCTGAGCCCTGAGCCTGAGGATGGGGTAGAGGGGAGAGCCTTCCAGGTGGACCCACT
TCAGGCAGTGTGACGCTGGGGGTGCTCCCACTCCGAGCAGGCCAGAACATCCTGCTTCTGGTGTGGCC
ATGGACCTGGCAGGCGCAGAGGGTGGCTTCAGCAGCACGTGTGAAGTGAAGTCGAGTCACAGATATC
AATGATCACGCCCTGAGTTCATCACTTCCCAGATTGGGCCTATAAGCCTCCCTGAGGATGTGGAGCCC
GGGACTCTGGTGGCCATGCTAACAGCCATTGATGCTGACCTCGAGCCCGCCTTCCGCCTCATGGATTTT
GCCATTGAGAGGGGAGACACAGAAGGGACTTTTGGCCTGGATTGGGAGCCAGACTCTGGGCATGTTAGA
CTCAGACTCTGCAAGAACCTCAGTTATGAGGCAGCTCCAAGTCATGAGGTGGTGGTGGTGCAGAGT
GTGGCGAAGCTGGTGGGGCCAGGCCAGGCCCTGGAGCCACCGCCACGGTGACTGTGCTAGTGGAGAGA
GTGATGCCACCCCCAAGTTGGACCAGGAGAGCTACGAGGCCAGTGTCCCATCAGTGCCCCAGCCGGC
TCTTCTGCTGACCATCCAGCCCTCCGACCCATCAGCCGAACCCTCAGGTTCTCCCTAGTCAATGAC
TCAGAGGGCTGGCTCTGCATTGAGAAATCTCCGGGGAGGTGCACACCGCCAGTCCCTGCAGGGCGCC
CAGCCTGGGGACACCTACACGGTGCTTGTGGAGGCCAGGATACAGATGAGCCGAGACTGAGCGCTTCT
GCACCCCTGGTGATCCACTTCTAAAGGCCCTCTGCCACAGCCCTGACTCTTGCCCTGTGCCCTCC
CAATACCTCTGCACACCCGCCAAGACCATGGCTTGATCGTGAGTGACCCAGCAAGGACCCCGATCTG
GCCAGTGGGCACGGTCCCTACAGCTTACCCTTGGTCCCAACCCACGGTGCAACGGGATTGGCGCTC
CAGACTCTCAATGGTTCCCATGCCTACCTACCTTGGCCCTGCATTGGGTGGAGCCAGTGAACACATA
ATCCCCGTGGTGGTCAGCCACAATGCCAGATGTGGCAGCTCCTGGTTCGAGTGATCGTGTGTGCTGC
AACGTGGAGGGGAGTGCATGCGCAAGGTGGGCCGATGAAGGGCATGCCACGAAGCTGTCGGCAGTG
GGCATCCTTGTAGGCACCCTGGTAGCAATAGGAATCTTCTCATCCTCATTTTACCCCACTGGACCATG
TCAAGGAAGAAGGACCCGGATCAACCAGCAGACAGCGTGGCCCTGAAGGCGACTGTCTGA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001204746

Insert Size: 2199 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001204746.1](#)

RefSeq Size: 2697 bp

RefSeq ORF: 2199 bp

Locus ID: 1014

UniProt ID: [O75309](#)

Cytogenetics: 16q22.1

Protein Families: Transmembrane

MW: 79.2 kDa

Gene Summary: This gene is a member of the cadherin superfamily, genes encoding calcium-dependent, membrane-associated glycoproteins. Mapped to a previously identified cluster of cadherin genes on chromosome 16q22.1, the gene localizes with superfamily members CDH1, CDH3, CDH5, CDH8 and CDH11. The protein consists of an extracellular domain containing 6 cadherin domains, a transmembrane region and a truncated cytoplasmic domain but lacks the prosequence and tripeptide HAV adhesion recognition sequence typical of most classical cadherins. Expression is exclusively in kidney, where the protein functions as the principal mediator of homotypic cellular recognition, playing a role in the morphogenic direction of tissue development. Alternatively spliced transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Mar 2011]

Transcript Variant: This variant (4) lacks an in-frame segment in the 5' coding region, as compared to variant 1. The resulting isoform (4) lacks an internal segment, as compared to isoform 1.