

## Product datasheet for **SC331638**

### TEX28 (NM\_001205201) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TEX28 (NM\_001205201) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** TEX28  
**Synonyms:** CXorf2; fTEX; MRX99; TEX28P1; TEX28P2  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331638 representing NM\_001205201.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGTTTTAAAGGCGGAACACACCAGGAGCCCCAGCGCAACCCTCCCCTCCAATGTGCCTTCATGCCGG
TCCCTGTCATCCAGCGAAGACGGCCCCAGTGGCCCTTCCAGCCTCGCAGATGGAGGCCTAGCCCACAAC
TTACAGGATAGTGTGAGGACCGCATCCTCTACCTCTCAGAGCAGCTGAGAGTGGAGAAGGCCAGTCGG
GATGGCAACACTGTGAGTACCTCAAGTGGTATCCAAAGCAGACCGGCACCAGGTGCCGCACATCCAG
CAGGCCTTTGAGAAGGTGAACCAGCGCCTCTGCCACCATCGCCAGATCGAGCACAGGCTCCACCAG
TGTCACCAGCAGCTCCAGGAGCTGGAGGAAGGCTGCAGGCCCGAGGGCTTACTGCTGATGGCAGAAAAGC
GACCCAGCCAACCTGCGAGCCACCAGTGAAGGCCCTGCTTTCAGAGCCCCCGAGCCAGGTGGGGAA
GACGGGCCGGTCAACCTGCCTCATGCCAGCAGGCCCTTCACTTGGAGAGTGCCTTCCAGAGCTTACAG
CAGGGGACGTGCTTAGAGACAGAGGATGTGGCCAGCAACAAAACCTGCTGTTGCAGAAGGTAAGGCA
GAGCTGGAAGAAGCCAAGAGGTTCCACATCAGCCTCCAGGAGTCTATCAGACCTAAAGGAGAGGTCT
CTGACTGACCTGCAGCTGTTGCTGGAGTCCCTCAGGAGGAGAAGTGTAGGCAAGCATTGATGGAAGAA
CAGGTGAATGGTCGCCTGCAGGGACAGCTGAATGAGATTTACAACCTCAAACACAATCTGGCCTGCAGC
GAAGAGAGAATGGCCTATCTATCCTATGAGAGAGCCAAGGAAATATGGGAGATCACGGAGACCTTCAAG
AGCCGAATATCCAAGCTGGAGATGCTACAGCAAGTCAACCAACTGGAGGCAGCGGAGCACCTCAAAGC
CGTCCCCCGCAGATGTTGTTCAAGTTCCTGAGTCCGCGCCTCTCACTGGCAACCGTCTCTTGGTCTTT
GTCTCCACCTTGTGCTGCTGCCCTCGTCACTGATCAGCTCACGCCTGTGCACCTGCACCATGCTGATG
CTGATCGGGCTTGGGGTCTGGCCTGGCAGAGGTGGCGGCCATCCCTGCCACAGACTGGCAGGAATGG
GTCCCCTCCAGGTGTAGACTGTACTCCAAGGACTCTGGGCCTCCAGCAGATGGACCTTAA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001205201  
**Insert Size:** 1233 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001205201.1</u>
<b>RefSeq Size:</b>	1564 bp
<b>RefSeq ORF:</b>	1233 bp
<b>Locus ID:</b>	1527
<b>UniProt ID:</b>	<u>O15482</u>
<b>Cytogenetics:</b>	Xq28
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	46.1 kDa