

Product datasheet for **SC331642**

ACSL6 (NM_001205247) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ACSL6 (NM_001205247) Human Untagged Clone
Tag:	Tag Free
Symbol:	ACSL6
Synonyms:	ACS2; FACL6; LACS2; LACS5; LACS 6
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331642 representing NM_001205247.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCTGACCTTCTCCTCGTGTCGGGGGCTCCCTCTGGCTATTCGTAGAGTTTGTCTCTCACTTCTG
GAGAAGATGCAGACACAGGAGATCCTGAGGATACTGCGACTGCCTGAGCTAGGTGACTGGGACAGTTT
TTCCGCAGCCTCTCGGCCACCACCCTCGACAGTGGCGGGGCACGGCGATCTGTGATTGGGTGGCCCT
CAGCTACTTACCCTACTATGATGATGCCCGGACCATGTACCAGGTGTTCCGCCGTGGGCTTAGCATC
TCAGGGAATGGGCCCTGTCTTGGTTTCAGGAAGCCTAAGCAGCCTTACCAGTGGCTGTCCTACCAGGAG
GTGGCCGACAGGGCTGAATTTCTGGGGTCCGGACTTCTCCAGCACAATTGTAAAGCATGCACTGATCAG
TTTATTGGTGTTTTTGCACAAAATCGGCCAGAGTGGATCATTGTGGAGCTGGCCTGCTACACATATTCC
ATGGTGGTGGTCCCGCTCTATGACACCCTGGGCCCTGGGGCTATCCGCTACATCATCAATACAGCGGAC
ATCAGCACCGTGATTGTGGACAAACCTCAGAAGGCTGTGCTTCTGCTAGAGCATGTGGAGAGGAAGGAG
ACTCCAGGCCCAAGCTGATCATCTCATGGACCATTGAAGAAGCCCTGAAAGAGAGAGGGCAGAAG
TGGGGGTGGTCATTAAGTCCATGCAGGCCGTGGAGGACTGTGGCAAGAGAATCACCAGGCTCCTGTG
CCCCCGCAGCCTGATGACCTCTCCATTGTGTGTTTCAAAAGCGGCACGACAGGGAACCCAAAAGGTGCG
ATGCTCACCCATGGGAAGTGGTGGCTGATTTCTCAGGCTTTCTGAAAGTGACAGAGAGTCAAGTGGGCT
CCCACTTGTGCGGATGTGCACATTTCTATTTGCCTTTAGCACACATGTTTGAGCGAATGGTGCAGTCT
GTCGTCTATTGCCACGGAGGGCGTGTGGCTTCTTCCAGGGAGATATCCGCCCTTCTCAGATGACATG
AAGGCTCTATGCCCCACCATCTTCCCTGTGGTCCCACGACTGCTGAACCGGATGTACGACAAGATCTTC
AGCCAGGCAAACACACCATTAAGCGCTGGCTCCTGGAGTTTGCAGCAAAGCGTAAGCAAGCCGAGGTC
CGGAGTGAATCATCAGGAATGATAGTATCTGGGATGAACTCTTCTTTAATAAGATTGAGCCAGTCTT
GGTGGGTGTGTGGGATGATTGTTACTGGAGCAGCCCCAGCATCACCAACAGTTCTGGGATTTCTCCGG
GCAGCTTAGGGTGCCAGGTTTATGAAGTTATGGCCAACTGAGTGCACAGCTGGATGTACCTTCAAC
ACTCCTGGCGACTGGACCTCAGGGCACGTAGGGGCCACTTCCCTGCAATCATATCAAGTCCGTTGAT
GTTGAGGAACTGAACTACTGGGCCGCAAGGAGAGGGAGAGATATGTGTGAGAGGACCAATGTGTTC
AAAGGCTACTTGAAGATCCAGACAGGACGAAGGAGGCCCTGGACAGCGATGGCTGGCTTACACTGGA
GACATCGGAAAATGGCTGCCGCGCAGGAACTCTTAAAATTATTGATCGGAAAAAGCATATATTTAACTT
GCTCAGGGAGAATATGTTGCACCCGAGAAGATTGAGAACATCTACATCCGGAGCCAACTGTGGCGCAA
ATCTATGTCCATGGGACAGCTTAAAGGCCTTTTTGGTAGGCATTGTTGTGCCTGACCCTGAAGTTATG
CCCTCCTGGGCCAGAAGAGAGGAATTGAAGGAACATATGCAGATCTCTGCACAAATAAGGATCTGAAG
AAAGCCATTTTGAAGATATGGTGAAGTTAGGAAAAGAAAGTGGACTCCATTCTTTTGGAGAGGTTAAA
GCCATTCACATCCATTCTGACATGTTCTCAGTTCAAAATGGCTTGCTGACACCAACTAAAAGCTAAG
AGACCTGAGCTGAGAGAGTACTTCAAAAAACAATAGAAGAGCTTTACTCAATCTCCATGTGA
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_001205247

Insert Size: 2064 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001205247.1
RefSeq Size:	6449 bp
RefSeq ORF:	2064 bp
Locus ID:	23305
UniProt ID:	Q9UKU0
Cytogenetics:	5q31.1
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Adipocytokine signaling pathway, Fatty acid metabolism, Metabolic pathways, PPAR signaling pathway
MW:	76.6 kDa
Gene Summary:	<p>The protein encoded by this gene catalyzes the formation of acyl-CoA from fatty acids, ATP, and CoA, using magnesium as a cofactor. The encoded protein plays a major role in fatty acid metabolism in the brain. Translocations with the ETV6 gene are causes of myelodysplastic syndrome with basophilia, acute myelogenous leukemia with eosinophilia, and acute eosinophilic leukemia. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Apr 2011]</p> <p>Transcript Variant: This variant (3) uses an alternate in-frame splice junction at the 3' end of an exon compared to variant 1. The resulting isoform (c) has the same N- and C-termini but lacks an internal segment compared to isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>