

Product datasheet for **SC331673**

SP100 (NM_001206702) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SP100 (NM_001206702) Human Untagged Clone
Tag: Tag Free
Symbol: SP100
Synonyms: lysp100b
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331673 representing NM_001206702.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGCAGGTGGGGCGGCGACCTGAGCACCAGGAGGCTGAATGAATGTATTTACCAGTAGCAAATGAG
ATGAACCATCTTCCTGCACACAGCCACGATTTGCAAAGGATGTTACGGAAGACCAGGGTGTAGATGAC
AGGCTGCTCTATGACATTGTATTCAAGCACTTCAAAGAATAAGGTGGAGATTTCAAATGCAATAAAA
AAGACATTTCCATTCTCGAGGGCCTCCGTGATCGTGATCTCATCAGAAATAAAATGTTTGAAGATTCT
CAAGATTCTGTAGAAACCTGGTCCCTGTACAGAGAGTGGTGTACAATGTTCTTAGTGAAGTGGAGAAG
ACATTTAACCTGCCAGTCTGGAAGCACTGTTACGCGATGTCAACATGCAGGAATACCCCGATTTAATT
CACATTTATAAAGGCTTTGAAAAATGTAATCCATGACAAATTCCTCTCCAAGAAAGTGAAGAAGAAGAG
AGGGAGGAGAGGTCTGGCCTCCAATAAGTCTTGAACAAGGAAGTGGTAAAACTCTTTTCGAAGCCTG
ACTTGGCCACCTTCGGGTTCCCATCTCATGCTGGTACAACCCCACTGAAAATGGACTCTCAGAGCAC
CCCTGTGAAACAGAACAGATAAATGCAAAGAGAAAAGATACAACCAAGTACAAAGATGATTCGCTAGGA
AGCCAACAACAAATGAACAATGTGCTCAAAGGCTGAGCCAACAGAGTCCCTGCGAACAAATGTGCTGC
CAAGTGAATAATGGGGATGCTGGAAGGGAGATGCCCTGCCGTTGCCCTGTGATGAAGAAAGCCAGAG
GCAGAGCTACACAACCATGGAATCCAAATTAATTCCTGTTCTGTGCGACTGGTGGATATAAAAAAGGAA
AAGCCATTTTCTAATTCAAAAGTTGAGTGCCAAGCCCAAGCAAGAACTCATATAACAGGCATCTGAC
ATAATAGTCATCAGCAGTGAGGACTCTGAAGGATCCACTGACGTTGATGAGCCCTTAGAAGTCTTCATC
TCAGCACCCGAGAAGTGAAGCTGTGATCAATAATGACAACCTTTAGAATCAAATGATGAAAAGGAGGGC
CAAGAAGCCACTTGCTCACGACCCAGATTGTACCAGAGCCATGGATTTAGAAAATTATCTACATTC
AGAGAAAGTTTTAAGAAAAGAGTGATAGGACAAGACCAGACTTTTCAGAATCCAGTGAAGGAGGGCG
CCCGCAGAAGCCTCGAGCGGGGCACTGAGAAGCAAGCATGGTGAAGGCTCCTATGACTTCTAGAAGT
ACATCTACTTGGAGAATACCCAGCAGGAAGAGACGTTTCAGCAGTAGTACTTTTCAGACCTGAGTAAT
GGAGAAGAGCTTCAGGAAACCTGCAGCTCATCCCTAAGAAGAGGGTCAGGTAAGAAGATTAG
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001206702
Insert Size: 1443 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001206702.1
RefSeq Size:	2015 bp
RefSeq ORF:	1443 bp
Locus ID:	6672
UniProt ID:	P23497
Cytogenetics:	2q37.1
Protein Families:	Transcription Factors
MW:	53.8 kDa
Gene Summary:	<p>This gene encodes a subnuclear organelle and major component of the PML (promyelocytic leukemia)-SP100 nuclear bodies. PML and SP100 are covalently modified by the SUMO-1 modifier, which is considered crucial to nuclear body interactions. The encoded protein binds heterochromatin proteins and is thought to play a role in tumorigenesis, immunity, and gene regulation. Alternatively spliced variants have been identified for this gene; one of which encodes a high-mobility group protein. [provided by RefSeq, Aug 2011]</p> <p>Transcript Variant: This variant (4) lacks multiple exons which represent the 3' coding and UTR regions of variant 1. Variant 4 represents use of an alternate transcription termination site resulting in a much shorter protein (isoform 4, also known as SP100A) with a distinct C-terminus.</p>