

## Product datasheet for **SC331678**

### Thyroid Peroxidase (TPO) (NM\_001206745) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Thyroid Peroxidase (TPO) (NM_001206745) Human Untagged Clone
Tag:	Tag Free
Symbol:	Thyroid Peroxidase
Synonyms:	MSA; TDH2A; TPX
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

**Fully Sequenced ORF:** >SC331678 representing NM\_001206745.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAGAGCGCTCGCTGTGCTGTCTGTACGCTGGTTATGGCTGCACAGAAGCCTTCTCCCTTCATC
TCGAGAGGGAAAGAACTCCTTTGGGAAAGCCTGAGGAGTCTCGTGTCTCTAGCGCTTTGGAGGAAAGC
AAGCGCTGGTGGACACCGCCATGTACGCCACGATGCAGAGAAACCTCAAGAAAAGAGGAATCCTTTCT
CCAGCTCAGCTTCTGTCTTTTCCAAACTCCTGAGCCAACAAGCGGAGTGATTGCCCGAGCAGCAGAG
ATAATGGAACATCAATACAAGCGATGAAAAGAAAAGTCAACCTGAAAACCTCAACAATCACAGCATCCA
ACGGATGCTTTATCAGAAGATCTGCTGAGCATCATTGCAAACATGTCTGGATGTCTCCCTTACATGCTG
CCCCAAAATGCCAAACACTTGCCTGGCGAACAATAACAGGCCCATCACAGGAGCTTGCAACAACAGA
GACCACCCAGATGGGGCGCCTCCAACACGGCCCTGGCACGATGGCTCCCTCCAGTCTATGAGGACGGC
TTCAGTCAGCCCCGAGGCTGGAACCCCGGCTTCTGTACAACGGGTTCCCACTGCCCCGGTCCGGGAG
GTGACAAGACATGTCATTCAAGTTTCAATGAGGTTGTACAGATGATGACCGCTATTCTGACCTCTG
ATGGCATGGGACAATACATCGACCACGACATCGCGTTACACCACAGAGCACCAGCAAAGCTGCCTTC
GGGGAGGGGCTGACTGCCAGATGACTTGTGAGAACAAAACCCATGTTTTCCATAACAACCTCCCGAG
GAGGCCCGGCCGCGCCGCGGCACCGCTGTCTGCCCTTCTACCGCTCTCGGCCGCTGCCGCACCGGG
GACCAAGGCGCGCTCTTTGGGAACCTGTCCACGGCCAACCCGCGCAGCAGATGAACGGGTTGACCTCG
TTCTGGACGCGTCCACCGTGTATGGCAGCTCCCCGGCCCTAGAGAGGCAGCTGCGGAACTGGACCACT
GCCGAAGGGTGTCTCCGCGTCCACGCGCGCCTCCGGGACTCCGGCCGCGCCTACCTGCCCTTCGTGCCG
CCACGCGCGCTGCGGCCTGTGCGCCGAGCCGGCATCCCCGGAGAGACCCGCGGGCCCTGCTTCTGTG
GCCGAGACGGCCGCGCCAGCGAGGTCCCTCCCTGACGGCACTGCACACGCTGTGGTGCAGGAGCAC
AACCGCTGGCCGCGCGCTCAAGGCCCTCAATGCGCACTGGAGCGCGGACGCGGTGTACCAGAGGGCG
CGCAAGTGTGGGCGCTCTGCACAGATCATCACCTGAGGGATTACATCCCCAGGACTCTGGGACCC
GAGGCCTTCCAGCAGTAGTGGTCCCTATGAAGGCTATGACTCCACCGCAACCCCACTGTGTCCAAC
GTGTTCTCCACAGCCGCTTCCGCTTCCGCCATGCCACGATCCACCCGCTGGTGAGGAGGCTGGACGCC
AGCTTCCAGGAGCACCCGACCTGCCCGGGCTGTGGTGCACCAAGGCTTTCTTACGCCATGGACATTA
CTCCGTGGAGGTTACAATGAGTGGAGGGAGTCTGCGGCCTGCCTCGCCTGGAGACCCCCGCTGACCTG
AGCACAGCCATCGCCAGCAGGAGCGTGGCCGACAAGATCTGGACTTGTACAAGCATCCTGACAACATC
GATGTCTGGTGGGAGGCTTAGCTGAAAACCTCCTCCCCAGGGCTCGGACAGGGCCCTGTTTGCCTGT
CTCATTGGGAAGCAGATGAAGGCTCTGCGGACGGTACTGGTTTTGGTGGGAGAACAGCCACGTCTTC
ACGGATGCACAGAGGCGTGAAGGCTCTGCGGACGGTACTGGTTTTGGTGGGAGAACAGCCACGTCTTC
ACCAGGGTGCCTATGGATGCCTTCCAAGTCGGCAAATCCCCGAAGACTTTGAGTCTTGACAGCATC
ACTGGCATGAACCTGGAGGCTGGAGGGAAACCTTCTCAAGACGACAAGTGTGGCTTCCAGAGAGC
GTGGAGAATGGGACTTTGTGCACTGTGAGGAGTCTGGGAGGCGCGTGTGGTGTATTCTGCCGGCAC
GGGTATGAGTCCAAGCCGGGAGCAGCTCACTTGACCCAGGAAGGATGGGATTTCCAGCTCCCTC
TGCAAAGATGTGAACGAGTGTGCAGACGGTGCCACCCCCCTGCCACGCTCTGCGAGGTGCAGAAAC
ACCAAAGCGGCTTCCAGTGTCTCTGCGCGGACCCCTACGAGTTAGGAGACGATGGGAGAACCTGCGTA
GACTCCGGGAGGCTCCCTCGGGTACTTGGATCTCCATGTGCTGGTGTCTGCTGATCGGAGGCTTC
GCAGGTCTCACCTCGACGGTGATTTGCAGGTGGACACGCACTGGCACTAAATCCACACTGCCATCTCG
GAGACAGGCGGAGGAACTCCCGAGCTGAGATGCGGAAAGCACCAAGCCGTAGGGACCTCACCGCAGCGG
GCCGACGCTCAGGACTCGGAGCAGGAGTGTGGATGGAAGCCGGGATACTCACAGGCTGCCGAGA
GCCCTGA
  
```

**Restriction Sites:** SgfI-NotI  
**ACCN:** NM\_001206745  
**Insert Size:** 2631 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001206745.1</a>
<b>RefSeq Size:</b>	2980 bp
<b>RefSeq ORF:</b>	2631 bp
<b>Locus ID:</b>	7173
<b>UniProt ID:</b>	<a href="#">P07202</a>
<b>Cytogenetics:</b>	2p25.3
<b>Protein Families:</b>	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
<b>Protein Pathways:</b>	Autoimmune thyroid disease, Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling pathway, Metabolic pathways, Tyrosine metabolism
<b>MW:</b>	96.7 kDa
<b>Gene Summary:</b>	<p>This gene encodes a membrane-bound glycoprotein. The encoded protein acts as an enzyme and plays a central role in thyroid gland function. The protein functions in the iodination of tyrosine residues in thyroglobulin and phenoxy-ester formation between pairs of iodinated tyrosines to generate the thyroid hormones, thyroxine and triiodothyronine. Mutations in this gene are associated with several disorders of thyroid hormonogenesis, including congenital hypothyroidism, congenital goiter, and thyroid hormone organification defect IIA. Multiple transcript variants encoding distinct isoforms have been identified for this gene, but the full-length nature of some variants has not been determined. [provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (7) lacks a segment in the 5' UTR and an in-frame exon in the central coding region, compared to variant 1, resulting in an isoform (b, also known as TPO2) that is shorter than isoform a.</p>