

## Product datasheet for SC331683

### PKM2 (PKM) (NM\_001206796) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PKM2 (PKM) (NM_001206796) Human Untagged Clone
Tag:	Tag Free
Symbol:	PKM2
Synonyms:	CTHBP; HEL-S-30; OIP3; p58; PK3; PKM2; TCB; THBP1
Vector:	pCMV6-Entry (PS100001)
Fully Sequenced ORF:	>SC331683 representing NM_001206796. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGCAGTGGAGCTCAGAGAGAGGAGAACGGCTCCTCACGCCTGGGGCCTGCTCTTCAGAAGTCCCCAGC
GCCGTTCTTCCAGATCAGGCGGCTCTCCAGGGCACACCGTATTCAGCTCTGAGCGGTCTTTGCTAGTG
AGGCCAAGGAGCCACCCTGAGCCAAAAGGGGAGCATTATGTACCAGGAAAGCCCAACCCAGAGAACCAA
AGGACCTCAGCAGCCATGTGCAAGCCCCATAGTGAAGCCGGGACTGCCTTCATTAGACCCAGCAGCTG
CACGCAGCCATGGCTGACACATTCCTGGAGCACATGTGCCGCTGGACATTGATTACCACCCATCACA
GCCCGAAACTGGCATCATCTGTACCATTGGCCAGCTTCCCGATCAGTGGAGACGTTGAAGGAGATG
ATTAAGTCTGGAATGAATGTGGCTCGTCTGAACCTCTCTCATGGAACCTATGAGTACCATGCGGAGACC
ATCAAGAATGTGCGCACAGCCACGGAAAGCTTTGCTTCTGACCCCATCCTCTACCGGCCCGTTGCTGTG
GCTCTAGACACTAAAGGACCTGAGATCCGAACTGGGCTCATCAAGGGCAGCGGCACTGCAGAGGTGGAG
CTGAAGAAGGGAGCCACTCTCAAATCACGCTGGATAACGCCTACATGGAAAAGTGTGACGAGAACATC
CTGTGGCTGGACTACAAGAACATCTGCAAGGTGGTGAAGTGGGAGCAAGATCTACGTGGATGATGGG
CTTATTTCTCTCCAGGTGAAGCAGAAAGGTGCCGACTTCTGGTGACGGAGGTGGAAAATGGTGGCTCC
TTGGGCAGCAAGAAGGGTGTGAACCTTCTGGGCTGCTGTGGACTTGCCTGCTGTGTCGGAGAAGGAC
ATCCAGGATCTGAAGTTTGGGGTCGAGCAGGATGTTGATATGGTGTTCGCTCATTATCCGCAAGGCA
TCTGATGTCCATGAAGTTAGGAAGTCTGGGAGAGAAGGAAAGAACATCAAGATTATCAGCAAAATC
GAGAATCATGAGGGGTTTCGGAGTTTGTGAAATCCTGGAGGCCAGTGTGGGATCATGGTGGCTCGT
GGTGTCTAGGCATTGAGATTCCTGCAGAGAAGGTCTTCTTGTCTCAGAAGATGATGATTGGACGGTGC
AACCGAGCTGGGAAGCCTGTCTGTGCTACTCAGATGCTGGAGAGCATGATCAAGAAGCCCCGCCCC
ACTCGGGCTGAAGGCAGTGTGTGGCAATGCAGTCTGGATGGAGCCGACTGCATCATGCTGTCTGGA
GAAACAGCCAAAGGGGACTATCCTCTGGAGGCTGTGCCATGCAGCACCTGATAGCTCGTAGGCTGAG
GCAGCCATGTTCCACCGCAAGCTGTTTGAAGAAGTGTGCGAGCCCTCAAGTCACTCCACAGACCTCATG
GAAGCCATGGCCATGGGAGCGTGGAGGCTTCTATAAGTGTTCAGCAGCAGCTTTGATAGTTCTGACG
GAGTCTGGCAGGTCTGCTACCAGGTGGCCAGATACCGCCACAGTGCACCCATCATTGCTGTGACCCGG
AATCCCCAGACAGCTCGTCAGGCCACCTGTACCGTGGCATCTTCCCTGTGCTGTGCAAGGACCCAGTC
CAGGAGGCTGGGCTGAGGACGTGGACCTCCGGGTGAACCTTGGCCATGAATGTTGGCAAGGCCCGAGGC
TTCTTCAAGAAGGGAGATGTGGTCATTGTGCTGACCGGATGGCGCCCTGGCTCCGGCTTACCAACACC
ATGCGTGTGTTCTGTGCCGTGA

```



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001206796
<b>Insert Size:</b>	1818 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001206796.1</a>
<b>RefSeq Size:</b>	2853 bp
<b>RefSeq ORF:</b>	1818 bp
<b>Locus ID:</b>	5315
<b>UniProt ID:</b>	<a href="#">P14618</a>
<b>Cytogenetics:</b>	15q23
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Glycolysis / Gluconeogenesis, Metabolic pathways, Purine metabolism, Pyruvate metabolism, Type II diabetes mellitus
<b>MW:</b>	65.9 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein involved in glycolysis. The encoded protein is a pyruvate kinase that catalyzes the transfer of a phosphoryl group from phosphoenolpyruvate to ADP, generating ATP and pyruvate. This protein has been shown to interact with thyroid hormone and may mediate cellular metabolic effects induced by thyroid hormones. This protein has been found to bind Opa protein, a bacterial outer membrane protein involved in gonococcal adherence to and invasion of human cells, suggesting a role of this protein in bacterial pathogenesis. Several alternatively spliced transcript variants encoding a few distinct isoforms have been reported. [provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (4) represents the longest transcript and encodes the longest isoform (c).</p>