

## Product datasheet for **SC331685**

### PKM2 (PKM) (NM\_001206798) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PKM2 (PKM) (NM\_001206798) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** PKM2  
**Synonyms:** CTHBP; HEL-S-30; OIP3; p58; PK3; PKM2; TCB; THBP1  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331685 representing NM\_001206798.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGTCACCGGAAGCCCAACCCAGAGAACCAAGGACCTCAGCAGCCATGTCGAAGCCCCATAGTGAAG
CCGGGACTGCCTTCATTAGACCCAGCAGCTGCACGCAGCCATGGCTGACACATTCCTGGAGCACATGT
GCCGCCTGGACATTGATTCACCACCCATCACAGCCCGAACACTGGCATCATCTGTACCATTGTACCAT
GCGGAGACCATCAAGAAATGTGCGCACAGCCACGGAAAGCTTTGCTTCTGACCCCATCTCTACCGGCC
GTTGCTGTGGCTCTAGACACTAAAGGACCTGAGATCCGAAGTGGGCTCATCAAGGCAGCGGCACCTGCA
GAGGTGGAGCTGAAGAAGGGAGCCACTCTCAAAATCACGCTGGATAACGCCTACATGGAAAAGTGTGAC
GAGAACATCCTGTGGCTGGACTACAAGAACATCTGCAAGGTGGTGAAGTGGGCAGCAAGATCTACGTG
GATGATGGGCTTATTTCTCTCCAGGTGAAGCAGAAAGGTGCCGACTTCTGGTGACGGAGGTGAAAAT
GGTGGCTCCTTGGGCAGCAAGAAGGGTGTGAACCTTCTGGGGCTGCTGTGGACTTGCTGCTGTGTGCG
GAGAAGGACATCCAGGATCTGAAGTTTGGGGTCGAGCAGGATGTTGATATGGTGTGGTTCGCTCATTTCATC
CGCAAGGCATCTGATGTCCATGAAGTTAGGAAGTCTGGGAGAGAAGGGAAAGAACATCAAGATTATC
AGCAAAATCGAGAATCATGAGGGGGTTCGGAGGTTTGATGAAATCCTGGAGGCCAGTGATGGGATCATG
GTGGCTCGTGGTATCTAGGCATTGAGATTCTGCAGAGAAGGCTTCTTCTGCTCAGAAGATGATGATT
GGACGGTGCAACCGAGCTGGGAAGCCTGCATCTGTGCTACTCAGATGCTGGAGAGCATGATCAAGAAG
CCCCGCCCACTCGGGCTGAAGGCAGTGATGTGGCCAATGCAGTCTGGATGGAGCCGACTGCATCATG
CTGTCTGGAGAAACAGCCAAAGGGGACTATCCTCTGGAGGCTGTGCGCATGCAGCACCTGATTGCCCGT
GAGGCAGAGGCTGCCATCTACCACTTGCAATTATTTGAGGAACTCCGCCCTGGCGCCATTACCAGC
GACCCACAGAAGCCACCGCCGTGGGTGCCGTGGAGGCCCTCTTCAAGTGTGAGTGGGCCATAATC
GTCCTACCAAGTCTGGCAGGTCTGCTCACCAGGTGGCCAGATACCGCCACGTGCCCCATCATTGCT
GTGACCCGGAATCCCAGACAGCTCGTCAGGCCACCTGTACCGTGGCATCTTCCCTGTGCTGTGCAAG
GACCCAGTCCAGGAGCCTGGGCTGAGGACGTGGACCTCCGGTGAAGTGGCCATGAATGTTGGCAAG
GCCCGAGGCTTCTCAAGAAGGGAGATGTGGTCATTGTGCTGACCGGATGGCGCCTGGCTCCGGCTTC
ACCAACACCATGCGTGTTTCTGTGCCGTGA
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001206798  
**Insert Size:** 1551 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001206798.1</u>
<b>RefSeq Size:</b>	2541 bp
<b>RefSeq ORF:</b>	1551 bp
<b>Locus ID:</b>	5315
<b>UniProt ID:</b>	<u>P14618</u>
<b>Cytogenetics:</b>	15q23
<b>Protein Families:</b>	Druggable Genome
<b>Protein Pathways:</b>	Glycolysis / Gluconeogenesis, Metabolic pathways, Purine metabolism, Pyruvate metabolism, Type II diabetes mellitus
<b>MW:</b>	56.3 kDa
<b>Gene Summary:</b>	<p>This gene encodes a protein involved in glycolysis. The encoded protein is a pyruvate kinase that catalyzes the transfer of a phosphoryl group from phosphoenolpyruvate to ADP, generating ATP and pyruvate. This protein has been shown to interact with thyroid hormone and may mediate cellular metabolic effects induced by thyroid hormones. This protein has been found to bind Opa protein, a bacterial outer membrane protein involved in gonococcal adherence to and invasion of human cells, suggesting a role of this protein in bacterial pathogenesis. Several alternatively spliced transcript variants encoding a few distinct isoforms have been reported. [provided by RefSeq, May 2011]</p> <p>Transcript Variant: This variant (6) differs in the 5' UTR and coding sequence, lacks an alternate exon, and has an alternate in-frame coding exon compared to variant 4. The resulting isoform (e) has a shorter and distinct N-terminus, lacks an alternate internal segment, and has a different internal segment compared to isoform c.</p>