

Product datasheet for **SC331700**

TAX1BP1 (NM_001206901) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAX1BP1 (NM_001206901) Human Untagged Clone
Tag:	Tag Free
Symbol:	TAX1BP1
Synonyms:	CALCOCO3; T6BP; TXBP151
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331700 representing NM_001206901.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGACATCCTTTCAAGAAGTCCCATTGCAGACTTCCAACCTTTGCCCATGTCATCTTTCAAATGTGGCC
AAGAGTTACCTTCCTAATGCACACCTGGAATGTCATTACACCTTAACCCATATATTCATCCACATCCA
AAAGATTGGGTTGGTATATTCAGGTTGGATGGAGTACTGCTCGTGATTATTACACGTTTTTATGGTCC
CCTATGCCTGAACATTATGTGGAAGGATCAACAGTCAATTGTGTACTAGCATTCCAAGGATATTACCTT
CCAAATGATGATGGAGAATTTTATCAGTTCTGTTACGTTACCCATAAGGGTGAAATTCGTGGAGCAAGT
ACACCTTTCCAGTTTCGAGCTTCTTCTCCAGTTGAAGAGCTGCTTACTATGGAAGATGAAGGAAATCT
GACATGTTAGTGGTGACCACAAAAGCAGGCCTTCTTGAGTTGAAAATTGAGAAAACCATGAAAGAAAA
GAAGAACTGTTAAAGTTAATTGCCGTTCTGAAAAAGAAACAGCACAACTTCGAGAACAAGTTGGGAGA
ATGGAAGAGAACTTAACCATGAGAAAGAAAGATGTGACCAACTGCAAGCAGAACAAAAGGGTCTTACT
GAAGTAACACAAAGCTTAAAAATGGAATGAAGAGTTTAAAGAGAGTTTCAAGTATGCTACATCCAAA
GCCCATCAGCTTGAGGAAGATATTGTGTCAGTAACACATAAAGCAATTGAAAAAGAAACCGAATTAGAC
AGTTTAAAGGACAACTCAAGAAGGCACAACATGAAAGAGAACAACTTGAATGTCAGTTGAAGACAGAG
AAGGATGAAAAGGAACCTTATAAGGTACATTTGAAGAATACAGAAATAGAAAATACCAAGCTTATGTCA
GAGGTCCAGACTTTAAAAATTTAGATGGGAACAAAGAAAGCGTGATTACTCATTTCAAAGAAGAGATT
GGCAGGCTGCAGTTATGTTTGGCTGAAAAGGAAAACTGCAAAGAACTTTCTGCTTACAACCTCAAGT
AAAGAAGATACTGTTTTTAAAGGAGCAACTTCGTAAGCAGAGGAACAGGTTCAAGCAACTCGGCAA
GAAGTTGTCTTTCTGGCTAAAGAACTCAGTGATGCTGTCAACGTACGAGACAGAACGATGGCAGACCTG
CATACTGCACGCTTGAAAAACGAGAAAGTAAAAAGCAGTTAGCTGATGCAGTGGCAGAACTTAACTA
AATGCTATGAAAAAGATCAGGACAAGACTGATACACTGGAACACGAACATAAGAAGAGAAGTTGAAGAT
CTGAAACTCCGTCTTCAGATGGCTGCAGACCATTATAAAGAAAAATTTAAGGAATGCCAAAGGCTCCAA
AAACAAATAAACAACTTTAGATCAATCAGCTAATAATAATAATGTCTTCACAAAGAAAAACGGGGAAT
CAGCAGAAAGTGAATGATGCTTCAGTAAACACAGACCCAGCCACTTCTGCCTCTACTGTAGATGTAAG
CCATCACCTTCTGCAGCAGAGGCAGATTTTGACATAGTAACAAAGGGGCAAGTCTGTGAAATGACAAA
GAAATTGCTGACAAAACAGAAAAGTATAATAATGTAAACAACTTTCAGGATGAGAAAGCAAAATGC
AATAAATATGCTGATGAACCTGCAAAAATGGAGCTGAAATGGAAGAACAAGTAAAAATTGCTGAAAT
GTAAACTTGAAGTACAGGACAATTATAAAGAACTTAAAGGAGTCTAGAAAATCCAGCA
GAAAGGAAATGGAAGATGGAGCAGATGGTGCTTTTACCCAGATGAAATACAAAGGCCACCTGTCAGA
GTCCCCTCTTGGGACTGGAAGACAATGTTGTCTGCAGCCAGCCTGCTCGAACTTTAGTCGGCCTGAT
GGCTTAGAGGACTCTGAGGATAGCAAGAAGATGAGAATGTGCCTACTGCTCCTGATCCTCAAGTCAA
CATTTACGTGGGCATGGGACAGGCTTTTGTCTTCCAGCTTTGATGTTACAGAAGGTGCCCTC
TGTGAGTTAATGTTTCCTCACTATGATCAGAGCAATTTGAAGAATGTTGAAAGTCACTGGAAG
GTGTGCCCGATGTGCAGCGAGCAGTTCCCTCCTGACTATGACCAGCAGGTGTTTGAAGGCATGTGCAG
ACCCATTTTGATCAGAATGTTCTAAATTTGACTAG
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001206901

Insert Size: 2244 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001206901.1</u>
RefSeq Size:	3383 bp
RefSeq ORF:	2244 bp
Locus ID:	8887
UniProt ID:	<u>Q86VP1</u>
Cytogenetics:	7p15.2
Protein Families:	Druggable Genome
MW:	86.3 kDa
Gene Summary:	This gene encodes a HTLV-1 tax1 binding protein. The encoded protein interacts with TNFAIP3, and inhibits TNF-induced apoptosis by mediating the TNFAIP3 anti-apoptotic activity. Degradation of this protein by caspase-3-like family proteins is associated with apoptosis induced by TNF. This protein may also have a role in the inhibition of inflammatory signaling pathways. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, May 2011] Transcript Variant: This variant (3) contains an alternate 5' non-coding exon, and uses an alternate in-frame donor splice site at an internal coding exon compared to variant 1. This results in a shorter isoform (2, also known as TXBP151-L) missing an internal protein segment compared to isoform 1. Variants 2 and 3 encode the same isoform.