

Product datasheet for **SC331708**

SLC6A12 (NM_001206931) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLC6A12 (NM_001206931) Human Untagged Clone
Tag: Tag Free
Symbol: SLC6A12
Synonyms: BGT-1; BGT1; GAT2
Vector: pCMV6-Entry (PS100001)
Fully Sequenced ORF: >SC331708 representing NM_001206931.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGACGGGAAGGTGGCAGTGCAAGAGTGTGGCCTCCTGCAGTCTCCTGGGTCCCCGAGGAGGGAGAG
AAGTTGGACCAGGAAGACGAGGACCAGGTGAAGGATCGGGGCCAATGGACCAACAAGATGGAGTTTGTG
CTGTCAGTGGCCGGGAGATCATTGGGCTGGGCAATGTCTGGAGGTTTCCCTATCTCTGTACAAAAAC
GGAGGTGGAGCCTTCTTCATCCCCTACTTCATCTTCTTTGTCTGCGGCATCCCGGTGTTCTTCTCTG
GAGGTGGCGTTGGCCAATACACCAGCCAAGGGAGTGTACAGCCTGGAGGAAGATCTGCCCCCTCTTC
CAGGGCATTGGTCTGGCATCTGTGGTCATCGAGTCATATTTGAATGTCTACTACATCATCATCCTTGCC
TGGGCTCTTCTACCTGTTTCAGCTCCTCACTTCTGAGCTGCCCTGGAGGACCTGCAACAACCTTTTGG
AACACAGAGCATTGCACGGACTTCTGAACCACTCAGGAGCCGGCACAGTGACCCCAATTTGAGAATTTT
ACCTCACCTGTGATGGAATCTGGGAGAGACGAGTTCTGGGCATCACCTCGGGCATCCATGACCTGGGC
TCCCTGCGCTGGGAGCTGGCCCTGTGCCTCCTGCTCGCCTGGGTGATCTGCTATTTCTGCATCTGGAAG
GGGTCAAGTCCACAGGCAAGGTGGTTTATTTACAGCCACGTTTCCGTACCTGATGCTTGTGATTTTG
CTGATCAGAGGTGTCACCTTCCCGGAGCCTACCAGGGCATCATCTACTACTTGAAGCCAGATTTGTTT
CGCCTCAAGGACCCTCAGGTGTGGATGGATGCGGGCACCCAGATCTTCTCTCCTTTGCCATCTGCCAG
GGGTGCTGACAGCCCTGGGCAGCTACAACAAGTATACAACAACACTGCTACAAGGACTGCATCGCCCTC
TGCTTCTGAACAGTGCCACCAGCTTTGTGGCTGGGTTTGTGCTTCTCCATCCTGGGCTTCATGTCC
CAAGAGCAAGGGGTGCCATTTCTGAAGTGGCCGAGTCAGGTCCTGGGCTGGCCTTCATCGCCTTCCCC
AAGGCTGTGACTATGATGCCCTATCCCAGCTGTGGTCTGCCTGTTCTTTATCATGCTCATATTCTTA
GGGCTGGACAGCCAGTTTGTCTGTGTGGAGTGCCTGGTGACAGCCTCCATAGACATGTTCCCCAGGCAG
CTCCGGAAGAGCGGGCGGCGGAGCTCCTCATCCTCACCATCGCCGTGATGTGCTACCTGATAGGGCTT
TTCTGGTCACCGAGGGCGGGATGTACATCTTCCAGCTGTTTACTACTATGCTTCCAGTGGCATATGC
CTGCTGTTCTGTGATTGTTTGAAGTGGTCTGCATAAGCTGGGTGATGGGGCGGACCGTTTCTATGAC
AACATTGAGGACATGATTGGCTACCGCCATGGCCCTGGTGAAGATCTCCTGGCTTCTCCTGACCCCT
GGACTTTGCCTGGCCACTTCTCTTCTCCTTGGCAAGTACACCCCCCTCAAGTACAACAACGCTAT
GTGATCCCGCCCTGGGATACTCCATTGGCTGGTCTCCTGGCTCTGTCTCCATGGTCTGTGTCCCACTC
TTGCTGTCATCACCTCCTGAAGACTCGGGTCTTTTTCAGGAAGCGTCTGCGTCAGCTCATCACCCCT
GACTCCAGTCTGCCACAGCCCAAGCAACATCCCTGCTTGGATGGCAGTGTGGCCGGAACCTTTGGGCC
TCCCCAACAGGGGAAGGACTGATAGCCGGGAGAAGGAGACCCATTTGTAG
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001206931
Insert Size:	1845 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001206931.1
RefSeq Size:	3177 bp
RefSeq ORF:	1845 bp
Locus ID:	6539
UniProt ID:	P48065
Cytogenetics:	12p13.33
Protein Families:	Druggable Genome, Transmembrane
MW:	69.4 kDa
Gene Summary:	<p>Transports betaine and GABA. May have a role in regulation of GABAergic transmission in the brain through the reuptake of GABA into presynaptic terminals, as well as in osmotic regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR compared to variant 1. Variants 1, 2, 3 and 4 encode the same protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>