

## Product datasheet for **SC331768**

### DHRS7C (NM\_001220493) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** DHRS7C (NM\_001220493) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** DHRS7C  
**Synonyms:** SDR32C2  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331768 representing NM\_001220493.  
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
ATGGGAGTCATGGCCATGCTGATGCTCCCCCTGCTGCTGCTGGGAATCAGCGGCCTCCTCTTCATTTAC
CAAGAGGTGTCCAGGCTGTGGTCAAAGTCAGCTGTGCAGAACAAAGTGGTGGTGATCACCGATGCCATC
TCAGGACTGGGCAAGGAGTGTGCTCGGGTGTTCACACAGGTGGGCAAGGCTGGTGCTGTGTGGAAAG
AACTGGGAGAGGCTAGAGAACCTATATGATGCCTTGATCAGCGTGGCTGACCCCAGCAAGCAGACATTC
ACCCCAAAGCTGGTCCTGTTGGACCTCTCAGACATCAGCTGTGTCCAGATGTGGCAAAGAAGTCTCTG
GATTGCTATGGCTGTGTGGACATCCTCATCAACAATGCCAGTGTGAAGGTGAAGGGGCCTGCCATAAG
ATTTCTCTGGAGCTCGACAAAAGATCATGGATGCCAATTAATTTGGCCCATCACATTGACGAAAGCC
CTGCTTCCCAACATGATCTCCCGGAGAACAGGCCAAATCGTGTTAGTGAATAATATCCAAGGGAAGTTT
GGAATCCCCTTCCGTACGACTTACGCTGCCTCCAAGCACGCAGCCCTGGGCTTCTTTGACTGCCTCCGA
GCCGAAGTGGAGGAATACGATGTTGTCATCAGCACCGTGAGCCCGACTTTCATCCGGTCGTACCACGTG
TATCCAGAGCAAGGAACTGGGAAGCTTCCATTTGGAAATTTCTTTTCAGGAAGCTGACCTACGGCGTG
CACCCAGTAGAGGTGGCGGAGGAGGTGATGCGCACCGTGCGGAGGAAGAAGCAAGAGGTGTTATGGCC
AACCCCATCCCAAGGCCGCCGTGACGTCCGCACCTTCTCCCGGAGTCTTTTTCGCCGTGGTGGCC
TGTGGGTGAAGGAGAAGCTCAATGTCCCGGAGGAGGGTAA
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001220493  
**Insert Size:** 939 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001220493.1](#)

**RefSeq Size:** 1014 bp

**RefSeq ORF:** 939 bp

**Locus ID:** 201140

**UniProt ID:** [A6NNS2](#)

**Cytogenetics:** 17p13.1

**Protein Families:** Druggable Genome, Transmembrane

**MW:** 34.9 kDa

**Gene Summary:** Putative oxidoreductase.[UniProtKB/Swiss-Prot Function]  
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).