

## Product datasheet for SC331810

### ZNF345 (NM\_001242475) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** ZNF345 (NM\_001242475) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** ZNF345  
**Synonyms:** HZF10  
**Vector:** pCMV6-Entry (PS100001)  
**Fully Sequenced ORF:** >SC331810 representing NM\_001242475.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGAAAAACCTTACAAAACACAGCATTGAGTGTTCAAGTTTCAGAGGTGATTGGGAATGTAAAAACCAG
TTTGAGAGAAAACAGGGATCTCAGGAAGGACATTTTCAGTGAAATGATATTTACTCCTGAAGACATGCC
ACTTTTCAGTATCCAGCATCAGAGAATTCATACTGATGAGAACTCCTTGAATGTAAGGAATGTGGGAAG
GATTTTAGTTTTGTATCAGTCCTTGTTCGACATCAGCGAATTCATACTGGTGAGAAACCTTATGAATGC
AAAGAATGTGGCAAGGCCTTTGGTAGTGGTGCAAACCTTGCTTACCATCAAAGAATTCATACTGGTGAG
AAGCCTTTTGAATGTAAGAATGTGGGAAGGCCTTTGGTAGTGGTCAAACCTTACTCACCATCAGAGA
ATTCATACTGGTGAGAAACCTATGAGTGTAAAGGAATGTGGGAAAGCCTTTAGTTTTGGATCAGGCCTT
ATTCGACATCAGATCATTACAGTGGTGAGAAGCCTTATGAGTGTAAAGGAATGTGGGAAGTCCTTTAGT
TTTGAATCAGCCCTTATTCGGCATCACAGAATTCACACAGGTGAGAAACCTTATGAATGTATAGATTGT
GGTAAAGCCTTTGGCAGTGGTTCAAACCTTACTCAACATCGGCGGATTACTACTGGTGAGAAACCTTAT
GAATGCAAAGCATGTGGAATGGCCTTTAGCAGTGGTTCGGCTTACTCGGCATCAGAGAATTCATACC
GGTGAGAAACCATATATGTAATGAATGTGGTAAGGCCTTTAGTTTTGGATCAGCCCTTACTCGACAT
CAAAGAATTCATACTGGTGAGAAACCTTATGATGTAAGGAATGTGGGAAGGCTTTTAAATAGTGGCTCA
GATCTCACTCAGCATCAGAGAATTCACACTGGTGAGAAACCTATGAGTGTAAAGGAGTGTGAGAAAGCC
TTTAGAAGTGGTTCAAACCTTATTCAGCATCAAAGAATGCATACTGGAGAGAAACCTTATGAATGTAAG
GAATGTGGGAAGACCTTTAGTAGTGGTTCAGACCTTACTCAACATCACAGAATTCATACTGGTGAGAAA
CCCTATGAATGTAAGGAATGTGGGAAGGCCTTTGGTAGTGGCTCAAACCTTATCCAACACCAGCTAATC
CATACTGGTGAAAGACCCTATGAATGTAAGAATGTGGAAAGTCCCTTTAGTAGTGGTTCAGCTCTTAAT
CGGCACCAGAGAATACACACTGGTGAGAAACCTATGAATGTAAGGAGTGTGGGAAGGCTTTTTATAGT
GGCTCAAGCCTTACTCAGCATCAGAGAATTCATACAGGTGAGAAACCTTTATGAATGTAAGAAGTGTGGG
AAGCCTTATGGGAGGGATTTCAGAGTTTCAGCAACATAAGAAAAGTCATAATGGTAAGAAACTCTGCGAA
TTGGAAACTATAAATGA
  
```

**Restriction Sites:** Sgfl-Mlul  
**ACCN:** NM\_001242475  
**Insert Size:** 1467 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001242475.1</a>
<b>RefSeq Size:</b>	2945 bp
<b>RefSeq ORF:</b>	1467 bp
<b>Locus ID:</b>	25850
<b>UniProt ID:</b>	<a href="#">Q14585</a>
<b>Cytogenetics:</b>	19q13.12
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	55.4 kDa
<b>Gene Summary:</b>	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (4) differs in the 5' UTR and encodes the same protein, compared to variant 1.