

Product datasheet for **SC331848**

DLGAP1 (NM_001242762) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DLGAP1 (NM_001242762) Human Untagged Clone
Tag:	Tag Free
Symbol:	DLGAP1
Synonyms:	DAP-1; DAP-1-ALPHA; DAP-1-BETA; DAP1; DLGAP1A; DLGAP1B; GKAP; SAPAP1
Vector:	pCMV6-Entry (PS100001)



[View online »](#)

Fully Sequenced ORF: >SC331848 representing NM_001242762.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

ATGGATTTGAAGACATTAACACTGTTTAAATTCAGTTGAGATGTGGTTGGCTCTTATACATTTGGAAT
AAAGCATTTATGGCACATCTCAGGGTCCACAAGATGAATGGACAGGGTACACCCACGAGGTAAGAT
GATGAAATTCATGCCGAAGAATGCGGAGTGGCAGTTATATCAAGGCCATGGGGGATGAAGACAGTGGA
GACTCAGACACGAGTCCTAAGCCTTCTCCAAAAGTTGCTGCGCGGAGAGAAAGCTATCTCAAGCTACT
CAGCCATCCCTTACAGAACTCACCACACTCAAAATCTCCAATGAACACTACCCAAACTCCAGATCCGG
AGTCATAGTTACCTGAGGGCAGTGAGTGAAGTCTCCATCAACCGAGCCTGGACAGCCTGGACCCTGCA
GGCTTGCTCACATCACCAAAGTCCGCTCCAGGAATGAGAGCTACATGCGAGCCATGAGCACCATCAGC
CAGGTGAGCGAGATGGAAGTGAACGGGCAGTTCGAGTCCGTGTGCGAGTCCGTGTTGAGCAGCTGGAG
TCGCAGGCCGTGGAAGCCTGGACCTGCCATGCCCGGCTGCTCCGCATGCGGAGCCACAGCTATGTG
CGGGCCATTGAGAAAGGCTGCTCCAGGACGACGAGTGCCTGTCCTGAGGTCGCTCGCCGCCGCGC
ACCACCACCCTTAGGACCATCCAGAGCAGCACGGCCAGAGTAGCACAGAGTCAGCCCAGGATGCC
TACATGGACGGACAGGGCCAGCAGGAGATATTATCAGCCAGTCTGGACTCAGCAACTCCACCGAGAGC
CTGGACAGTATGAAGGCTCTGACAGCCGCCATCGAAGCTGCAAACGCCAGATCCATGGCCCTGCCAGT
CAACACATGGGAATAACACTGCCACCCTCACCACCACGACTACCATAGCCACCCTCACCACGGAGGAC
AGGAAGAAGGACCACTTTAAGAAAAATCGATGCCTGTCTATCGGGATACAGGTGGATGATGCTGAAGAA
CCTGACAAAACAGGGGAGAATAAAGCACCCAGTAAAGTTCAGTCCGTGGGAGTGAAGTAGAAGAAGAG
AAGTGCTTCCGCAGGTTCACTCGATCCAACAGTGTGACGACAGCAGTACAGGCCGACCTGGACTTCCAT
GATAATCTGAAAATTCTCTGGAATCTATAGAGGACAATTCGTGTCTGGCCCCATGGCCAGACAGTTC
TCCCGCATGCCAGCACCTCCACAGTCAAGATTCAGGGCTCAGGAAACCATACCATGCTGTGCCGCC
GATGATGACTTTGACACGGATTTGACCCTCTATTCTGCCTCCTCCGGACCCCTGGATTGACTATC
ACTGAAGACCCTCTGGAGCCGTGCAAAGGTCAAGTGTGCCACCGGATGGCCACTGGTTCCTGAAGCTT
CTCCAGGCAGAGCGAGACCCATGGAGGGTGGTCAACAGATGGAGCGGGAAGAACGGAAAAACAAC
CTGCCGAAGACATTCTAGGAAAAATCCGAACCGCAGTGGGCAGTGCCCAACTTCTCATGGCCAGAAA
TTCTACCAGTTCAGAGAACTGTGTGAAGAAAACCTGAATCCTAATGCTCATCCAAGACCCACCTCCAG
GATTTGGCGGGTTTTGGGACATGCTGCAGTTGCCATAGAAAATATTAGTATGAAATTTGATGAACCT
CATCAGTTAAAGCCAATAATTGAAACAGATGGATCCTCTTGACAAGAAGGAGAGAAGGGCCCTCCT
CCAGTGCCAAAGAAGCCGCGAAGGGCCCCGCGCGCTGATCCGGGAGCGCTCGCTGGAGAGCTCGCAG
CGCCAGGAGGCCGCAAGCGCTGATGGCCGCCAAGCGCCGCGCTCCGTCCGCCAGAACTCGGCCACC
GAGAGCGCCGAGAGCATCGAGATCTACATCCCGAGGCGCAGACCCGGCTGA
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001242762

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242762.1](#)

RefSeq Size: 5214 bp

RefSeq ORF: 1986 bp

Locus ID: 9229

UniProt ID: [O14490](#)

Cytogenetics: 18p11.31

MW: 74.1 kDa

Gene Summary: Part of the postsynaptic scaffold in neuronal cells.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (4) lacks a portion of the 5' UTR and 5' coding region, initiates translation at an alternate start codon, and has multiple coding region differences, compared to variant 1. The encoded isoform (4) has a distinct N-terminus and is shorter than isoform 1.
Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.